

mercymedical 
EQUIPMENT COMPANY

enovis™



AIRCAST®

PROCARE®



Federal Supply Schedule (FSS)

Product Catalog

Mercy Medical Equipment Company
1303 S Comal, San Antonio TX 78207
SAM UEI: SN3TFMPA2TH6 | FSS Contract # V797D-50479

To Order:
Email: info@emergymedical.com
Phone: 210-224-9714 | Fax: 210-224-9812

TABLE OF CONTENTS



SPECIALTY

Exos Bracing Systems	02
Cryo Compression Therapy	11
Dynamic Compression Venaflow	17
Anti-Embolism Stockings	23

UPPER EXTREMITY

Cervical Collars	26
Clavicle Supports	30
Shoulder Immobilizers	31
Arm Sling	33
Elbow Braces and Supports	34
Wrist and Hand Supports	37
Wrist and Thumb Supports	41
Lumbro-Sacral Supports	48
Abdominal & Rib Supports	51

LOWER EXTREMITY

Knee Immobilizers	54
Post-Op Knee Braces	56
Thigh and Calf Supports	58
Patellofemoral Supports	60
Hinged Knee Braces	63
Ankle Braces and Supports	71
Plantar Fasciitis	74
Post-Op Shoes	75
Walking Braces	80

PATIENT CARE

Pressure Care and Bandages	93
Fracture Bracing and Splints	102

EXOS BRACING SYSTEMS

**Short Arm Fracture Brace Open Thumb**

The open thumb design allows for enhanced ease of fit and patient comfort while the lycra top cover allows for faster heating and cooling to speed the application process. The ability to heat form and custom fit ensures proper uniform support and comfortable fit. Provides support and protection for major injuries to the wrist. Used primarily for fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
312-31-XXXX	312-32-XXXX	XS
312-41-XXXX	312-42-XXXX	S
312-51-XXXX	312-52-XXXX	M
312-61-XXXX	312-62-XXXX	L
312-71-XXXX	312-72-XXXX	XL

See color options and sizing chart on page 8.

**Pediatric Short Arm Fracture Brace Open Thumb**

The smallest of patients can benefit from Exos technology. The open thumb design allows for enhanced ease of fit and patient comfort while the lycra top cover allows for faster heating and cooling to speed the application process. The ability to heat form and custom fit ensures proper uniform support and comfortable fit. Exos braces are radiolucent and do not need to be removed for x-ray imaging. Provides support and protection for major injuries to the wrist. Used primarily for fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Does not contain latex.

Left	Right	Size
311-11-XXXX	311-12-XXXX	XXXS
311-21-XXXX	311-22-XXXX	XXS

See color options and sizing chart on page 8.

**Wrist Brace with Boa**

Intended for carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapho-lunate dislocations, triangular fibro-cartilage complex (TFCC) tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures. May be used to support stable wrist fractures. May be used to control wrist motion for other injuries that require stabilization. Does not contain latex.

Left	Right	Size
221-21-XXXX	221-22-XXXX	XXS
221-31-XXXX	221-32-XXXX	XS
221-41-XXXX	221-42-XXXX	S
221-51-XXXX	221-52-XXXX	M
221-61-XXXX	221-62-XXXX	L
221-71-XXXX	221-72-XXXX	XL

See color options and sizing chart on page 8.

**Wrist Brace**

Intended for carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapho-lunate dislocations, triangular fibro-cartilage complex (TFCC) tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures. May be used to support stable wrist fractures. May be used to control wrist motion for other injuries that require stabilization. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
220-31-1111	220-32-1111	XS
220-41-1111	220-42-1111	S
220-51-1111	220-52-1111	M
220-61-1111	220-62-1111	L
220-71-1111	220-72-1111	XL

Available in black only. See sizing chart on page 8.



Thumb Spica Fracture Brace

Designed to provide stabilization to the thumb, first metacarpophalangeal (MCP) joint, the first carpometacarpal (CMC) joint, wrist joint, distal radius and distal ulna. Useful for all injuries to the thumb, such as Bennett's or complex Rolando's fractures, suspected or confirmed navicular (scaphoid) fractures, and for other carpal bone injuries. Also useful for fractures to the distal radius and/or distal ulna, Colles fractures, Salter-Harris pediatric wrist and forearm injuries. Used for acute injuries as well as for post-op support. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
315-31-1111	315-32-1111	XS
315-41-1111	315-42-1111	S
315-51-1111	315-52-1111	M
315-61-1111	315-62-1111	L
315-71-1111	315-72-1111	XL

Available in black only. See sizing chart on page 8.



Long Thumb Spica

Provides primary stabilization for the first metacarpophalangeal (MCP) joint and first carpometacarpal (CMC) joints, as well as secondary support for the radiocarpal and radioulnar wrist joints. Provides stability for injuries to the base of the thumb at the CMC joint, such as Bennett's or Rolando's fractures, or to stabilize a painful CMC joint. Also useful for confirmed or suspected navicular (scaphoid) fractures, and for injuries or arthritis involving the trapezium or trapezoid bones of the wrist. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
230-31-1111	230-32-1111	XS
230-41-1111	230-42-1111	S
230-51-1111	230-52-1111	M
230-61-1111	230-62-1111	L
230-71-1111	230-72-1111	XL

Available in black only. See sizing chart on page 8.



Long Thumb Spica with Boa

Provides primary stabilization for the first metacarpophalangeal (MCP) joint and first carpometacarpal (CMC) joints, as well as secondary support for the radiocarpal and radioulnar wrist joints. Provides stability for injuries to the base of the thumb at the CMC joint, such as Bennett's or Rolando's fractures, or to stabilize a painful CMC joint. Also useful for confirmed or suspected navicular (scaphoid) fractures, and for injuries or arthritis involving the trapezium or trapezoid bones of the wrist. Does not contain latex.

Left	Right	Size
231-21-XXXX	231-22-XXXX	XXS
231-31-XXXX	231-32-XXXX	XS
231-41-XXXX	231-42-XXXX	S
231-51-XXXX	231-52-XXXX	M
231-61-XXXX	231-62-XXXX	L
231-71-XXXX	231-72-XXXX	XL

See color options and sizing chart on page 8.



Extended Short Thumb Spica

Immobilizes the base of the thumb and the interphalangeal joint, to stabilize injuries to the ulnar collateral ligament (UCL) commonly known as Gamekeeper's thumb, fractures or avulsion injuries to the proximal first phalanx or distal first metacarpal, severe sprains of the thumb, or for arthritic conditions involving the metacarpophalangeal (MCP) joint, or the interphalangeal (I-P joint). May also be used to immobilize the thumb joints following surgical repair of soft tissue injuries. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
211-31-1111	211-32-1111	XS
211-41-1111	211-42-1111	S
211-51-1111	211-52-1111	M
211-61-1111	211-62-1111	L
211-71-1111	211-72-1111	XL

Available in black only. See sizing chart on page 8.



Boxer's Fracture Brace

Used to stabilize fractures to the 4th and/or 5th metacarpal (MC) bones, and also provides stabilization to the 4th and 5th interphalangeal and wrist joints. May be used for metacarpal (MC) head/neck fracture of the 4th and/or 5th metacarpal. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
325-31-1111	325-32-1111	XS
325-41-1111	325-42-1111	S
325-51-1111	325-52-1111	M
325-61-1111	325-62-1111	L
325-71-1111	325-72-1111	XL

Available in black only. See sizing chart on page 8.




Radial Gutter Fracture Brace

Used to stabilize fractures to the 2nd and/or 3rd metacarpal (MC) bones and the respective metacarpal-phalangeal (MCP) joints of the hand. This pre-fabricated, custom-molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The radial gutter portion of the Radial Gutter Fracture Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 2nd and 3rd digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
326-31-1111	326-32-1111	XS
326-41-1111	326-42-1111	S
326-51-1111	326-52-1111	M
326-61-1111	326-62-1111	L
326-71-1111	326-72-1111	XL

Available in black only. See sizing chart on page 8.



Hand Based Ulnar Gutter Brace

The Hand Based Ulnar Gutter Brace is used to stabilize fractures to the 4th and/or 5th metacarpals or phalanges. This molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The ulnar gutter portion of the Hand Based Ulnar Gutter Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 4th and 5th digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
125-31-1111	125-32-1111	XS
125-41-1111	125-42-1111	S
125-51-1111	125-52-1111	M
125-61-1111	125-62-1111	L
125-71-1111	125-72-1111	XL

Available in black only. See sizing chart on page 8.



Hand Based Radial Gutter Brace

The Hand Based Radial Gutter Brace is used to stabilize fractures to the 2nd and/or 3rd metacarpals or phalanges. This molded orthosis is designed to provide superior fit and support compared to plaster, fiberglass, or pre-made braces. The radial gutter portion of the Hand Based Radial Gutter Brace allows the caregiver to achieve greater than 70° of flexion at the MCP joints of the 2nd and 3rd digits. Does not contain latex.

Left	Right	Size
N/A	N/A	XXS
126-31-1111	126-32-1111	XS
126-41-1111	126-42-1111	S
126-51-1111	126-52-1111	M
126-61-1111	126-62-1111	L
126-71-1111	126-72-1111	XL

Available in black only. See sizing chart on page 8.



Humeral Fracture Brace

Used to stabilize fractures to the mid-shaft of the humerus. The Humeral Fracture Brace may be appropriate for some injuries to the more proximal humerus. Does not contain latex.

Left	Size
330-40-1111	S/Univ.
330-50-1111	M/Univ.
330-60-1111	L/Univ.

Available in black only. See sizing chart on page 8.

Use the H measurement to select the appropriate size using the size chart below.
H=Height

SMALL 4¼"
MEDIUM 5"
LARGE 6"



Exos Top Loading Oven

Convenient sized oven used to heat the Exos Brace in order to gently mold the brace to the patient's limb.

	Description
800-OVEN3P	Exos Top Loading Oven
Accessories	
800-ER	Oven Extender Ring
800-OVENBAG	Oven Bag

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
125-31-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, LT, XS	\$85.63
125-32-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, RT, XS	\$85.63
125-41-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, LT, SM	\$85.63
125-42-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, RT, SM	\$85.63
125-51-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, LT, MD	\$85.63
125-52-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, RT, MD	\$85.63
125-61-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, LT, LG	\$85.63
125-62-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, RT, LG	\$85.63
125-71-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, LT, XL	\$85.63
125-72-1111	Enovis Exos ULNAR GUTTER BRACE, HAND BASED, BLK, RT, XL	\$85.63
126-31-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, LT, XS	\$85.63
126-32-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, RT, XS	\$85.63
126-41-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, LT, SM	\$85.63
126-42-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, RT, SM	\$85.63
126-51-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, LT, MD	\$85.63
126-52-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, RT, MD	\$85.63
126-61-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, LT, LG	\$85.63
126-62-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, RT, LG	\$85.63
126-71-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, LT, XL	\$85.63
126-72-1111	Enovis Exos RADIAL GUTTER BRACE, HAND BASED, BLK, RT, XL	\$85.63
211-31-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, LT, XS	\$90.51
211-32-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, RT, XS	\$90.51
211-41-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, LT, S	\$90.51
211-42-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, RT, S	\$90.51
211-51-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, LT, M	\$90.51
211-52-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, RT, M	\$90.51
211-61-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, LT, L	\$90.51
211-62-1111	Enovis Exos EXTENDED SHRT THMB SPICA UBL, BLK, RT, L	\$90.51
211-71-1111	Enovis Exos EXTENDED SHRT THMB SPICA, BLK, LT, XL	\$90.51
211-72-1111	Enovis Exos EXTENDED SHRT THMB SPICA, BLK, RT, XL	\$90.51
220-31-1111	Enovis Exos WRIST BRACE, BLK, LT, XS	\$118.91
220-32-1111	Enovis Exos WRIST BRACE, BLK, RT, XS	\$118.91
220-41-1111	Enovis Exos WRIST BRACE, BLK, LT, S	\$118.91
220-42-1111	Enovis Exos WRIST BRACE, BLK, RT, S	\$118.91
220-51-1111	Enovis Exos WRIST BRACE, BLK, LT, M	\$118.91
220-52-1111	Enovis Exos WRIST BRACE, BLK, RT, M	\$118.91
220-61-1111	Enovis Exos WRIST BRACE, BLK, LT, L	\$118.91
220-62-1111	Enovis Exos WRIST BRACE, BLK, RT, L	\$118.91

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
220-71-1111	Enovis Exos WRIST BRACE, BLK, LT, XL	\$118.91
220-72-1111	Enovis Exos WRIST BRACE, BLK, RT, XL	\$118.91
221-21-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, XXS	\$123.67
221-21-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, LT, XXS	\$123.67
221-21-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, LT, XXS	\$123.67
221-21-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, LT, XXS	\$123.67
221-21-4444	Enovis Exos WRIST BRACE W/BOA, PINK, LT, XXS	\$123.67
221-22-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, XXS	\$123.67
221-22-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, RT, XXS	\$123.67
221-22-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, RT, XXS	\$123.67
221-22-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, RT, XXS	\$123.67
221-22-4444	Enovis Exos WRIST BRACE W/BOA, PINK, RT, XXS	\$123.67
221-31-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, XS	\$123.67
221-31-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, LT, XS	\$123.67
221-31-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, LT, XS	\$123.67
221-31-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, LT, XS	\$123.67
221-31-4444	Enovis Exos WRIST BRACE W/BOA, PINK, LT, XS	\$123.67
221-32-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, XS	\$123.67
221-32-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, RT, XS	\$123.67
221-32-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, RT, XS	\$123.67
221-32-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, RT, XS	\$123.67
221-32-4444	Enovis Exos WRIST BRACE W/BOA, PINK, RT, XS	\$123.67
221-41-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, S	\$123.67
221-41-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, LT, S	\$123.67
221-41-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, LT, S	\$123.67
221-41-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, LT, S	\$123.67
221-41-4444	Enovis Exos WRIST BRACE W/BOA, PINK, LT, S	\$123.67
221-42-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, S	\$123.67
221-42-2293	Enovis Exos WRIST BRACE W/BOA, SPIDER, RT, S	\$123.67
221-42-2388	Enovis Exos WRIST BRACE W/BOA, CAMO, RT, S	\$123.67
221-42-3285	Enovis Exos WRIST BRACE W/BOA, POLKA-DOT, RT, S	\$123.67
221-42-4444	Enovis Exos WRIST BRACE W/BOA, PINK, RT, S	\$123.67
221-51-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, M	\$123.67
221-52-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, M	\$123.67
221-61-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, L	\$123.67
221-62-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, L	\$123.67
221-71-1111	Enovis Exos WRIST BRACE W/BOA, BLK, LT, XL	\$123.67
221-72-1111	Enovis Exos WRIST BRACE W/BOA, BLK, RT, XL	\$123.67

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
230-31-1111	Enovis Exos LNG THMB SPICA, BLK, LT, XS	\$123.67
230-32-1111	Enovis Exos LNG THMB SPICA, BLK, RT, XS	\$123.67
230-41-1111	Enovis Exos LNG THMB SPICA, BLK, LT, S	\$123.67
230-42-1111	Enovis Exos LNG THMB SPICA, BLK, RT, S	\$123.67
230-51-1111	Enovis Exos LNG THMB SPICA, BLK, LT, M	\$123.67
230-52-1111	Enovis Exos LNG THMB SPICA, BLK, RT, M	\$123.67
230-61-1111	Enovis Exos LNG THMB SPICA, BLK, LT, L	\$123.67
230-62-1111	Enovis Exos LNG THMB SPICA, BLK, RT, L	\$123.67
230-71-1111	Enovis Exos LNG THMB SPICA , BLK, LT, XL	\$123.67
230-72-1111	Enovis Exos LNG THMB SPICA , BLK, RT, XL	\$123.67
231-21-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, LT, XXS	\$123.93
231-22-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, XXS	\$123.93
231-31-1111	Enovis Exos LONG THUMB SPICA W/BOA, BLK, LT, XS	\$123.93
231-32-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, XS	\$123.93
231-41-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, LT, S	\$123.93
231-42-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, S	\$123.93
231-51-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, LT, M	\$123.93
231-52-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, M	\$123.93
231-61-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, LT, L	\$123.93
231-62-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, L	\$123.93
231-71-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, LT, XL	\$123.93
231-72-1111	Enovis Exos LNG THMB SPICA W/BOA, BLK, RT, XL	\$123.93
311-11-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, LT, 3XS/PEDS	\$148.48
311-11-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, LT, 3XS/PEDS	\$148.48
311-11-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, LT, 3XS/PEDS	\$148.48
311-11-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, LT, 3XS/PEDS	\$148.48
311-11-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, LT, 3XS/PEDS	\$148.48
311-12-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, RT, 3XS/PEDS	\$148.48
311-12-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, RT, 3XS/PEDS	\$148.48
311-12-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, RT, 3XS/PEDS	\$148.48
311-12-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, RT, 3XS/PEDS	\$148.48
311-12-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, RT, 3XS/PEDS	\$148.48
311-21-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, LT, 2XS/PEDS	\$148.48
311-21-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, LT, 2XS/PEDS	\$148.48
311-21-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, LT, 2XS/PEDS	\$148.48
311-21-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, LT, 2XS/PEDS	\$148.48
311-21-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, LT, 2XS/PEDS	\$148.48
311-22-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, RT, 2XS/PEDS	\$148.48

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
311-22-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, RT, 2XS/PEDS	\$148.48
311-22-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, RT, 2XS/PEDS	\$148.48
311-22-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, RT, 2XS/PEDS	\$148.48
311-22-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, RT, 2XS/PEDS	\$148.48
311-31-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, LT, 1XS/PEDS	\$148.48
311-31-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, LT, 1XS/PEDS	\$148.48
311-31-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, LT, 1XS/PEDS	\$148.48
311-31-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, LT, 1XS/PEDS	\$148.48
311-31-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, LT, 1XS/PEDS	\$148.48
311-32-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, RT, 1XS/PEDS	\$148.48
311-32-2293	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, SPIDER, RT, 1XS/PEDS	\$148.48
311-32-2388	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, CAMO, RT, 1XS/PEDS	\$148.48
311-32-3285	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, POLKA-DOT, RT, 1XS/PEDS	\$148.48
311-32-4444	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, PINK, RT, 1XS/PEDS	\$148.48
312-31-1111	Enovis Exos SHORT ARM FRACTURE BRACE, OPEN THUMB, XS, LT, BLK	\$144.14
312-31-2293	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL LEFT SPIDER	\$148.48
312-31-2388	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL LEFT CAMOU	\$148.48
312-31-3285	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL LEFT POLKA	\$148.48
312-31-4444	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL LEFT PINK	\$148.48
312-32-1111	Enovis Exos SHORT ARM FRACTURE BRACE, OPEN THUMB, XS,RT, BLK	\$144.14
312-32-2293	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL RIGHT SPIDER	\$148.48
312-32-2388	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL RIGHT CAMOU	\$148.48
312-32-3285	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL RIGHT POLKA	\$148.48
312-32-4444	Enovis Exos SHORT ARM FRACTURE BRACE-XSMALL RIGHT PINK	\$148.48
312-41-1111	Enovis Exos SHRT ARM FRACTURE BRACE, OPEN THUMB, BLK, LT, SM	\$144.14
312-41-2293	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL LEFT SPIDER	\$148.48
312-41-2388	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL LEFT CAMOU	\$148.48
312-41-3285	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL LEFT POLKA	\$148.48
312-41-4444	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL LEFT PINK	\$148.48
312-42-1111	Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB, SM,RT, BLK	\$144.14
312-42-2293	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL RIGHT SPIDER	\$148.48
312-42-2388	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL RIGHT CAMOU	\$148.48
312-42-3285	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL RIGHT POLKA	\$148.48
312-42-4444	Enovis Exos SHORT ARM FRACTURE BRACE-SMALL RIGHT PINK	\$148.48
312-51-1111	Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB, M,LT, BLK	\$144.14
312-52-1111	Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB, M,RT, BLK	\$144.14
312-61-1111	Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB,L,LT, BLK	\$144.14
312-62-1111	Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB,L,RT, BLK	\$144.14

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
312-71-1111	312-71-1111 Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB,XL,LT, BLK	\$144.14
312-72-1111	312-72-1111 Enovis Exos SHORT ARM FRACTURE BRACE , OPEN THUMB,XL,RT, BLK	\$144.14
315-31-1111	315-31-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, LT, XS	\$145.81
315-32-1111	315-32-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, RT, XS	\$145.81
315-41-1111	315-41-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, LT, S	\$145.81
315-42-1111	315-42-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, RT, S	\$145.81
315-51-1111	315-51-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, LT, M	\$145.81
315-52-1111	315-52-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, RT, M	\$145.81
315-61-1111	315-61-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, LT, L	\$145.81
315-62-1111	315-62-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, RT, L	\$145.81
315-71-1111	315-71-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, LT, XL	\$145.81
315-72-1111	315-72-1111 Enovis Exos THMB SPICA FRACTURE BRACE, BLK, RT, XL	\$145.81
325-31-1111	325-31-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, LT, XS	\$144.14
325-32-1111	325-32-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, RT, XS	\$144.14
325-41-1111	325-41-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, LT, S	\$144.14
325-42-1111	325-42-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, RT, S	\$144.14
325-51-1111	325-51-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, LT, M	\$144.14
325-52-1111	325-52-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, RT, M	\$144.14
325-61-1111	325-61-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, LT, L	\$144.14
325-62-1111	325-62-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, RT, L	\$144.14
325-71-1111	325-71-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, LT, XL	\$144.14
325-72-1111	325-72-1111 Enovis Exos BOXER'S FRACTURE BRACE, BLK, RT, XL	\$144.14
326-31-1111	326-31-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, LT, XS	\$148.48
326-32-1111	326-32-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, RT, XS	\$148.48
326-41-1111	326-41-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, LT, SM	\$148.48
326-42-1111	326-42-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, RT, SM	\$148.48
326-51-1111	326-51-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, LT, MD	\$148.48
326-52-1111	326-52-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, RT, MD	\$148.48
326-61-1111	326-61-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, LT, LG	\$148.48
326-62-1111	326-62-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, RT, LG	\$148.48
326-71-1111	326-71-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, LT, XL	\$148.48
326-72-1111	326-72-1111 Enovis Exos RADIAL GUTTER FRACTURE BRACE, BLK, RT, XL	\$148.48
330-40-1111	330-40-1111 Enovis Exos HUMERAL FRACTURE BRACE, UNIV, S, BLK	\$158.72
330-50-1111	330-50-1111 Enovis Exos HUMERAL FRACTURE BRACE, UNIV, M, BLK	\$158.72
330-60-1111	330-60-1111 Enovis Exos HUMERAL FRACTURE BRACE, UNIV, L, BLK	\$158.72

CRYO-COMPRESSION THERAPY



Cryo/Cuff™ IC Cooler

Integrated cold and compression is clinically proven to help reduce post-operative swelling and regain range of motion. The Cryo/Cuff IC provides both gravity and motorized cold and compression therapy for flexibility of care. An integrated pneumatic pump within the cooler lid provides automated cold and compression therapy, with simplicity of design and ease-of-use operation.

Indication: Ideal for post-operative recovery, trauma, athletic training rooms, and home use.

	Description	Size
51A	IC Cooler Only	Univ.
51A10A	IC Cooler with Ankle Cryo/Cuff	Univ.
51A11A	IC Cooler with Knee Cryo/Cuff	M
51A11B	Cryo/Cuff IC Cooler with Knee Cryo/Cuff	L
51A12A	Cryo/Cuff IC Cooler with Shoulder	Univ.
51A12AXL	Cryo/Cuff IC Cooler with Shoulder, XLarge Strap	Univ.
51A14A	Cryo/Cuff IC Cooler with Back/Hip/Rib	Univ.
Accessories		
25-4882	Cryo/Cuff IC Cooler Power Supply	Univ.
25-0238	Cryo/Cuff IC Lid	Univ.
2130	Tube Assembly	Univ.
2002	Insulation Disk	Univ.





AIRCAST®

Cryo/Cuff™ Gravity Cooler

An integrated approach to cold therapy that combines cold and focal compression delivered using a simple, gravity-fed system. Designed with safety and ease-of-use in mind.

Indication: Ideal for post-operative recovery, trauma, athletic training rooms, and home-use.

	Description	Size
2125	Gravity Cooler Only	Univ.
10A	Gravity Cooler w/Ankle Cryo/Cuff	Univ.
10C	Gravity Cooler w/Foot Cryo/Cuff	S/M
10B	Gravity Cooler w/Foot Cryo/Cuff	L
10P	Gravity Cooler w/Pediatric Ankle Cryo/Cuff	Ped.
11C	Gravity Cooler w/Knee Cryo/Cuff	S
11A	Gravity Cooler w/Knee Cryo/Cuff	M
11B	Gravity Cooler w/Knee Cryo/Cuff	L
12A	Gravity Cooler w/Shoulder Cryo/Cuff	Univ.
12AXL	Gravity Cooler w/Shoulder XL Strap Cryo/Cuff	Univ.
14A	Gravity Cooler w/Back/Hip/Rib Cryo/Cuff	Univ.
15A	Gravity Cooler w/Elbow Cryo/Cuff	Univ.
16A	Gravity Cooler w/Hand/Wrist Cryo/Cuff	Univ.



AIRCAST®

Cooler Hanger

For use with any Aircast Cryo/Cuff coolers. Enables cooler to be hung from a bed rail.

	Description
2140	Cooler Hanger



AIRCAST®

Ankle Cryo/Cuff

Anatomically designed to completely fit the ankle. Minimizes hemarthrosis, swelling and pain. Universal adult size. Also available in pediatric size.

Indication: Ideal use for acute sprains, trauma, post-op, and rehabilitation.

	Description	Length x Height	Size
10A01	Ankle Cryo/Cuff Only	11.5" x 10"	Univ.
10P01	Ankle Cryo/Cuff Only	8" x 7"	Ped.



AIRCAST®

Foot Cryo/Cuff

Anatomically designed to provide complete foot coverage for optimal management of swelling and pain.

Indication: Ideal use for turf toe and post-op.

	Description	Circumference*	Size
10C01	Foot Cryo/Cuff Only	9" – 13"	M
10B01	Foot Cryo/Cuff Only	10" – 17"	L

*Measure forefoot circumference.



AIRCAST®

Knee Cryo/Cuff

Anatomically designed to completely fit the knee. Minimizes hemarthrosis, swelling, and pain. Available in three adult sizes. Also available in pediatric size.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Circumference*	Size
11C01	Knee Cryo/Cuff Only	10" – 19"	S
11A01	Knee Cryo/Cuff Only	18" – 23"	M
11B01	Knee Cryo/Cuff Only	20" – 31"	L
11P01	Knee/Elbow Cryo/Cuff	Pediatric	P

*Measure 6" above the patella.



AIRCAST®

Self-Contained (SC) Knee Cryo/Cuff

A totally self-contained knee cuff that can be filled with ice and water - no cooler required. The SC Knee Cryo/Cuff remains cold for over one hour. A hand bulb is included to add compression as needed.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Size
11ASC	Self-Contained Knee Cryo/Cuff	Univ.



AIRCAST®

Shoulder Cryo/Cuff

Anatomically designed to cover the shoulder area completely. Minimizes hemarthrosis, swelling, and pain. Can be ordered with an XL strap to accommodate larger torso/chest circumferences.

Indication: Ideal use for trauma, post-op, rehabilitation, and sports injuries.

	Description	Circumference*	Size
12B01	Shoulder Cryo/Cuff Only	Up to 34" and below	S
12A01	Shoulder Cryo/Cuff Only	34"+	Univ.
12AXL01	Shoulder Cryo/Cuff with X-Long Strap Only	42" – 54"	Univ.

Accessories

1200XL	Shoulder Cryo/Cuff XL Strap	Univ.
--------	-----------------------------	-------

*Measure torso circumference. (Chest circumference)



AIRCAST®

Calf Cryo/Cuff

Anatomically designed to completely fit the calf. Minimizes swelling and pain. 10" long.

Indication: Ideal use for acute injury, sports injuries, and soft tissue injury.

	Description	Circumference*	Size
13C01	Calf Cryo/Cuff Only	14" – 20"	Univ.

*Measure mid-calf circumference.



AIRCAST®

Thigh Cryo/Cuff

Designed to completely fit the thigh. Minimizes swelling and pain.

Indication: Ideal use for stress fractures, sports injuries, and soft tissue injury.

Description	Circumference*	Size
13A01 Thigh Cryo/Cuff Only	17" – 23"	S/M
13B01 Thigh Cryo/Cuff Only	20" – 27"	L/XL

*Measure thigh circumference.



AIRCAST®

Back/Hip/Rib Cryo/Cuff

Universal fit. The versatile design may be applied either to the back, hip, or rib area. Minimizes swelling and pain.

Indication: Ideal use for chronic/acute pain; post-op and sports injuries.

Description	Size
14A01 Back/Hip/Rib Cryo/Cuff Only	Univ. up to 45" waist



AIRCAST®

Elbow Cryo/Cuff

Anatomically designed to fit the elbow. Minimizes hemarthrosis, swelling, and pain. One size for adults. Also available in a pediatric size.

Indication: Ideal use for epicondylitis (tennis elbow), post-op, and throwing injuries.

Description	Size
15A01 Elbow Cryo/Cuff Only	Univ.
11P01 Elbow/Knee Cryo/Cuff Only	Ped.



AIRCAST®

Hand/Wrist Cryo/Cuff

Anatomically designed to cover the hand and wrist. Versatile design allows repositioning for desired hand coverage. Minimizes swelling and pain. For added support, aluminum support bars are incorporated and can be removed if increased range of motion is desired. 12" long.

Indication: Ideal use for carpal tunnel syndrome, tendonitis, post-op, and sports injuries.

Description	Circumference*	Size
16A01 Hand/Wrist Cryo/Cuff Only	12"	Univ.

*Measure forearm circumference.

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
10A	Enovis AIRCAST ANKLE CRYO/CUFF WITH COOLER	\$151.63
10A01	Enovis AIRCAST ANKLE CUFF ONLY	\$72.55
10B	Enovis AIRCAST LARGE FOOT CRYO/CUFF W/COOLER	\$151.63
10B01	Enovis AIRCAST LARGE FOOT CUFF ONLY	\$72.55
10C	Enovis AIRCAST MED FOOT CRYO/CUFF WITH COOLER	\$151.63
10C01	Enovis AIRCAST MED FOOT CUFF ONLY	\$72.55
10P01	Enovis AIRCAST PEDIATRIC ANKLE CUFF ONLY	\$101.54
11A	Enovis AIRCAST MEDIUM KNEE CRYO/CUFF W/COOLER	\$151.63
11A01	Enovis AIRCAST MEDIUM KNEE CUFF ONLY	\$70.45
11AL314959	Enovis AIRCAST UNIV KY MED KNEE C/C W/COOLER	\$135.56
11B	Enovis AIRCAST LARGE KNEE CRYO/CUFF W/COOLER	\$151.63
11B01	Enovis AIRCAST LARGE KNEE CUFF ONLY	\$70.45
11C	Enovis AIRCAST SMALL KNEE CRYO/CUFF W/COOLER	\$148.50
11C01	Enovis AIRCAST SMALL KNEE CUFF ONLY	\$72.55
11P01	Enovis AIRCAST PEDIATRIC KNEE/ELBOW CUFF ONLY	\$101.54
1220	Enovis AIRCAST STRAP,TORSO SHOULDER CRYO/CUFF	\$10.95
1220XL	Enovis AIRCAST STRAP,TORSO SHLDR CRYO/CUFF,XL	\$10.19
12A	Enovis AIRCAST SHOULDER CRYO/CUFF W/COOLER	\$166.62
12A01	Enovis AIRCAST SHOULDER CRYO/CUFF PACKAGING	\$112.59
12AXL	Enovis AIRCAST SHOULDER CUFF W/COOLER,XLSTRAP	\$166.62
12AXL01	Enovis AIRCAST SHOULDER CUFF ONLY W/XL STRAP	\$115.17
12B01	Enovis AIRCAST SMALL SHOULDER CRYO/CUFF	\$112.59
13A01	Enovis AIRCAST THIGH CUFF ONLY	\$72.55
13B01	Enovis AIRCAST THIGH CUFF ONLY, XL	\$72.55
13C01	Enovis AIRCAST CALF CUFF ONLY	\$72.55
14A	Enovis AIRCAST BACK/HIP/RIB CRYO/CUFF W/COOLR	\$151.63
14A01	Enovis AIRCAST BACK/HIP/RIB CUFF ONLY	\$72.55
15A	Enovis AIRCAST ELBOW CRYO/CUFF WITH COOLER	\$151.63
15A01	Enovis AIRCAST ELBOW CUFF ONLY	\$72.55
16A	Enovis AIRCAST HAND & WRIST CRYO/CUFF W/COOLR	\$151.63
16A01	Enovis AIRCAST HAND/WRIST CUFF ONLY	\$72.55
2002	Enovis AIRCAST INSULATION DISK AUTOCHILL	\$3.95
2125	Enovis AIRCAST COOLER AND TUBE ASSEMBLY	\$58.34
2130	Enovis AIRCAST CRYO/CUFF TUBE ASSEMBLY	\$13.78
2140	Enovis AIRCAST COOLER HANGER	\$40.73
51A	Enovis AIRCAST CRYO/CUFF IC COOLER	\$130.18
51A10A	Enovis AIRCAST CRYO/CUFF IC W/ANKLE	\$188.21
51A11A	Enovis AIRCAST CRYO/CUFF IC W/KNEE, MED	\$188.21

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
51A11B	Enovis AIRCAST CRYO/CUFF IC W/KNEE, LG	\$188.21
51A12A	Enovis AIRCAST CRYO/CUFF IC W/SHOULDER	\$201.11
51A14A	Enovis AIRCAST CRYO/CUFF IC W/HIP,RIB, BACK	\$188.21

DYNAMIC COMPRESSION - VENAFLOW

Vascular Systems

While intermittent pneumatic compression is clinically proven to reduce the risks of deep vein thrombosis, The VenaFlow system is the only device on the market proven to mimic ambulation. Through its unique technology, it has been shown to reduce DVT in comparison to other pneumatic compression systems.^{1, 2, 3, 4}

VENAFLOW®

VenaFlow Elite System

Designed as a prophylaxis for Deep Vein Thrombosis, the VenaFlow Elite System combines two proven technologies, rapid inflation and graduated sequential compression that work to mimic ambulation and accelerate venous velocity. For the superior emptying of veins, VenaFlow Elite also incorporates asymmetric compression. The breathable and comfortable calf, foot and thigh cuffs enhance the patient experience and are compatible with one universal pump.

Indication: For the prevention of thrombus formation.

	Description
30B	Elite System
30B-B	Elite System with Battery (additional charge)

References

¹Eisele R, Kinzl L, Koelsch T: "Rapid-Inflation Pneumatic Compression For Prevention of Deep Venous Thrombosis." JBJS, 2008.

²Dai G, Gertler GP, Kamm RD: "The Effects of External Compression on Venous Blood Flow and Tissue Deformation in the Lower Leg." Journal of Biomechanical Engineering, 1999.

³Gardner AMN, Fox RH: "The Return of Blood to the Heart, Venous Pumps in Health and Disease." Second Edition, London, John Libbey & Company Ltd, 1993.

⁴Lachiewicz PF, Kelley SS, Haden LR: "Two Mechanical Devices for Prophylaxis of Thromboembolism After Total Knee Arthroplasty." JBJS, 2004.





VENAFLOW®

VenaFlow Elite Calf Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite calf cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. Available in sterile variation. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3040	Calf Cuff	19"	11.5"	2
3041	Sterile Calf Cuff	19"	11.5"	1
3042	Calf Cuff XL	22"	11.5"	2
3043	Bariatric Calf Cuff	30"	11.5"	2
3044	Sterile Calf Cuff XL	22"	11.5"	1
3050	Foam Calf Cuff	19"	11.5"	2



VENAFLOW®

VenaFlow Elite System Tubing (30B, 30B-B Models)

Designed for use with the VenaFlow Elite System. The tubing connects the VenaFlow Elite System to a VenaFlow cuff and is available in a variety of lengths.

Indication: To be used in conjunction with the VenaFlow Elite System for the prevention of thrombus formation.

	Description	Length
3008	Regular	5.5'
3008XL	X-Long	8.5'
3008XXL	XX-Long	10.5'
3008XXXL	XXX-Long	12.5'



VENAFLOW®

VenaFlow Elite Foot Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite foot cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3046	Foot Cuff	16"	9"	2



VENAFLOW®

VenaFlow Elite Thigh Cuff (30B, 30B-B Models)

Used exclusively with the VenaFlow Elite System, the VenaFlow Elite thigh cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circumference	Length	Qty.
3045	Thigh Cuff	29"	21"	2



VENAFLOW®

VenaFlow Elite System Repair Kits (30B, 30B-B Models)

Repair kits for use with the VenaFlow Elite tubing system.

Indication: To be used in conjunction with the VenaFlow Elite System for the prevention of thrombus formation.

	Description
3070	Power cord (U.S.)
3071	Bed hanger
3072	Tube attachment tag
3073	Fuse



VENAFLOW®

VenaFlow® System

Designed as a prophylaxis for Deep Vein Thrombosis, the VenaFlow System combines two proven technologies, rapid inflation and graduated sequential compression, to accelerate venous velocity and enhance fibrinolysis. For the superior emptying of veins, VenaFlow also incorporates asymmetric compression. The breathable and comfortable calf, foot and thigh cuffs enhance the patient experience and are compatible with one universal pump.

Indication: For the prevention of thrombus formation.

	Description
30A	VenaFlow System



VENAFLOW®

VenaFlow Foot Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3016-PL	Foot Cuff	16"	9"	2



VENAFLOW®

VenaFlow Calf Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. Available in sterile variation. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3010-PL	Calf Cuff	19"	11.5"	2
3011-PL	Sterile Calf Cuff	19"	11.5"	1
3012-PL	Calf Cuff XL	22"	11.5"	2
3013-PL	Bariatric Calf Cuff	30"	11.5"	2
3014-PL	Sterile Calf Cuff XL	22"	11.5"	1
3030-PL	Foam Calf Cuff	19"	11.5"	2



VENAFLW®

VenaFlow Thigh Cuff (30A Model)

Used exclusively with the VenaFlow System, the VenaFlow cuffs are latex-free and may be placed directly against the skin. The cuff material is light, cool and comfortable to promote patient compliance. With slight adjustment, the patient may ambulate while wearing cuffs. For single patient use only.

Indication: For the prevention of thrombus formation.

	Description	Circum.	Length	Qty.
3015-PL	Thigh Cuff	29"	21"	2



VENAFLW®

VenaFlow System Repair Kits (30A Model)

Repair kits for use with the VenaFlow system.

Indication: To be used in conjunction with the VenaFlow System for the prevention of thrombus formation.

	Description
11-0201	Tube Storage Strap Repair Kit
11-0202	Power Cord Repair Kit
11-0203	Pump/Tubing Connector Repair Kit
11-0204	Pump Case (Back Half) Repair Kit
11-0205	Bed Hanger Foam Repair Kit
11-0206	Fuse Repair Kit



VENAFLW®

VenaFlow System Tubing (30A Model)

Designed for use with the VenaFlow System. The tubing connects the VenaFlow System to a VenaFlow cuff and is available in a variety of lengths.

Indication: To be used in conjunction with the VenaFlow System for the prevention of thrombus formation.

	Description	Length
3007-PL	Regular	5.5'
3007XL-PL	X-Long	8.5'
3007XXL-PL	XX-Long	10.5'

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
3004	Enovis AIRCAST PRESSURE VERIFICATION KIT	\$144.38
3007-PL	Enovis AIRCAST VF TUBE ASSY W/POS LOCK	\$55.37
3007-PL-V	Enovis AIRCAST 12.5FT-PL TUBE ASSEMBLY	\$64.81
3007XL	Enovis AIRCAST VENAFLOW TUBE ASSEMBLY, X LONG	\$66.44
3007XL-PL	Enovis AIRCAST VF XL TUBE ASSY W/POS LOCK	\$66.44
3007XXL	Enovis AIRCAST VENAFLOW TUBE ASSEMBLY,XX LONG	\$66.44
3007XXL-PL	Enovis AIRCAST VF XXL TUBE ASSY W/POS LOCK	\$66.44
3008	Enovis AIRCAST TUBE ASSY, VENAFLOW ELITE, STD	\$46.19
3008XL	Enovis AIRCAST TUBE ASSY, VENAFLOW ELITE, XL	\$55.41
3008XXL	Enovis AIRCAST TUBE ASSY, VENAFLOW ELITE, XXL	\$55.41
3008XXXL	Enovis AIRCAST TUBE ASSY, VENAFLOW ELITE,XXXL	\$55.41
3010	Enovis AIRCAST VENAFLOW DSP CALF CUFF ASSY PR	\$72.18
3010-PL	Enovis AIRCAST VF DSP CALF CUFF W/POS LOCK	\$72.18
3011	Enovis AIRCAST STERILE VENAFLOW DSP CUFF	\$101.32
3011-PL	Enovis AIRCAST VF-PL STERILE CUFF PKG	\$108.76
3012	Enovis AIRCAST VF DISP CALF CUFF ASSY (PR),XL	\$72.18
3012-PL	Enovis AIRCAST XL CALF CUFF W/POS LOCK	\$72.18
3013	Enovis AIRCAST BARIATRC VF DISP CLF CF ASY PR	\$102.89
3013-PL	Enovis AIRCAST BARIATRIC C/CUFF W/POS LOCK	\$70.69
3014	Enovis AIRCAST STERILE VF DISP CUFF, XL	\$101.32
3014-PL	Enovis AIRCAST VF XL-PL STERILE CUFF PKG	\$108.76
3015	Enovis AIRCAST VENAFLOW DISP THI CUFF ASSY PR	\$116.12
3015-PL	Enovis AIRCAST VF THIGH CUFF W/POS LOCK	\$116.12
3016	Enovis AIRCAST VENAFLOW DSP FOOT CUFF ASSY PR	\$127.15
3016-PL	Enovis AIRCAST VF FOOT CUFF W/POS LOCK	\$127.15
3030-PL	Enovis AIRCAST VENAFLOW FOAM CALF CUFF	\$72.18
3040	Enovis AIRCAST VENAFLOW ELITE CALF CUFF	\$39.37
3041	Enovis AIRCAST VF ELITE STERILE CALF CUFF	\$101.32
3042	Enovis AIRCAST VENAFLOW ELITE XL CALF CUFF	\$39.37
3043	Enovis AIRCAST VF ELITE BARIATRIC CALF CUFF	\$41.52
3044	Enovis AIRCAST VF ELITE STERILE XL CALF CUFF	\$101.32
3045	Enovis AIRCAST VENAFLOW ELITE THIGH CUFF	\$108.02
3046	Enovis AIRCAST VENAFLOW ELITE FOOT CUFF	\$66.33
3050	Enovis AIRCAST VENAFLOW ELITE FOAM CALF CUFF	\$39.37
3073	Enovis AIRCAST REPL FUSES,2 COUNT,VF ELITE	\$6.82
30A	Enovis AIRCAST VENAFLOW SYSTEM	\$1,941.00
30B	Enovis AIRCAST VENAFLOW ELITE SYSTEM	\$2,427.00
30B-B	Enovis AIRCAST VF ELITE SYSTEM W/BATTERY	\$2,547.00

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
30B-S	Enovis AIRCAST VENAFLOW ELITE SYSTEM, WITH S-MODE	\$2,427.00
30B-SB	Enovis AIRCAST VENAFLOW ELITE SYSTEM, WITH S-MODE & BATTERY	\$2,547.00
30M	Enovis AIRCAST VENAPRO SYSTEM	\$359.79

STATIC COMPRESSION- ANTI EMBOLISM STOCKINGS

VENAFLow[®]AES

VenaFlow Knee High Anti Embolism Stockings

Designed with quality, fit and safety in mind, the VenaFlow Anti Embolism stockings are key in preventing VTE. They enhance venous flow velocity while keeping patients comfortable during and after their hospital stay.

Indication: For the prevention of thrombus formation.

	Description	Qty.	Size
Regular [®]			
500R-2	Knee High Stocking	12 pair/pk.	S
500R-3	Knee High Stocking	12 pair/pk.	M
500R-4	Knee High Stocking	12 pair/pk.	L
500R-5	Knee High Stocking	12 pair/pk.	XL
500R-6	Knee High Stocking	12 pair/pk.	XXL
500R-7	Knee High Stocking	12 pair/pk.	XXXL
Long [®]			
500L-2	Knee High Stocking	12 pair/pk.	S
500L-3	Knee High Stocking	12 pair/pk.	M
500L-4	Knee High Stocking	12 pair/pk.	L
500L-5	Knee High Stocking	12 pair/pk.	XL
500L-6	Knee High Stocking	12 pair/pk.	XXL
500L-7	Knee High Stocking	12 pair/pk.	XXXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.





VENAFLOW[®] AES

VenaFlow Thigh High Anti Embolism Stockings

Designed with quality, fit and safety in mind, the VenaFlow Anti Embolism stockings are key in preventing VTE. They enhance venous flow velocity while keeping patients comfortable during and after their hospital stay.

Indication: For the prevention of thrombus formation.

Thigh High	Thigh High w/Belt	Qty.	Size
Regular[®]			
501R-2	502R-2	6 pair/pk.	S
501R-3	502R-3	6 pair/pk.	M
501R-4	502R-4	6 pair/pk.	L
501R-5	502R-5	6 pair/pk.	XL
501R-6	502R-6	6 pair/pk.	XXL
Long[®]			
501L-2	502L-2	6 pair/pk.	S
501L-3	502L-3	6 pair/pk.	M
501L-4	502L-4	6 pair/pk.	L
501L-5	502L-5	6 pair/pk.	XL
501L-6	502L-6	6 pair/pk.	XXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.

Thigh High	Thigh High w/Belt	Qty.	Size
Short[®]			
501S-2	502S-2	6 pair/pk.	S
501S-3	502S-3	6 pair/pk.	M
501S-4	502S-4	6 pair/pk.	L
501S-5	502S-5	6 pair/pk.	XL
501S-6	502S-6	6 pair/pk.	XXL

[®]All stockings provide compression in mmHg: Compression at ankle 18 mmHg, Compression at calf 14 mmHg, Compression at knee 8 mmHg, Compression at lower thigh 10 mmHg, Compression at upper thigh 8 mmHg.

VENAFLOW ANTI EMBOLISM KNEE LENGTH STOCKINGS SIZING CHART [®]					
	Size	1 - Calf Circumference	2 - Length		Color
			Regular	Long	
	S	< 30.5 cm / 12"	< 41 cm / < 16"	> 41 cm / < 16"	Yellow
	M	30.5 – 38 cm / 12" – 15"	< 43 cm / < 17"	> 43 cm / < 17"	White
	L	38 – 44.5 cm / 15" – 17.5"	< 46 cm / < 18"	> 46 cm / < 18"	Blue
	XL	44.5 – 51 cm / 17.5" – 20"	< 46 cm / < 18"	> 46 cm / < 18"	Green
	XXL	51 – 58.4 cm / 20" – 23"	< 46 cm / < 18"	> 46 cm / < 18"	Purple
	XXXL	58.5 – 66 cm / 23" – 26"	< 46 cm / < 18"	> 46 cm / < 18"	Red
*Knee Length Sizing: 1 - Measure calf circumference, 2 - Measure length.					

VENAFLOW ANTI EMBOLISM THIGH LENGTH STOCKINGS SIZING CHART*							
	Size	1 - Thigh Circumference	2 - Calf Circumference	3 - Length			Color
				Short	Regular	Long	
	S	< 63.5 cm / 25"	< 30.5 cm / < 12"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Yellow
	M	< 63.5 cm / 25"	30.5 – 38 cm / 12" – 15"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	White
	L	< 63.5 cm / 25"	38 – 44.5 cm / 15" – 17.5"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Blue
	XL	63.5 – 81.5 cm / 25" – 31.8"	44.5 – 55 cm / 17.5" – 21.6"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Green
	XXL	81 – 91 cm / 31.8" – 35.8"	55 – 66 cm / 21.6" – 26"	< 74 cm / < 29"	74 – 84 cm / < 29" – 33"	> 84 cm / > 33"	Purple
*Thigh High Sizing: 1 - Measure Thigh Circumference, 2 - Measure Calf Circumference, 3 - Measure Length.							

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
500L-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, S, 12/BOX	\$103.58
500L-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, M, 12/BOX	\$103.58
500L-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, L, 12/BOX	\$103.58
500L-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, XL, 12/BOX	\$103.58
500L-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, XXL, 12/BOX	\$124.55
500L-7	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, LONG, XXXL, 12/BOX	\$124.55
500R-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, S, 12/BOX	\$103.58
500R-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, M, 12/BOX	\$103.58
500R-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, L, 12/BOX	\$103.58
500R-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, XL, 12/BOX	\$103.58
500R-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, XXL, 12/BOX	\$124.55
500R-7	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, BELOW-KNEE, W/INSPECTION TOE, REGULAR, XXXL, 12/BOX	\$124.55
501L-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, LONG, S, 6/BOX	\$55.57
501L-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, LONG, M, 6/BOX	\$55.57
501L-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, LONG, L, 6/BOX	\$55.57
501L-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, LONG, XL, 6/BOX	\$55.57
501L-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, LONG, XXL, 6/BOX	\$55.57
501R-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, REGULAR, S, 6/BOX	\$55.57
501R-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, REGULAR, M, 6/BOX	\$55.57
501R-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, REGULAR, L, 6/BOX	\$55.57
501R-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, REGULAR, XL, 6/BOX	\$55.57
501R-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, REGULAR, XXL, 6/BOX	\$55.57
501S-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, SHORT, S, 6/BOX	\$55.57
501S-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, SHORT, M, 6/BOX	\$55.57
501S-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, SHORT, L, 6/BOX	\$55.57
501S-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, SHORT, XL, 6/BOX	\$55.57
501S-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/INSPECTION TOE, SHORT, XXL, 6/BOX6PK	\$55.57
502L-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, LONG, L, 6/BOX	\$72.07
502L-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, LONG, XL, 6/BOX	\$72.07
502L-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, LONG, XXL, 6/BOX	\$72.07
502R-2	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, REGULAR, S, 6/BOX	\$72.07
502R-4	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, REGULAR, L, 6/BOX	\$72.07
502R-5	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, REGULAR, XL, 6/BOX	\$72.07
502R-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, REGULAR, XXL, 6/BOX	\$72.07
502S-3	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, SHORT, M, 6/BOX	\$72.07
502S-6	Enovis AIRCAST VENAFLOW ANTI-EMBOLISM STOCKINGS, THIGH LENGTH, W/WAIST ATTACHMENT & INSPECTION TOE, SHORT, XXL, 6/BOX	\$72.07

CERVICAL COLLARS

PROCARE®**ProCare® Transitional 172 Cervical Collar**

The ProCare® Transitional 172 Cervical Collars were designed to provide optimal support and comfort. The collar provides substantial motion restriction while minimizing painful pressure points that can lead to skin maceration and/or breakdown. Enhanced comfort using soft, breathable foam padding helps promote and encourage patient compliance. Removable open cell foam pads wick moisture away from the skin to ensure proper patient hygiene and comfort while minimizing skin irritation. X-ray, CT and MRI lucent.

Indication: Restricts cervical spine flexion, extension and rotation to promote patient recovery.

Adult Collar Only	Adult Collar w/ Repl. Pads	Height	Neck Circ.	Size
Adult Collars				
79-83283	79-83363	2.25"	13" – 22"	Adult Short
79-83285	79-83365	3.0"	13" – 22"	Adult Regular
79-83287	79-83367	3.75"	13" – 22"	Adult Tall
79-83288	79-83368	4.5"	13" – 22"	Adult X-Tall
Adult Collar Accessories				

Pediatric Collar Only	Pediatric Collar w/ Repl. Pads	Age	Height	Neck Circ.	Size
Pediatric Collars					
79-83272	79-83351	1-18 mo.	0.75"	8" – 12"	PED-1, Infant Short
79-83273	79-83352	9-24 mo.	1.25"	11" – 13"	PED-2, Infant
79-83275	79-83353	1-3 yrs.	2.0"	11.5" – 14"	PED-3, Child Toddler
79-83277	79-83354	2-5 yrs.	2.25"	12" – 14.5"	PED-4, Child Short
79-83278	79-83355	3-6 yrs.	2.5"	13" – 16"	PED-5, Child Regular
Pediatric Collar Accessories					





PROCARE®

Aspen® Vista® TX Collar

Vista TX encompasses all of the features and benefits of the Vista collar, but with a thoracic extension. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. x-ray, CT and MRI lucent.

Indication: Motion control of cervical spine flexion, extension and rotation to promote patient recovery.

Collar Only	Collar w/Repl. Pads	Neck Circum.	Size
79-83390	79-83360	13" - 19"	Adult/Univ.
Accessories			
79-83372	Vista TX Repl. Pad Set	Adult	
79-83377	Back Only	Adult	



PROCARE®

Universal Clinic Collar

Medium density foam collar provides consistent fit and support. Fits a wide range of neck sizes. Additional stockinette cover provided for better patient hygiene. Latex free.

Indication: Ideal for supporting the cervical spine in neutral position.

	Height	Length	Size
79-83520	2.5"	24"	Univ.
79-83500	3"	24"	Univ.
79-83510	4"	24"	Univ.



PROCARE®

Aspen® Vista® Colla

Innovative height adjustment technology offers six (6) collars in one. Two-piece polyethylene design enhances patient comfort and is easily adjusted with Velcro™ tabs. Anterior and posterior flex tabs on the front and back panels provide improved rigidity and durability. The large trachea opening provides functional access for the airway management and other essential procedures. Removable, foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation. x-ray, CT and MRI lucent.

Indication: Motion control of cervical spine flexion, extension and rotation to promote patient recovery.

Collar Only	Collar w/Repl. Pads	Neck Circum.	Size
79-83370	79-83380	13" - 19"	Adult/Univ.
Accessories			
79-83371	Vista Replacement Pad Set	Adult	
79-83377	Back Only	Adult	



PROCARE®

Form Fit Cervical Collar

Firm or medium density foam collar helps support cervical spine in neutral position. Contour design for use in either flexed or extended position. Cotton stockinette cover for better patient hygiene. Sold as each or in a pack of six (6). Latex free.

Indication: Ideal for post-operative rehabilitation; cervical sprains and strains.

Firm Density	Med. Density	Med. Dens. (6-pk.)	H" x L	Neck Circum.	Size
			2" x 14"	7" - 12"	Ped./2XS
			2" x 16.5"	9" - 13"	XS
79-83003	79-83013	79-83043	3" x 18.5"	11" - 16"	S
79-83004	79-83014	79-83044	3" x 22"	15" - 20"	S/Long
79-83005	79-83015	79-83045	4" x 20"	13" - 18"	M
79-83006	79-83016	79-83046	4" x 23"	16" - 21"	M/Long
79-83007	79-83017	79-83047	4.5" x 22.5"	15" - 20"	L
79-83008	79-83018	79-83048	4.5" x 24"	16" - 22"	XL
79-83009	79-83019		2" x 24"	17" - 22"	Long/Nar.

*Height = measurement from sternal notch to chin.



PROCARE®

Patriot® Collar

One-piece injection molded collar provides proper support and helps restrict cervical spine flexion, extension and rotation to promote patient recovery following trauma. Constructed with durable high density polyethylene lined with a soft closed cell foam for maximum patient comfort. Both universal sizes have multiple heights to fit most adults and pediatrics. Oversized trachea opening provides quick access to carotid pulse monitoring and emergency tracheotomies. Stacks flat to preserve space and is x-ray, CT and MRI lucent. Latex free.

Indication: Ideal for cervical spine immobilization and extraction.

	Height Adjustment*	Neck Circum.	Size
79-83130	1.5", 1.75", 2"	8" – 18"	Pediatric
79-83131	1.75", 2.5", 3", 3.5"	11" – 23"	Adult

*Height = measurement from sternal notch to chin.



PROCARE®

Miami Jr.® Collar

Designed to meet anatomical needs of children up to 12 years old. Padding helps retard skin breakdown.

Indication: Ideal for motion control of cervical spine flexion, extension and rotation to promote patient recovery.

	Age	Height*	Neck Circum.
Ped 0			
79-83281 MJ P0	0-6 mo.		
79-83206 MJR P0 w/repl. Pad	0-6 mo.		
79-83209 Repl. Pad Set Only			
Ped 1			
79-83251 MJ P1	6 mo. – 2 yrs.	1.125"	≤ 9"
79-83227 MJR P1 w/repl. pad	6 mo. – 2 yrs.	1.125"	≤ 9"
79-83253 Repl. Pad Set Only			
Ped 2			
79-83252 MJ P2	2 – 6 yrs.	1.25"	8" – 11"
79-83228 MJR P2 w/repl. pad	2 – 6 years	1.25"	8" – 11"
79-83254 Repl. Pad Set Only			
Ped 3			
79-83202 MJ P3	6 – 12 yrs.	1.75"	9" – 13"
79-83222 MJR P3 w/ repl. pad	6 – 12 yrs.	1.75"	9" – 13"

*Height = measurement from sternal notch to chin.



PROCARE®

Miami J® Collar

Designed to meet the needs of the adult ambulatory patient. Padding helps retard skin breakdown.

Indication: Ideal for motion control of cervical spine flexion, extension and rotation to promote patient recovery.

		Height*	Neck Circum.	Size
79-83201	MJ 200S	1"	10" – 20"	Super Short
79-83203	MJ 300	2"	10" – 20"	S
79-83205	MJ 400	2.5"	10" – 20"	M
79-83207	MJ 500	3"	10" – 20"	L
79-83208	MJ 200L	1.75"	14" – 24"	Stout
Standard Collar With Replacement Pad				
79-83221	MJR200S	1"	10" – 20"	Super Short
79-83232	MJR 250	1.25"	8" – 14"	XS
79-83233	MJR 300	2"	10" – 20"	S
79-83235	MJR 400	2.5"	0" – 20"	M
79-83237	MJR 500	3"	10" – 20"	L
79-83238	MJR 200L	1.75"	14" – 24"	Stout
Replacement Pads Only				
79-83224	MJP 100			Adult
79-83225	MJP 200L			Stout

*Height = measurement from sternal notch to chin

*Height = measurement from sternal notch to chin.



PROCARE®

One-Piece Flat Collar

Contoured plastic/foam design provides immobilization for extrication and rehabilitation. One piece construction with large trachea opening stacks flat to preserve space. Color coded closure straps for quick size reference. CT and MRI compatible.

Indication: Ideal for cervical spine immobilization and extraction.

	Color	Height*	Size
79-83021	Pink	1.5"	Infant
79-83022	Aqua	1.75"	Pediatric
79-83028	Purple	1.75"	Low
79-83023	Blue	2.5"	Short
79-83025	Orange	3"	Regular
79-83027	Green	3.5"	Tall

*Height = measurement from sternal notch to chin.

SHOULDER IMMOBILIZERS

**PROCARE®****Elastic Shoulder Immobilizer**

Elastic waist band construction for compressive support. Durable yet comfortable foam wrist and nylon humeral cuff with contact closure for ease of application. Fits right or left arm.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

		Rib Measurement	Size
79-84032	Male	24" – 32"	XS
79-84033	Male	32" – 36"	S
79-84035	Male	36" – 40"	M
79-84037	Male	40" – 44"	L
79-84038	Male	44" – 48"	XL
79-84039	Male	48" – 52"	2XL
79-84043	Female	24" – 28"	S
79-84045	Female	28" – 32"	M
79-84047	Female	32" – 36"	L
79-84048	Female	36" – 42"	XL
79-84049	Female	42" – 48"	2XL

**PROCARE®****Deluxe Sling and Swathe**

Cool, lightweight pile/foam construction provides patient comfort during rehabilitation following surgery or dislocation. Two over-the-shoulder criss-cross straps with a slide buckle for proper strap adjustment. Shoulder straps attach to the sling with contact closure to ensure proper weight distribution and ease of use. One size fits up to 54" waist. Fits right and left.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

		Size
79-84230		Univ.
79-99600	Swathe only	54" length

**PROCARE®****Universal Clavicle**

Foam/nylon construction with contact closure for easy adjustment. Movable D-rings allow for proper adjustment.

Indication: Ideal for clavicular fractures and postural problems.

	Circum.*	Size
79-85100	Up to 46"	Univ.

*Chest circumference measurement taken at bottom of



PROCARE®

Hot/Cold Therapy Wrap

Safe, non-toxic pack is microwave and freezer safe. Ideal for use on shoulder, wrist, elbow, back, knee and ankle.

Indication: Therapeutic wrap for muscle strains and sprains.

	Size
79-97210	Univ.
Replacement Hot/Cold Pack Only	
79-97220	



PROCARE®

Deluxe Arm Sling w/Pad

Generous size envelope to support the arm. Cotton/poly construction with web strap and contact closure. Thumb loop helps prevent sling from migrating. Shoulder pad included.

Indication: Ideal for cast support; injuries to the arm, wrist or hand.

	Depth	Length*	Size
79-840		8.25"	Pediatric/2XS
79-84002	7"	11"	Pediatric/XS
79-84003	7.5"	12.5"	S
79-84		15"	M
79-840		17.5"	L
79-84		20"	XL

*Measure length from elbow crease to base of index finger.



PROCARE®

Sling and Swathe

Sling and swathe constructed with durable cotton/duck material. Slide buckle closure on soft cotton woven shoulder strap. Available in four (4) sizes.

Indication: Ideal for immobilization of the shoulder after surgery, injury or cast removal.

	Swathe Size	Envelope	Size
79-84243	4" x 44"	9" x 14"	S
79-84245	4" x 49"	10" x 16"	M
79-84247	4" x 54"	11" x 17"	L
79-84248	4" x 59"	12" x 19"	XL



PROCARE®

Universal Surround® Elbow w/FLOAM™

Comfortable neoprene strap features FLOAM bladder and wraparound universal design with loop-lock closure. FLOAM is intended to reduce impact stress to the epicondyle and helps dissipate muscle force over a greater area. Durable to withstand wear during exercise, sporting activities, and repetitive work conditions.

Indication: Ideal for bursitis, sprains, strains, tendonitis, Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.*	Size
79-82570	≤ 18"	Univ.

*Forearm circumference.



PROCARE®

ElbowRANGER® Motion Control Splint

For non-operative application, collateral ligament strains and ruptures, and ulnar nerve transposition. May be used following early cast removal to stabilize fractures, and for post-operative immobilization. Removable hand grip assists in the prevention of pronation/supination. Aluminum R.O.M. bilateral hinges allow adjustments in 30° increments with formed Kydex® shells. Universal application fits right or left arm.

Indication: Ideal for post-operative immobilization and protected range-of-motion.

	Circum.*	Size
79-94200	8" – 11.5"	S
79-94210	11.5" – 14.5"	M
79-94220	14.5" – 19"	L

*Measurement at mid-bicep.



AIRCAST®

Pneumatic Armband

A single, pre-inflated aircell concentrates compression on the extensor muscle, not around the arm—providing more support, less constriction. Available in black and beige. Universal fit.

Indication: Ideal for Tennis elbow; Golfers elbow (epicondylitis).

	Circum.*	Color	Size
05A-B	8" – 14"	Black	Univ.
05A	8" – 14"	Beige	Univ.

PROCARE®

Tennis Elbow Support

Combines soothing warmth of neoprene elbow sleeve with a counterforce strap that restricts full extension of muscles and tendons. Contact closure strap with D-ring provides exact fit and adjustment.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.*	Length	Size
79-82322	8" – 9"	8"	XS
79-82323	9" – 10"	8"	S
79-82325	10" – 12"	8"	M
79-82327	12" – 14"	8"	L
79-82328	14" +	8"	XL

*Forearm circumference.

*Forearm circumference.



PROCARE®

Arm Band w/Compression Pad

Quality nylon construction with reinforced neoprene pad to alleviate stress over lateral/medial epicondyle to address tendonitis, bursitis, and strains. D-ring with contact closure for adjustable compression. One size fits all.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.*	Size
79-97000	≤ 18"	Univ.

*Forearm circumference.



PROCARE®

Elbow Sleeve

Neoprene sleeve provides maximum warmth and comfort without roll or ride-up. Offers support and relief to elbow joint, muscles and tendons. 8.5" length.

Indication: Ideal for mild support and compression of the elbow.

	Circum.*	Size
79-82312	8" – 9"	XS
79-82313	9" – 10"	S
79-82315	10" – 12"	M
79-82317	12" – 14"	L
79-82318	14" – 16"	XL
79-82319	16" – 18"	2XL

*Forearm circumference.



PROCARE®

Tennis Elbow Support w/FLOAM™

Combines neoprene band with FLOAM bladder for compression over tendon insert point. FLOAM has no memory, which prohibits pressure from being displaced to bony prominences.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Circum.*	Size
79-81183	6" – 8"	S
79-81185	8" – 11"	M
79-81187	11" – 14"	L
79-81188	14" – 16"	XL

*Forearm circumference.



PROCARE®

Clinic Tennis Elbow

Two-inch wide nylon strap with D-ring closure to allow comfortable fit and compression to address lateral/medial epicondylitis and tendonitis. One size fits most patients.

Indication: Ideal for Medial epicondylitis (Golfer's elbow) and Lateral epicondylitis (Tennis elbow).

	Qty.	Size
79-81031	6/pk.	Univ.

WRIST AND HAND SUPPORTS



PROCARE®

ComfortFORM™ Wrist

Durable lightweight foam laminate is Lycra® lined for breathability and patient comfort. Preformed aluminum stay and loop/lock closure helps provide anatomically correct fit and proper support.

Indication: Ideal for sprains, strains and control of Carpal Tunnel syndrome symptoms.

Right	Wrist Circum.	Length	Size
79-87281 79-87291	≤ 4.5"	7.0"	Ped./2XS



PROCARE®

CTS Wrist Support

Contoured cotton/elastic construction provides comfortable compression for proper support following mild sprains, strains and cast removal. Circumferential strap for improved wrist stabilization. Preformed aluminum stay provides anatomically correct fit. Contact closure.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right	Left	Wrist Circum.	Size
79-87152	79-87162	≤ 5.5"	Ped./XS
79-87153	79-87163	5.5" – 6.5"	S
79-87155	79-87165	6.5" – 7.5"	M



PROCARE®

Elastic Wrist Brace

Soft, durable cotton/elastic construction provides compressive support. Contoured elastic body and preformed aluminum stay for anatomically correct fit and proper support following mild sprain, strain or cast removal. Contact closure straps for ease of application. Designed to prevent brace from extending beyond palmar crease allowing full finger range-of-motion.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

Right		Circum.	Length	Size
79-87072	79-87082	≤ 5.5"	7"	Ped./XS
79-87073	79-87083	5.5" – 6.5"	7"	S
79-87075	79-87085	6.5" – 7.5"	7"	M
79-87077	79-87087	7.5" – 8.5"	7"	L
79-87078	79-87088	8.5"+	7"	XL

**PROCARE®****Universal Elastic Wrist Brace**

Cotton/elastic construction provides compressive support. Contoured elastic body and aluminum stay intended to provide anatomically correct fit. Converts from right to left by changing location of aluminum stay. Contact closure straps for easy adjustment and application.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

	Wrist Circum.	Size
79-87092	≤ 5.5"	Ped./XS
79-87093	5.5" – 6.5	S
79-87095	6.5" – 7.5	M
79-87097	7.5" – 8.5"	L
79-87098	8.5"+	XL

**PROCARE®****Universal CTS Wrist Brace**

Panel elastic construction. Converts from right to left by moving aluminum palmar stay. Circumferential wrist strap provides improved stabilization. Contact closure.

Indication: Ideal for Carpal Tunnel Syndrome and immobilization of the wrist.

	Wrist Circum.	Size
79-97012	≤ 5.5"	Ped./XS
79-97013	5.5" – 6.5"	S
79-97015	6.5" – 7.5"	M
79-97017	7.5" – 8.5"	L
79-97018	8.5" – 9.5"	XL
79-97019	9.5" – 10.5"	2XL

**PROCARE®****Canvas B.A.T.H.™ Wrist Splint**

Constructed of heavy duty canvas for durability. Splint features a double contour stay for precise fit. Slip on design for easy application. Does not extend past palmar crease which does not limit finger dexterity.

Indication: Ideal for wrist sprains or strains; Carpal Tunnel Syndrome and Tendonitis.

Right		Wrist Circum.	Length	Size
79-87441	79-87481	3.5" – 4.5"	5.25"	Ped.
79-87444	79-87484	4" – 5"	5.25"	Child
79-87442	79-87482	4.5" – 5.5"	7.5"	XS
79-87443	79-87483	5.5" – 6.5"	8"	S
79-87445	79-87485	6.5" – 7.5"	8"	M
79-87447	79-87487	7.5" – 8.5"	8"	L
79-87448	79 488	8.5"+	8"	XL

**PROCARE®****Universal Wrist-O-Prene™**

Wraparound neoprene brace provides compression and support for treatment of Carpal Tunnel Syndrome (CTS) and wrist sprains and strains. Removable contoured aluminum palmar stay. Inside contact closure allows for ease of application and proper fit. Durable construction withstands wear during exercise, sporting activities and repetitive work conditions.

Indication: Ideal for compression and support for treatment of Carpal Tunnel Syndrome (CTS) and wrist sprains and strains.

Right	Wrist Circum.	Size
79-82	≤ 11"	Univ.



PROCARE®

Universal Wrist & Forearm Supports

Breathable nylon/fiber with flannelette liner. Loop/lock closure. Adjustable aluminum palmar and dorsal stays for support. Available in 7" and 10" length.

Indication: Ideal for post surgical immobilization, post fracture casting, wrist sprains and strains.

Right	Left	Wrist Circum.	Length	Size
79-87050	79-87060	≤ 9"	10"	Univ.
79-87070	79-87080	≤ 9"	7"	Univ.



PROCARE®

Safety Wrist

Ergonomically designed brace allows for full hand dexterity during repetitive daily activities. Durable, cotton/elastic construction, contact closure straps and plastic palmer stay provide support to help address wrist sprains, tendonitis, and Carpal Tunnel Syndrome.

Indication: Ideal for wrist sprains, tendonitis and Carpal Tunnel Syndrome.

Right	Left	Wrist Circum.	Length	Size
79-82223	79-82233	5.5" – 6.5"	6"	S
79-82225	79-82235	6.5" – 7.5"	6"	M
79-82227	79-82237	7.5" – 8.5"	6"	L
79-82228	79-82238	8.5" – 9.5"	6"	XL
79-82229	79-82239	9.5"+	6"	2XL



PROCARE®

ThumbSPICA™

Lightweight, perforated outer foam shell with a cool, cotton-terry liner. Elastic, circumferential, contact closure straps for proper positioning of splint.

Indication: Ideal for soft tissue injury to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

Right	Left	Length	Size®
79-87113	79-87114	9"	S/M
79-87117	79-87118	9"	L/XL

®S/M - Fits most women and smaller men, L/XL - Fits most men and larger women.



PROCARE®

Abducted ThumbSPICA™

Volara™ foam laminate construction with elastic contact closure straps.

Indication: Ideal for soft tissue injury to the thumb, arthritis, tendonitis, deQuervain's Syndrome.

Right	Left	Size
79-87100	79-87101	Univ.



PROCARE®

Universal Cock-Up Splint

Two piece design constructed of durable perforated vinyl, foam, and tricot. Available in 7" and 10" length. May be worn on either hand. Dorsal and palmar aluminum stays provide proper patient support.

Indication: Ideal for immobilization of the wrist and post-fracture casting.

7"	10"	Size
79-87010	79-87020	Univ.
79-87021		Univ./XL



PROCARE®

Wrist Wrap

Adjustable neoprene wraparound design provides a secure compressive support. Fits right or left. Two styles.

Indication: Ideal for compressive support for weak or injured wrist.

		Size
79-82050	Wistlet	Univ.
79-82300	Nylon 2	Univ.



PROCARE®

ComfortFORM™ Wrist w/Abducted Thumb

Durable lightweight fabric is Lycra® lined for breathability and patient comfort. Dual aluminum stays support palmar surface of wrist and extensor surface of thumb. Stays can be adjusted for proper angulation. Contact closure straps.

Indication: Ideal for wrist sprains, strains, scaphoid injuries, Carpal Tunnel, Gamekeeper's thumb, and deQuervain's Syndrome symptoms.

Right		Wrist Circum.	Length	Size
79-87302	79-87312	≤ 5.5"	7.5"	Ped./XS
79-87303	79-87313	5.5" – 6.5"	7.5"	S
79-87305	79-87315	6.5" – 7.5"	8"	M
79-87307	79-87317	7.5" – 8.5"	8.5"	L
79-87308	79-87318	8.5"+	8.5"	XL



PROCARE®

Wrist/Thumb Wrap

Neoprene construction provides compressive wrist support with thumb abduction. One size fits most patients. Right or left.

Indication: Ideal for soft tissue injuries to the thumb and wrist.

	Size
79-82600	Univ.

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
81-05A	Enovis AIRCAST PNEUMATIC ARMBAND	\$17.91
81-05AB	Enovis AIRCAST PNEUMATIC ARMBAND-BLACK	\$17.91
79-81183	Enovis PROCARE TENNIS ELBOW SUPPORT,S	\$18.34
79-81185	Enovis PROCARE TENNIS ELBOW SUPPORT,M	\$18.34
79-81187	Enovis PROCARE TENNIS ELBOW SUPPORT,L	\$18.34
79-81188	Enovis PROCARE TENNIS ELBOW SUPPORT,XL	\$18.84
79-82223	Enovis PROCARE SAFETY WRIST,RT,S	\$18.48
79-82300	Enovis PROCARE UNIVERSAL WRIST WRAP	\$6.52
79-82312	Enovis PROCARE ELBOW SLEEVE,XS	\$11.82
79-82313	Enovis PROCARE ELBOW SLEEVE,S	\$11.82
79-82315	Enovis PROCARE ELBOW SLEEVE,M	\$11.82
79-82317	Enovis PROCARE ELBOW SLEEVE,L	\$11.82
79-82318	Enovis PROCARE ELBOW SLEEVE,XL	\$11.82
79-82319	Enovis PROCARE ELBOW SLEEVE,XXL	\$12.01
79-82322	Enovis PROCARE TENNIS ELBOW SUPPORT,XS	\$15.13
79-82323	Enovis PROCARE TENNIS ELBOW SUPPORT,S	\$15.13
79-82325	Enovis PROCARE TENNIS ELBOW SUPPORT,M	\$15.13
79-82327	Enovis PROCARE TENNIS ELBOW SUPPORT,L	\$15.13
79-82328	Enovis PROCARE TENNIS ELBOW SUPPORT,XL	\$15.13
79-82329	Enovis PROCARE TENNIS ELBOW SUPPORT,XXL	\$15.13
79-82570	Enovis PROCARE SURROUND ELBOW, UNIV.	\$14.85
79-82600	Enovis PROCARE WRIST/THUMB WRAP,UNIV	\$16.00
79-83003	Enovis PROCARE FORM FIT CERVICAL COLLAR,S	\$9.55
79-83004	Enovis PROCARE FORM FIT CERVICAL COLLAR,SM/LONG	\$9.55
79-83005	Enovis PROCARE FORM FIT CERVICAL COLLAR,M	\$9.55
79-83006	Enovis PROCARE FORM FIT CERVICAL COLLAR,M/LONG	\$9.55
79-83007	Enovis PROCARE FORM FIT CERVICAL COLLAR,L	\$9.55
79-83008	Enovis PROCARE FORM FIT CERVICAL COLLAR,XL	\$9.55
79-83011	Enovis PROCARE FORM FIT CERVICAL COLLAR, XXS	\$7.49
79-83012	Enovis PROCARE FORM FIT CERVICAL COLL (XS)	\$7.49
79-83013	Enovis PROCARE FORM FIT CERVICAL COLL (S)	\$7.49
79-83014	Enovis PROCARE FORM FIT CERV COLL (SM/LONG)	\$7.49
79-83015	Enovis PROCARE FORM FIT CERVICAL COLL (M)	\$7.49
79-83016	Enovis PROCARE FORM FIT CERV COLL (M-LONG)	\$7.49
79-83017	Enovis PROCARE FORM FIT CERVICAL COLLAR (L)	\$7.49
79-83018	Enovis PROCARE FORM FIT CERV COLL (XL)	\$7.49
79-83019	Enovis PROCARE FORM FIT,CERV COLL,(X-LNG,NAR)	\$7.49
79-83022	Enovis PROCARE ONE-PIECE COLLAR,PEDIATRIC	\$25.63

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-83025	Enovis PROCARE ONE-PIECE COLLAR,REGULAR	\$25.63
79-83027	Enovis PROCARE ONE-PIECE COLLAR,TALL	\$25.63
79-83130	Enovis PROCARE ONE PIECE CERVICAL COLLAR, PEDIATRIC	\$9.84
79-83131	Enovis PROCARE ONE PIECE CERVICAL COLLAR, ADULT	\$9.84
79-83201	Enovis PROCARE XTEND 174, SUPER SHORT	\$72.25
79-83202	Enovis PROCARE XTEND 174 PED 3	\$72.25
79-83203	Enovis PROCARE XTEND 174, SHORT	\$72.25
79-83204	Enovis PROCARE XTEND 174, PED 3 REPL PAD SET	\$27.47
79-83205	Enovis PROCARE XTEND 174, REGULAR	\$72.25
79-83206	Enovis PROCARE XTEND 174, PED 0 W/REPL PAD	\$101.06
79-83207	Enovis PROCARE XTEND 174, TALL	\$72.25
79-83208	Enovis PROCARE XTEND 174, STOUT	\$72.25
79-83214	Enovis PROCARE ACC TRANSITIONAL 172 POSTERIOR PANEL, PED-1/PED-2	\$21.73
79-83221	Enovis PROCARE XTEND 174, SUPER SHORT W/REPL PAD	\$101.06
79-83222	Enovis PROCARE XTEND 174, PED 3 W/REPL PAD	\$101.06
79-83224	Enovis PROCARE XTEND 174 STANDARD REPL PAD SET	\$27.47
79-83225	Enovis PROCARE XTEND 174 STOUT REPL PAD SET	\$27.47
79-83227	Enovis PROCARE XTEND 174, PED 1 W/REPL PAD	\$101.06
79-83228	Enovis PROCARE XTEND 174, PED 2 W/REPL PAD	\$101.06
79-83232	Enovis PROCARE XTEND 174, XS W/REPL PAD	\$101.06
79-83233	Enovis PROCARE XTEND 174, SHORT W/REPL PAD	\$101.06
79-83235	Enovis PROCARE XTEND 174, REGULAR W/REPL PAD	\$101.06
79-83237	Enovis PROCARE XTEND 174, TALL W/REPL PAD	\$101.06
79-83238	Enovis PROCARE XTEND 174, STOUT W/REPL PAD	\$101.06
79-83251	Enovis PROCARE XTEND 174, PED 1	\$72.25
79-83252	Enovis PROCARE XTEND 174, PED 2	\$72.25
79-83253	Enovis PROCARE XTEND 174, PED 1 REPL PAD SET	\$27.47
79-83254	Enovis PROCARE XTEND 174, PED 2 REPL PAD SET	\$27.47
79-83272	Enovis PROCARE TRANSITIONAL 172 COLLAR, PED-1	\$45.18
79-83273	Enovis PROCARE TRANSITIONAL 172 COLLAR, PED-2	\$45.18
79-83275	Enovis PROCARE TRANSITIONAL 172 COLLAR, PED-3	\$45.18
79-83277	Enovis PROCARE TRANSITIONAL 172 COLLAR, PED-4	\$45.18
79-83278	Enovis PROCARE TRANSITIONAL 172 COLLAR, PED-5	\$45.18
79-83281	Enovis PROCARE XTEND 174, PED 0	\$72.25
79-83283	Enovis PROCARE TRANSITIONAL 172 COLLAR, ADULT SHORT	\$45.30
79-83285	Enovis PROCARE TRANSITIONAL 172 COLLAR, ADULT REGULAR	\$45.30
79-83287	Enovis PROCARE TRANSITIONAL 172 COLLAR, ADULT TALL	\$45.30
79-83288	Enovis PROCARE TRANSITIONAL 172 COLLAR, ADULT X-TALL	\$45.30

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-83325	Enovis PROCARE ACC TRANSITIONAL 172 COLLAR PAD SET, ADULT	\$19.24
79-83370	Enovis PROCARE VISTA COLLAR, ADULT	\$72.35
79-83371	Enovis PROCARE VISTA COLLAR REPL PADS, ADULT	\$34.23
79-83372	Enovis PROCARE VISTATX REPL PADS, ADULT	\$34.08
79-83377	Enovis PROCARE VISTA COLLAR BACK, LARGE	\$42.38
79-83380	Enovis PROCARE VISTA COLLAR SET, ADULT	\$101.30
79-83390	Enovis PROCARE VISTA COLLAR TX, ADULT	\$70.19
79-83500	Enovis PROCARE 3" UNIVERSAL CLINIC COLLAR	\$4.80
79-83510	Enovis PROCARE 4" UNIVERSAL CLINIC COLLAR	\$5.07
79-83520	Enovis PROCARE 2.5" UNIVERSAL CERVICAL COLLAR	\$6.68
79-84001-900	Enovis PROCARE DELUXE ARM SLING, XXS	\$9.69
79-84002-900	Enovis PROCARE DELUXE ARM SLING, XS	\$9.69
79-84003-900	Enovis PROCARE DELUXE ARM SLING, S	\$9.69
79-84005-900	Enovis PROCARE DELUXE ARM SLING, M	\$9.69
79-84007-900	Enovis PROCARE DELUXE ARM SLING, L	\$9.69
79-84008-900	Enovis PROCARE DELUXE ARM SLING, XL	\$9.69
79-84020	Enovis PROCARE FOAM ARM SLING PAD	\$3.12
79-84032	Enovis PROCARE ELASTIC SHOULDER IMMOB MALE,XS	\$14.98
79-84033	Enovis PROCARE ELASTIC SHOULDER IMMOB-MALE,S	\$14.98
79-84035	Enovis PROCARE ELASTIC SHOULDER IMMOB-MALE,M	\$14.98
79-84037	Enovis PROCARE ELASTIC SHOULDER IMMOB-MALE,L	\$14.98
79-84038	Enovis PROCARE ELASTIC SHOULDER IMMOB-MALE,XL	\$14.98
79-84039	Enovis PROCARE ELASTIC SHLDER IMMOB-MALE,XXL	\$14.98
79-84043	Enovis PROCARE ELASTIC SHOULDER IMMOB-FEM,S	\$14.95
79-84045	Enovis PROCARE ELASTIC SHOULDER IMMOB-FEM,M	\$14.95
79-84047	Enovis PROCARE ELASTIC SHOULDER IMMOB-FEM,L	\$14.95
79-84048	Enovis PROCARE ELASTIC SHOULDER IMMOB-FEM,XL	\$14.95
79-84049	Enovis PROCARE ELASTIC SHOULDER IMMOB-FEM,XXL	\$14.95
79-84230	Enovis PROCARE DELUXE SLING & SWATHE - UNIV	\$14.48
79-84243	Enovis PROCARE SLING & SWATHE (SM)	\$28.63
79-84245	Enovis PROCARE SLING & SWATHE (MED)	\$28.63
79-84247	Enovis PROCARE SLING & SWATHE (LG)	\$28.63
79-84248	Enovis PROCARE SLING & SWATHE (XL)	\$28.63
79-85100	Enovis PROCARE UNIVERSAL CLAVICLE SPLINT	\$14.14
79-87010	Enovis PROCARE 7 UNIVERSAL WRIST	\$21.53
79-87050	Enovis PROCARE UNIV WRIST/FOREARM SUPP-RIGHT	\$14.04
79-87070	Enovis PROCARE 7 UNIV. WRIST SUPPORT/RIGHT	\$19.03
79-87072-004	Enovis PROCARE ELASTIC WRIST, RT XS	\$5.88

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-87073	Enovis PROCARE ELASTIC WRIST,RT,S	\$7.41
79-87073-004	Enovis PROCARE ELASTIC WRIST,RT,S	\$10.29
79-87075	Enovis PROCARE ELASTIC WRIST,RT,M	\$7.41
79-87075-004	Enovis PROCARE ELASTIC WRIST,RT,M	\$10.29
79-87077	Enovis PROCARE ELASTIC WRIST,RT,L	\$7.41
79-87077-004	Enovis PROCARE ELASTIC WRIST,RT,L	\$10.29
79-87078	Enovis PROCARE ELASTIC WRIST,RT,XL	\$7.41
79-87078-004	Enovis PROCARE ELASTIC WRIST, RT XL	\$10.99
79-87082	Enovis PROCARE ELASTIC WRIST,LT,XS	\$7.41
79-87082-004	Enovis PROCARE ELASTIC WRIST, LT XS	\$5.88
79-87083	Enovis PROCARE ELASTIC WRIST,LT,S	\$7.41
79-87083-004	Enovis PROCARE ELASTIC WRIST,LT,S	\$10.29
79-87085	Enovis PROCARE ELASTIC WRIST,LT,M	\$7.41
79-87085-004	Enovis PROCARE ELASTIC WRIST,LT,M	\$10.29
79-87087	Enovis PROCARE ELASTIC WRIST,LT,L	\$7.41
79-87087-004	Enovis PROCARE ELASTIC WRIST,LT,L	\$10.29
79-87088	Enovis PROCARE ELASTIC WRIST,LT,XL	\$7.41
79-87088-004	Enovis PROCARE ELASTIC WRIST, LT XL	\$10.99
79-87100	Enovis PROCARE ABDUCTED THUMB SPICA ,UNIV,RT	\$14.08
79-87101	Enovis PROCARE ABDUCTED THUMB SPICA,UNIV,LEFT	\$14.08
79-87113	Enovis PROCARE THUMB SPICA,RT,S/M	\$20.44
79-87114	Enovis PROCARE THUMB SPICA,LT,S/M	\$20.44
79-87117	Enovis PROCARE THUMB SPICA,RT,L/XL	\$20.44
79-871108	Enovis PROCARE THUMB SPICA,LT,L/XL	\$20.44
79-87152	Enovis PROCARE CTS WRIST SUPPORT (XS, RIGHT)	\$9.83
79-87153	Enovis PROCARE CTS WRIST SUPPORT (SM, RIGHT)	\$9.83
79-87155	Enovis PROCARE CTS WRIST SUPPORT (MED, RIGHT)	\$9.83
79-87157	Enovis PROCARE CTS WRIST SUPPORT (LG, RIGHT)	\$9.83
79-87158	Enovis PROCARE CTS WRIST SUPPORT (XL, RIGHT)	\$9.83
79-87162	Enovis PROCARE CTS WRIST SUPPORT (XS, LEFT)	\$9.83
79-87163	Enovis PROCARE CTS WRIST SUPPORT (SM, LEFT)	\$9.83
79-87165	Enovis PROCARE CTS WRIST SUPPORT (MED, LEFT)	\$9.83
79-87167	Enovis PROCARE CTS WRIST SUPPORT (LG, LEFT)	\$9.83
79-87168	Enovis PROCARE CTS WRIST SUPPORT (XL, LEFT)	\$9.83
79-87302	Enovis PROCARE COMFORT FORM WRIST/THUMB,RT,XS	\$27.57
79-87303	Enovis PROCARE COMFORT FORM WRIST/THUMB RT,S	\$27.57
79-87305	Enovis PROCARE COMFORT FORM WRIST/THUMB,RT,M	\$27.57
79-87307	Enovis PROCARE COMFORT FORM WRIST/THUMB,RT,L	\$27.57

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-87308	Enovis PROCARE COMFORT FORM WRIST/THUMB,RT,XL	\$27.57
79-87312	Enovis PROCARE COMFORT FORM WRIST/THUMB,LT,XS	\$27.57
79-87313	Enovis PROCARE COMFORT FORM WRIST/THUMB,LT,S	\$27.57
79-87315	Enovis PROCARE COMFORT FORM WRIST/THUMB,LT,M	\$27.57
79-87317	Enovis PROCARE COMFORT FORM WRIST/THUMB,LT,L	\$27.57
79-87318	Enovis PROCARE COMFORT FORM WRIST/THUMB,LT,XL	\$27.57
79-87441	Enovis PROCARE B.A.T.H. WRIST SPLINT, R, PED	\$16.60
79-87442	Enovis PROCARE B.A.T.H. WRIST SPLINT, R, XS	\$16.60
79-87443	Enovis PROCARE B.A.T.H. WRIST SPLINT, R, S	\$16.60
79-87444	Enovis PROCARE B.A.T.H. WRIST SPLINT, R,CHILD	\$16.60
79-87481	Enovis PROCARE B.A.T.H. WRIST SPLINT, L, PED	\$16.60
79-87482	Enovis PROCARE B.A.T.H. WRIST SPLINT, L, XS	\$16.60
79-87483	Enovis PROCARE B.A.T.H. WRIST SPLINT, LEFT, S	\$16.60
79-87484	Enovis PROCARE B.A.T.H. WRIST SPLINT, L,CHILD	\$16.60
79-87488	Enovis PROCARE B.A.T.H. WRIST SPLINT, L, XL	\$25.34
79-94200	Enovis PROCARE ELBOW RANGER ,S	\$105.94
79-94210	Enovis PROCARE ELBOW RANGER,M	\$105.94
79-94220	Enovis PROCARE ELBOW RANGER,L	\$105.94
79-97000	Enovis PROCARE ARMBAND - UNIV	\$7.50
79-97012	Enovis PROCARE BLACK ELASTIC WRIST BRACE-XS	\$8.90
79-97013	Enovis PROCARE BLACK ELASTIC WRIST BRACE-S	\$8.90
79-97015	Enovis PROCARE BLACK ELASTIC WRIST BRACE-M	\$8.90
79-97017	Enovis PROCARE BLACK ELASTIC WRIST BRACE-L	\$8.90
79-97018	Enovis PROCARE BLACK ELASTIC WRIST BRACE-XL	\$8.90
79-97019	Enovis PROCARE BLACK ELASTIC WRIST BRACE-XXL	\$8.90
79-97210	Enovis PROCARE HOT/COLD THERAPY PACK	\$41.50

LUMBROSACRAL SUPPORTS

**PROCARE®****Sacro-Lumbar Support w/Compression Straps**

Low profile 9" elastic support with double side pull compression straps. Neoprene posterior panel accommodates ComfortFORM™ Moldable Insert #79-89120. Foam lumbar pad is included with 79-8250X style. No stays.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

Black	White	Waist Meas.	Size
79-82502		18" – 24"	XS
79-82503	79-89263	24" – 33"	S
79-82505	79-89265	33" – 39"	M
79-82507	79-89267	39" – 45"	L
79-82508	79-89268	45" – 53"	XL
79-82509	79-89269	53" – 59"	2XL
79-82509-10		59" – 65"	3XL
79-82509-11		65" – 72"	4XL

**PROCARE®****Clinic Retention Support with Compression Straps**

Foam/pile, with 9" panel elastic construction and double side pull elastic compression straps. Two sizes. Accommodates ComfortFORM™ moldable insert #79-89120.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89240	24" – 40"	S/M
79-89241	40" – 54"	L/XL

**PROCARE®****10" Criss-Cross Support**

Universal style criss-cross elastic support with two plastic posterior stays.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89040	28" – 50"	Univ.

**PROCARE®****Criss-Cross Support with Compression Straps**

10" elastic criss-cross support with compression straps. Four 9" stays—two aluminum that can be easily formed to help ensure proper support and two plastic to maintain position of support.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89183	24" – 30"	S
79-89185	30" – 36"	M
79-89187	36" – 42"	L
79-89188	42" – 48"	XL
79-89189	48" – 52"	2XL

**PROCARE®****10" Double-Pull Sacro-Lumbar Support**

Premium all-elastic construction with compression straps. Four 9" posterior plastic stays for stability and support.

Indication: Ideal for providing support and compression for strains, sprains of the low back and abdominal area.

	Waist Measurement	Size
79-89003	31" – 34"	S
79-89005	35" – 38"	M
79-89006	36" – 42"	M/L
79-89007	39" – 42"	L
79-89008	43" – 46"	XL
79-89009	47" – 54"	2XL
79-89009-10	55" – 60"	3XL



PRO CARE®

Industrial Back Support

8" black, mesh back support with contact closures. Double-pull elastic compression suspenders. Choice between sewn-on or optional attachable shoulder suspension straps.

Indication: Ideal for helping to prevent and protect back injuries.

Black w/o Suspenders	Blk. w/Suspenders	Waist	Size
	79-89142	22" – 24"	XS
79-89083	79-89143	24" – 30"	S
79-89085	79-89145	30" – 36"	M
79-89087	79-89147	36" – 42"	L
79-89088	79-89148	42" – 50"	XL
79-89089	79-89149	50" – 56"	2XL
	79-89149L	56" – 64"	3XL
	79-89149X	64" – 72"	4XL
79-89190	Optional Attachable Suspenders		Univ.



PRO CARE®

Industrial Back Support w/Compression Pad

Elastic construction with 8" molded back compression pad. Two 6" plastic stays to ensure proper support. Contact closure.

Indication: Ideal for helping to prevent and protect back injuries.

	Waist Measurement	Size
79-97033	24" – 30"	S
79-97035	30" – 36"	M
79-97037	36" – 42"	L
79-97038	42" – 50"	XL
79-89350	Suspenders	Univ.



PRO CARE®

Clinic Industrial Back Support

Lightweight, low profile 8" back support with double-pull elastic compression straps and four sewn-in 7" plastic stays. Contact closure. Sewn-on adjustable suspenders.

Indication: Ideal for helping to prevent and protect back injuries.

Right	Left	Size
79-89172	18" – 24"	XS
79-89173	24" – 30"	S
79-89175	30" – 36"	M
79-89177	36" – 42"	L
79-89178	42" – 50"	XL
79-89179	50" – 56"	2XL



PRO CARE®

Maternity Belt

8" plush elastic belt helps relieve discomfort of lower abdomen and back during all stages of pregnancy. Abdominal strap supports stomach while elastic back provides compression and support for lumbar area. Optional 8" ComfortFORM™ moldable back insert #79-89120 provides additional support.

Indication: Ideal for helping to relieve lower back discomfort during pregnancy.

	Hip Measurement*	Size
79-89293	32" – 42"	S
79-89295	42" – 52"	M
79-89297	52" – 62"	L
79-89298	62" – 72"	XL

*Measurement taken during pregnancy.

ABDOMINAL & RIB SUPPORTS



PRO-CARE®

Abdominal Binder

Elastic binder with four sewn-in plastic stays provides abdominal and lumbar support. Abdominal, post-natal or post-surgical application. Easy adjustment with contact closure. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

10"		14"	Waist Meas.	Size
79-89043	79-89323	79-89333	24" – 30"	S
79-89045	79-89325	79-89335	30" – 36"	M
79-89047	79-89327	79-89337	36" – 42"	L
79-89048	79-89328	79-89338	42" – 48"	XL
79-89049	79-89329	79-89339	48" – 54"	2XL



PRO-CARE®

Personal Abdominal Binder

Constructed of multi-panel, vented elastic. Multi-panel helps binder conform to patient's body. 10" height is an all elastic panel. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

	Height	Waist Meas.	Size
79-89364	9"	20" – 42"	S/M
79-89366	9"	36" – 65"	M/L
79-89376	10"	32" – 65"	M/L
79-89384	12"	20" – 42"	S/M
79-89386	12"	36" – 65"	M/L



PRO-CARE®

10" Universal Abdominal Support

10" elastic and pile/foam laminate with three 9" removable plastic stays. Contact closure to ensure ease of application and removal. Latex free.

Indication: Ideal for providing compression and support for strains and weakness of the abdominal area.

	Waist Meas.*	Size
79-89080	28" – 50"	Univ.

*Measurements taken from male waist or female hip.



PRO-CARE®

Universal Deluxe Rib Belt

6" two-panel elastic belt with flannel lining and contact closure. Choice of male or female styles. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

Male	Female	Rib Meas.	Qty.	Size
79-89150	79-89160	28" – 50"		Univ.
79-99040	79-99030	28" – 50"	6/pk.	Univ.
79-89151		50" – 65"		XL



PRO-CARE®

Universal Rib Belt

6" elastic and pile/foam construction. Contact closure. Choice of male or female styles. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

		Rib Meas.	Size
79-89050	Male	24" – 50"	Univ.
79-89060	Female	24" – 50"	Univ.



PRO-CARE®

Deluxe Rib Belt

6" elastic and foam/flannel construction contoured for male or female styles. Contact closure for ease of application. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

Male	Rib Meas.	Female	Rib Meas.	Size
79-89063	32" – 36"	79-89073	24" – 27"	S
79-89065	36" – 40"	79-89075	27" – 33"	M
79-89067	40" – 46"	79-89077	33" – 37"	L
79-89068	46" – 50"	79-89078	37" – 44"	XL
79-89069	50" – 54"	79-89079	44" – 50"	2XL



PRO-CARE®

8" Universal Rib Belt

Panel elastic with foam/flannel lining. Contact closure for proper adjustment. Latex free.

Indication: Ideal for providing compression and support for fractures and strains to the rib cage.

		Rib Meas.	Size
79-89100	Male	28" – 50"	Univ.
79-89110	Female	28" – 50"	Univ.

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-82502	Enovis PROCARE SAC LUMB W/COMPRESSION ST, XS	\$36.55
79-82503	Enovis PROCARE SAC LUMB W/COMPRESSION ST,S	\$36.55
79-82505	Enovis PROCARE SAC LUMB W/COMPRESSION ST,M	\$36.55
79-82507	Enovis PROCARE SAC LUMB W/COMPRESSION ST,L	\$36.55
79-82508	Enovis PROCARE SAC LUMB W/COMPRESSION ST,XL	\$36.55
79-82509	Enovis PROCARE SAC LUMB W/COMPRESSION,ST,XXL	\$36.74
79-82509-10	Enovis PROCARE SAC LUM W/COMPRESS.STRAP XXXL	\$38.71
79-82509-11	Enovis PROCARE SAC LUM W/COMPRESS.STRAP XXXXL	\$37.32
79-89003	Enovis PROCARE 10 DBL-PULL SACRO-LUMBAR (S)	\$22.96
79-89005	Enovis PROCARE 10 DBL PULL SACRO-LUMBAR (M)	\$24.57
79-89006	Enovis PROCARE 10 DBL PULL SACRO-LUMB(MED/LG)	\$24.57
79-89007	Enovis PROCARE 10 DBL PULL SACRO-LUMBAR (L)	\$24.57
79-89008	Enovis PROCARE 10 DBL PULL SACRO-LUMBAR (XL)	\$24.57
79-89009	Enovis PROCARE 10*DBL PULL SACRO LUMBAR XXL	\$24.57
79-89009-10	Enovis PROCARE 10* DBL PULL SACRO LUM XXXL	\$29.22
79-89040	Enovis PROCARE 10 CRISS-CROSS BINDER (UNIV)	\$27.09
79-89083	Enovis PROCARE INDUSTRIAL BACK SPRT,SM	\$47.80
79-89085	Enovis PROCARE INDUSTRIAL BACK SPRT,M	\$47.80
79-89087	Enovis PROCARE INDUSTRIAL BACK SUPPORT - LG	\$47.80
79-89088	Enovis PROCARE INDUSTRIAL BACK SUPPORT - XL	\$47.80
79-89089	Enovis PROCARE INDUSTRIAL BACK SUPPORT - XXL	\$47.80
79-89143	Enovis PROCARE IND BACK SUPP W/SUSPENDERS-SM	\$27.83
79-89145	Enovis PROCARE IND BACK SUPP W/SUSPENDERS-M	\$27.83
79-89147	Enovis PROCARE IND BACK SUPP W/SUSPENDERS-LG	\$27.83
79-89148	Enovis PROCARE IND BACK SUPP W/SUSPENDERS-XL	\$27.83
79-89149	Enovis PROCARE IND BACK SUPP W/SUSPENDERS-XXL	\$27.83
79-89149L	Enovis PROCARE IND BACK SUPP W/SUSP XXXL	\$31.91
79-89149X	Enovis PROCARE IND BACK SUPP W/SUSP - XXXXL	\$29.22
79-89183	Enovis PROCARE 10 CRISS-CROSS W/SIDE PULL-SM	\$35.45
79-89185	Enovis PROCARE 10 CRISS-CROSS W/SIDE PULL-M	\$35.45
79-89187	Enovis PROCARE 10 CRISS-CROSS W/SIDE PULL-LG	\$35.45
79-89188	Enovis PROCARE 10 CRISS-CROSS W/SIDE PULL-XL	\$35.45
79-89189	Enovis PROCARE 10 CRISS-CROSS BINDER - XXL	\$35.45
79-89240	Enovis PROCARE CLINIC RET. SUP W/COMP STR-S/M	\$40.15
79-89241	Enovis PROCARE CLINIC RET.SUP W/COMP STR-L/XL	\$44.34

KNEE IMMOBILIZERS

**PROCARE®****Quick-Fit® Basic Knee Splint**

Foam laminate construction with 3" elastic contact closure straps provides compressive support and immobilization. Dual contoured, I-beam posterior stays and moveable medial/lateral double stay casings for correct fit.

Indication: Ideal for immobilization of the knee joint.

	Length	Max Thigh Circum.®	Size
79-96011	10"	29"	Univ.
79-96012	12"	29"	Univ.
79-96014	14"	29"	Univ.
79-96016	16"	29"	Univ.
79-96018	18"	36"	Univ.
79-96019	20"	36"	Univ.
79-96021	22"	36"	Univ.
79-96024	24"	36"	Univ.
79-96026	26"	36"	Univ.

®Circumference 4" above patella.

**PROCARE®****Clinic 3-Panel Splint**

Universal three-piece design for ease of adjustment and correct size adjustment. Foam posterior panel may be trimmed for a customized fit. Foam/mesh side panels with I-beam stays and encircling loop-lock straps for maximum stability. Straight stay version features popliteal pad.

Indication: Ideal for post-operative immobilization; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

	Length	Max Thigh Circum.	Size
Straight Stays			
79-80810	16"		
79-80820	20"		
79-80830	24"		



PRO CARE®

Universal Foam Knee Splint

Constructed of open cell foam and detachable medial/lateral aluminum I-beam stay panels for easy adjustment and custom fit. Three double bend posterior aluminum stays. Alternating web straps provide dual adjustability, using a D-ring buckle and a hook tab.

Indication: Ideal for post-operative immobilization, grade 1 and grade 2 collateral ligament sprain, and patella subluxations.

	Length	Size
79-80600	8"	Univ.
79-96870	10"	Univ.
79-96871	12"	Univ.
79-96914	14"	Univ.
79-96880	16"	Univ.
79-96881	18"	Univ.
79-96882	20"	Univ.
79-96890	22"	Univ.
79-96891	24"	Univ.
79-96900	26"	Univ.



PRO CARE®

Universal Foam Knee with Loop Lock Closure

Constructed of durable open cell foam with moveable medial/lateral stay casings for proper positioning. Available with contoured or straight posterior stays. Immobilizer can be trimmed for custom fit. 4" elastic patella strap for compressive support. Straight stay version features open patella cut-out and popliteal pad.

Indication: Ideal for post-operative immobilization; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

Contour Stays	Straight Stays	Length	Max Thigh Circ.	Size	Strap Color
79-80611	79-80610	12"	22.5"	Univ.	White
79-80360		14"	29"	Univ.	Black
79-80621	79-80620	16"	24"	Univ.	White
79-80370		19"	29"	Univ.	Black
79-80631	79-80630	20"	29"	Univ.	White
79-80380		22"	29"	Univ.	Black
79-80641	79-80640	24"	29"	Univ.	White
79-80390		26"	29"	Univ.	Black
79-80651	79-80650	28"	29"	Univ.	White



PRO CARE®

Compression Knee Dressing

Trimable foam laminate body with 6" wide multi-panel plush elastic straps for full compression to aid in prevention of edema. Moveable medial/lateral stay casings for proper positioning. Adjustable foam pop pad included.

Indication: Ideal for providing compression and immobilization of the knee following arthroscopic and reconstructive knee surgery; grade 1 and grade 2 collateral ligament sprain; patellar subluxations.

Contoured Stays	Straight Stays	Without Stays	Length	Max Thigh Circ.	Size
79-80090	79-80320	79-80430	16"	29"	Univ.
79-80190	79-80330	79-80440	20"	30"	Univ.
79-80260	79-80340	79-80450	24"	32"	Univ.
79-80280			28"	32"	Univ.



PROCARE®

KneeRANGER® II - Universal

Post-operative immobilization with protected range-of-motion. Dual axis, polycentric hinge allows for convenient dial-in adjustment with flexion and extension stops at 0°, 15°, 30°, 60° and 90°. Available in two (2) lengths. Trimmable foam leg wraps for customized fit. Medial/lateral plastic shells riveted to hinge bars provide additional support following post-operative procedures. Accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for providing protected range-of-motion for post-operative immobilization.



PROCARE®

KneeRANGER® II - Sized

One-piece, 24" length design allows for easy application and adjustment. Dual axis hinge provides protected, full range-of-motion control with flexion and extension stops at 0°, 15°, 30°, 60°, and 90°. Breathable nylon/fiber/flannelette construction. Integrated strap/hinge bar design accommodates Ankle/Shoe Insert for added rotational control.

Indication: Ideal for immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs.

	Circum.*	Size
79-94253	14" – 16"	S
79-94255	16" – 18"	M
79-94257	18" – 20"	L
79-94258	22" – 24"	XL

*Circumference measurement taken 4" above mid-patella.



PROCARE®

KneeRANGER® Lite

Economical, lightweight, universal leg brace for warmer climates. Narrow band circumferential cuffs replace full cut leg wraps. Integrated strap/hinge bar design. Provides immobilization and protected range-of-motion in 7.5° increments from 0° – 120° of flexion and 0° – 70° of extension.

Indication: Ideal for immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs.

Length°	Qty.	Size
79-94		Univ.
Accessories		
79-94041	Repl wraps	Set
79-96170	Foam Condyle Pads	Set

*Maximum circumference: thigh 31", calf 22".



PROCARE®

Dual Axis Hinge

3/16" of set polycentric design provides stability and full range-of-motion control with convenient dial-in flexion/extension stops at 0°, 15°, 30°, 60° and 90°. Patent No. 5, 443,444. Sold in pairs.

Indication: Ideal for casting, fabrication and range-of-motion control.

	Length
79-94400	Dual Axis Hinge 10.75"
79-94410	Dual Axis Hinge 14.5"
79-94550	w/rubber condyle pad 10.75"
79-94560	w/rubber condyle pad and cap cover 10.75"
79-94580	w/rubber condyle pad 14.5"
79-94590	w/rubber condyle pad and cap cover 14.5"
79-93820	w/rubber condyle pad and cap cover 21.5"

THIGH, KNEE & CALF SUPPORTS

PROCARE®**Sport Knee Sleeve**

1/8" neoprene sleeve designed for durability and greater comfort. Extra length pad at the knee of offers localized heat retention, compressive support, and improved protection. Choice between open or closed patella. 13" Length.

Indication: Ideal for general compression and support.

Closed Patella	Open Patella	Open Pop.	Circum.*	Size
79-82002	79-82012	79-82122	13.5" – 15.5"	XS
79-82003	79-82013	79-82123	15.5" – 18"	S
79-82005	79-82015	79-82125	18" – 20.5"	M
79-82007	79-82017	79-82127	20.5" – 23"	L
79-82008	79-82018	79-82128	23" – 25.5"	XL
79-82009	79-82019	79-82129	25.5" – 28"	2XL
79-82009-10	79-82019-10	79-82129-10	28" – 31"	3XL
79-82019-11	79-82129-11		31" – 34"	4XL

*Measurement taken 6" above mid-patella.





PROCARE®

Clinic Knee Sleeve

1/8" single-sided nylon neoprene with contoured design for controlled, uniform compression. 10" length.

Indication: Ideal for general compression and support.

Open Patella	Circum.	Size
79-82612	13.5" – 15.5"	XS
79-82613	15.5" – 18"	S
79-82615	18" – 20.5"	M
79-82617	20.5" – 23"	L
79-82618	23" – 25.5"	XL
79-82619	25.5" – 28"	2XL



PROCARE®

Knee Support with Reinforced Patella

1/8" neoprene provides compressive support and warmth. Reversible. Blue/black neoprene. 10" length.

Indication: Ideal for Patellofemoral (knee cap) disorders.

	Circum.*	Size
79-82633	15.5" – 18"	S
79-82635	18" – 20.5"	M
79-82637	20.5" – 23"	L
79-82638	23" – 25.5"	XL
79-82639	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®


Thigh Sleeve

1/8" neoprene thigh sleeve provides firm compression and soothing warmth for hamstring and quadriceps strain. 12" length.

Indication: Ideal for strains or contusions to the quadriceps, hamstring, or groin muscles.

	Circum.*	Size
79-82332	18" – 19"	XS
79-82333	19" – 20"	S
79-82335	20" – 22"	M
79-82337	22" – 25"	L
79-82338	25" – 28"	XL
79-82339	28" – 31"	2XL

*Circumference measurement taken at mid-thigh.



PROCARE®

Calf Sleeve

1/8" neoprene sleeve contours to assure proper fit and correct degree of compression and support for strains and sprains. 12" Length.

Indication: Ideal for strains; shin splints.

	Circum.*	Size
79-82342	12" – 13"	XS
79-82343	13" – 14"	S
79-82345	14" – 15"	M
79-82347	15" – 17"	L
79-82348	17" – 19"	XL

*Circumference measurement taken at mid-calf.

PATELLOREMORAL

PROCARE®**Stabilized Knee Support**

1/8" neoprene with medial/lateral spiral stays for added support. Superior and inferior compression straps help prevent migration and provide additional support. Universal neoprene buttress may be adjusted to allow for proper patella management and stabilization.

Indication: Ideal for mild patella tendonitis, chondromalacia or subluxations.

Closed Pop.	Open Pop.	Circum.*	Size
79-82722	79-82752	13.5" – 15.5"	XS
79-82723	79-82753	15.5" – 18"	S
79-82725	79-82755	18" – 20.5"	M
79-82727	79-82757	20.5" – 23"	L
79-82728	79-82758	23" – 25.5"	XL
79-82729	79-82759	25.5" – 28"	2XL
79-82729-10		28" – 31"	3XL
79-82729-11		31" – 34"	4XL
79-82491	Horseshoe Buttress		Univ.

*Measurement taken 6" above mid-patella.





PROCARE®

Patella Stabilizer with Buttress

1/8" neoprene support with choice of adjustable felt horseshoe-buttress or sewn-in donut buttress to accommodate mild patella tendonitis, chondromalacia, or subluxation. Donut buttress may be trimmed to allow for medial, lateral, inferior, or superior buttress. Medial/lateral spiral stays and superior/inferior compression straps provide additional stability. Open popliteal available with donut buttress only.

Indication: For Patello Femoral (knee cap) disorders.

Horseshoe Buttress	Done Buttress w/Open Pop.	Circum.*	Size
79-92852	79-92842	13.5" – 15.5"	XS
79-92853	79-92843	15.5" – 18"	S
79-92855	79-92845	18" – 20.5"	M
79-92857	79-92847	20.5" – 23"	L
79-92858	79-92848	23" – 25.5"	XL
79-92859	79-92849	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Sport Knee Reinforced

1/8" neoprene sleeve with universal neoprene buttress offers full patellar control for chondromalacia, tendonitis, and patella tracking abnormalities. Neoprene buttress may be adjusted to allow for medial, lateral, inferior or superior patella stabilization. Choice between closed or open popliteal. Anterior oval compression pad with patella cutout provides additional warmth and compression.

Indication: Ideal for patella femoral tracking dysfunction; chondromalacia.

Closed Pop.	Open Pop.	Circum.*	Size
79-82702	79-82052	13.5" – 15.5"	XS
79-82703	79-82053	15.5" – 18"	S
79-82705	79-82055	18" – 20.5"	M
79-82707	79-82057	20.5" – 23"	L
79-82708	79-82058	23" – 25.5"	XL
79-82709	79-82059	25.5" – 28"	2XL
79-82709-10		28" – 31"	3XL
79-82491	Horseshoe Buttress		Univ.

*Measurement taken 6" above mid-patella.



PROCARE®

Reinforced Patella Stabilizer

1/8" neoprene with sewn-in horseshoe tubular buttress helps stabilize and improve patella tracking. Adjustable superior/inferior straps with contact closure for additional compression and stability.

Indication: Ideal for patella femoral tracking dysfunction; chondromalacia, patella tendonitis.

	Circum.*	Size
79-94433	15.5" – 18"	S
79-94435	18" – 20.5"	M
79-94437	20.5" – 23"	L
79-94438	23" – 25.5"	XL
79-94439	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Hinged Patella Stabilizer

1/8" neoprene construction with sewn-in tubular patella buttress and removable dual-axis polycentric hinges for additional medial/lateral support. Adjustable superior/inferior straps for compression and stability.

Indication: Ideal for medial/lateral instabilities with patella tendonitis, chondromalacia, subluxation, grade 1 collateral ligament sprains.

	Circum.*	Size
79-94483	15.5" – 18"	S
79-94485	18" – 20.5"	M
79-94487	20.5" – 23"	L
79-94488	23" – 25.5"	XL
79-94489	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PROCARE®

Lateral Patella Stabilizer

1/8" neoprene with tubular sewn-in "J" buttress and dynamic side-pull compression strap supports patella malalignment, dislocation, and subluxation. Sold in right and left configuration.

Indication: Ideal for supporting lateral patella subluxation, dislocation or lateral Patellofemoral malalignment.

Right	Left	Circum.*	Size
79-94463	79-94473	15.5" – 18"	S
79-94465	79-94475	18" – 20.5"	M
79-94467	79-94477	20.5" – 23"	L
79-94468	79-94478	23" – 25.5"	XL
79-94469	79-94479	25.5" – 28"	2XL
79-94469-10	79-94479-10	28" – 31"	3XL
79-94469-11	79-94479-11	31" – 34"	4XL

**Measurement taken 6" above mid-patella.*



PROCARE®

Universal Patella Knee with Buttress

1/8" neoprene wrap with open patella and spiral stays. Sleeve can accommodate buttress pad for additional patellar support.

Indication: Ideal for mild support and compression; mild patella bursitis.

	Circum.*	Size
79-82460	≤ 21"	Univ.
79-82461	≤ 25"	Univ./XL

**Circumference of knee fully extended.*



PROCARE®

Lites Visco Knee Support

Constructed of lightweight four-way stretch elastic, and two medial/lateral stays for compression and support. Silicone buttress provides patella stabilization.

Indication: Ideal for chondromalacia; osteoarthritis; patellofemoral tracking dysfunction; grade I collateral ligament sprain.

	Circum.*	Size
79-80162	14.5" – 16"	XS
79-80163	16" – 17"	S
79-80165	17" – 18"	M
79-80167	18" – 19.25"	L
79-80168	19.25" – 20.5"	XL
79-80169	20.5" – 22"	2XL
79-80169-10	22" – 23"	3XL

**Measurement taken 6" above mid-patella.*



PROCARE®

Surround® Patella Strap

The primary body of the strap is a soft, latex free nylon/foam laminate with contact closures. A unique cylindrical FLOAM™ bladder is encased inside the body of the primary patella strap and may be adjusted for added compression. Floam provides uniform compression and a customized fit while prohibiting pressure from being displaced to bony prominences. Patent pending.

Indication: Helps relieve anterior knee pain caused by patellar tendonitis, chondromalacia and Osgood Schlatter by applying pressure to the patella tendon.

	Circum. Just Below Kneecap	Size
79-80322	≤10"	XS
79-80323	10.5" – 12.5"	S
79-80325	13" – 15"	M
79-80327	15.5" – 17.5"	L
79-80328	18"+	XL

HINGED KNEE BRACES



Reaction Knee® Brace

The Reaction Web™ knee brace is a responsive, webbed approach to anterior knee pain that gives a distinct alternative to the basic knee sleeve. The elastomeric web is an innovative solution that provides energy dispersion to the painful knee, allowing comfortable motion in a brace that actually stays in place.

Indication: The elastomeric web design helps to reduce pain by dispersing energy across the knee. The web absorbs shock and shifts the peak loads away from the painful area of the knee. The elastomeric web acts to dynamically stabilize the patella on all sides, bringing the patella into proper tracking position to reduce pain for patients with general patellofemoral instabilities.

	Measurement*	Size
11-0215-2-XXXXX	13"-18.5" (33-47 cm)	XS/S
11-0215-3-XXXXX	18.5"-23.5" (47-60 cm)	M/L
11-0215-4-XXXXX	23.5"-29.5" (60-75 cm)	XL/XXL
11-0215-5-XXXXX	29.5"-31" (75-79 cm)	XXXL
Replacement Under Sleeve		
11-0541-2	13" -18.5" (33-47 cm)	XS/S
11-0541-3	18.5"-21" (47-53.25 cm)	M
11-0541-4	21"-23.5" (53.25-59.5 cm)	L
11-0541-5	23.5"-26.5" (59.5-67.25 cm)	XL
11-0541-6	26.5"-29.5" (67.25-75 cm)	XXL
11-0541-7	29.5"-32" (75-81.25 cm)	XXXL
Color Information = XXXXX		
Orange 14000	Red 01000	Green 11000
Pink 09000	Yellow 08000	Blue 02000
Grey	n/a (If color is not specified web will be gray)	

*Measurements listed below are the circumference taken around the thigh, 6" above mid-patella.





PRO CARE®

Hinged Knee Support

Compressive neoprene knee support is available with and without buttress. Universal buttress may be adjusted to allow for medial, lateral, inferior, or superior patella stabilization. Sturdy removable dual axis hinges and superior/inferior contact closure straps provide medial/lateral ligament support.

Indication: Ideal for moderate medial/lateral knee instability, mild ACL/PCL sprain, osteoarthritis, patellofemoral dysfunction.

3/16" w/Open Pop.	1/8" w/Reinf. Pad	Circum.*	Size
79-82152	79-82162	13.5" – 15.5"	XS
79-82153	79-82163	15.5" – 18"	S
79-82155	79-82165	18" – 20.5"	M
79-82157	79-82167	20.5" – 23"	L
79-82158	79-82168	23" – 25.5"	XL
79-82159	79-82169	25.5" – 28"	2XL
79-82159-10	79-82739-10	28" – 30.5"	3XL
	79-82739-11	30.5" – 33"	4XL

*Measurement taken 6" above mid-patella.

1/8" w/Univ. Butt.	1/8" w/Open Pop.	Circum.*	Size
79-82732	79-82742	13.5" – 15.5"	XS
79-82733	79-82743	15.5" – 18"	S
79-82735	79-82745	18" – 20.5"	M
79-82737	79-82747	20.5" – 23"	L
79-82739	79-82749	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.



PRO CARE®

WeekENDER™ Recreational Activity Brace

Lightweight, rigid, low profile. May be used for post-op, rehabilitation or functional application. Dual axis motion control hinges with foam/chamois condyle pads provide full and limited range-of-motion with flexion/extension stops at 0°, 15°, 30°, 60° and 90°.

Indication: Recommended for mild medial/lateral instability, including ACL/PCL laxity.

Right	Left	Circum.*	Size
79-94353	79-94363	14.5" – 18.5"	S
79-94355	79-94365	18.5" – 21"	M
79-94357	79-94367	21" – 23.5"	L
79-94358	79-94368	23.5" – 26.5"	XL
79-94359	79-94369	26.5" – 29.5"	2XL



PRO CARE®

Reddie® Brace

1/8" neoprene wraparound configuration and strap tabs provide ease of application. Posterior strap adjustment allows for one-time adjustment and proper positioning of hinges. Removable dual axis polycentric hinges provide medial/lateral support. Available with open popliteal only.

Indication: Ideal for mild MCL and/or LCL sprains, mild arthritis, post arthroscopic meniscectomy, or mild global instabilities.

	Circum.*	Size
79-82392	13.5" – 15.5"	XS
79-82393	15.5" – 18"	S
79-82395	18" – 20.5"	M
79-82397	20.5" – 23"	L
79-82398	23" – 25.5"	XL
79-82399	25.5" – 28"	2XL
79-82399-10	28" – 30.5"	3XL

*Measurement taken 6" above mid-patella.



PRO CARE®

ACTION™ Neoprene Brace/Wrap - 1/8"

Available in black 1/8" neoprene only. Front opening for easy on/off and custom fitting. Dual axis motion control hinges provide convenient dial-in adjustments with flexion/extension stops at 0°, 15°, 30°, 60° and 90°. Fully adjustable superior/inferior contact closure straps. Available in 13" or 17" length.

Indication: Ideal for mild to moderate ACL/PCL sprains; osteoarthritis, collateral ligament sprains and strains.

13"	17"	Circum.*	Size
79-94402	79-94412	13.5" – 15.5"	XS
79-94403	79-94413	15.5" – 18"	S
79-94405	79-94415	18" – 20.5"	M
79-94407	79-94417	20.5" – 23"	L
79-94408	79-94418	23" – 25.5"	XL
79-94409	79-94419	25.5" – 28"	2XL

*Measurement taken 6" above mid-patella.

Replacement Strap with Standard Foam Pads

79-94373	Thigh, Calf & Condyle Pads	S
79-94375	Thigh, Calf & Condyle Pads	M
79-94377	Thigh, Calf & Condyle Pads	L
79-94378	Thigh, Calf & Condyle Pads	XL
79-94379	Thigh, Calf & Condyle Pads	2XL

Chamois Replacement Pad Set

79-94323	Thigh, Calf & Condyle Pads	S
79-94325	Thigh, Calf & Condyle Pads	M
79-94327	Thigh, Calf & Condyle Pads	L
79-94328	Thigh, Calf & Condyle Pads	XL
79-94329	Thigh, Calf & Condyle Pads	2XL

*Measurement taken 6" above mid-patella.

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
81-02JL	Enovis AIRCAST AIR-STIRRUP ANKLE BRCE,PED,LT	\$50.72
81-02JR	Enovis AIRCAST AIR-STIRRUP ANKLE BRCE,PED,RT	\$50.72
79-80090	Enovis PROCARE 16 CKD KNEE IMMOBILIZER	\$48.43
79-80162	Enovis PROCARE LITES VISCO KNEE SUPPORT, XS	\$45.84
79-80163	Enovis PROCARE LITES VISCO KNEE SUPPORT, S	\$45.84
79-80165	Enovis PROCARE LITES VISCO KNEE SUPPORT, M	\$45.84
79-80167	Enovis PROCARE LITES VISCO KNEE SUPPORT, L	\$45.84
79-80168	Enovis PROCARE LITES VISCO KNEE SUPPORT, XL	\$45.84
79-80169	Enovis PROCARE LITES VISCO KNEE SUPPORT, XXL	\$45.84
79-80169-10	Enovis PROCARE LITES VISCO KNEE SUPPORT, XXXL	\$45.84
79-80190	Enovis PROCARE 20 CKD KNEE IMMOBILIZER	\$57.71
79-80260	Enovis PROCARE 24 CKD KNEE IMMOBILIZER	\$65.62
79-80280	Enovis PROCARE 28" CKD KNEE IMMOBILIZER	\$100.03
79-80322	Enovis PROCARE SURROUND PATELLA STRAP,XS	\$18.66
79-80323	Enovis PROCARE SURROUND PATELLA STRAP,S	\$18.66
79-80325	Enovis PROCARE SURROUND PATELLA STRAP,M	\$18.66
79-80327	Enovis PROCARE SURROUND PATELLA STRAP,L	\$18.66
79-80328	Enovis PROCARE SURROUND PATELLA STRAP,XL	\$18.66
79-80410	Enovis PROCARE 12"CLINIC 3-PANEL KNEE SPLIN,U	\$28.28
79-80420	Enovis PROCARE 16"CLINIC 3-PANEL KNEE SPL,UNI	\$26.33
79-80600	Enovis PROCARE 8" KNEE IMMOBILIZER	\$19.68
79-80610	Enovis PROCARE 12" UNIV KNEE IMMOB.W/ ST.STAY	\$53.52
79-80611	Enovis PROCARE 12" UNIV KNEE IMMOB.W/ CTR.ST.	\$54.05
79-80620	Enovis PROCARE 16" UNIV KNEE IMMOB.W/ ST.STAY	\$60.52
79-80621	Enovis PROCARE 16" UNIV KNEE IMMOB.W/ CTR.ST.	\$58.94
79-80630	Enovis PROCARE 20" UNIV KNEE IMMOB.W/ ST.STAY	\$66.68
79-80631	Enovis PROCARE 20" UNIV KNEE IMMOB.W/ CTR.ST.	\$61.03
79-80640	Enovis PROCARE 24" UNIV KNEE IMMOB.W/ ST.STAY	\$73.93
79-80641	Enovis PROCARE 24" UNIV KNEE IMMOB.W/ CTR.ST.	\$68.87
79-80650	Enovis PROCARE 28" UNIV KNEE IMMOB.W/ ST.STAY	\$70.63
79-80651	Enovis PROCARE 28" UNIV KNEE IMMOB.W/ CTR.ST.	\$71.57
79-80810	Enovis PROCARE 16" CLINIC 3-PANEL KNEE STR ST	\$32.42
79-80820	Enovis PROCARE 20" CLINIC 3-PANEL KNEE STR ST	\$35.15
79-80830	Enovis PROCARE 24" CLINIC 3-PANEL KNEE STR ST	\$39.65
79-82002	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,XS	\$31.68
79-82003	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,S	\$31.68
79-82005	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,M	\$31.68
79-82007	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,L	\$31.68

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-82007	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,L	\$31.68
79-82008	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,XL	\$31.68
79-82009	Enovis PROCARE SPORT KNEE SLV CLOSED PAT,XXL	\$31.68
79-82009-10	Enovis PROCARE SPORT KNEE SLV CLOSED PAT, 3XL	\$31.68
79-82019-11	Enovis PROCARE SPORT KNEE SLV OPEN PAT, XXXXL	\$31.18
79-82152	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,XS	\$44.17
79-82153	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,S	\$44.17
79-82155	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,M	\$44.17
79-82157	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,L	\$44.17
79-82158	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,XL	\$44.17
79-82159	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,XXL	\$44.17
79-82159-10	Enovis PROCARE HINGD KN SPPRT,W/OPEN POP,XXXL	\$44.17
79-82332	Enovis PROCARE THIGH SLEEVE,XS	\$15.72
79-82333	Enovis PROCARE THIGH SLEEVE,S	\$15.72
79-82335	Enovis PROCARE THIGH SLEEVE,M	\$15.72
79-82337	Enovis PROCARE THIGH SLEEVE,L	\$15.72
79-82338	Enovis PROCARE THIGH SLEEVE,XL	\$15.72
79-82339	Enovis PROCARE THIGH SLEEVE,XXL	\$16.00
79-82342	Enovis PROCARE CALF SLEEVE,XS	\$12.25
79-82343	Enovis PROCARE CALF SLEEVE,S	\$12.25
79-82345	Enovis PROCARE CALF SLEEVE,M	\$12.25
79-82347	Enovis PROCARE CALF SLEEVE,L	\$12.25
79-82348	Enovis PROCARE CALF SLEEVE,XL	\$12.25
79-82392	Enovis PROCARE REDDIE BRACE,XS	\$48.29
79-82393	Enovis PROCARE REDDIE BRACE,S	\$48.29
79-82395	Enovis PROCARE REDDIE BRACE,M	\$48.29
79-82397	Enovis PROCARE REDDIE BRACE,L	\$48.29
79-82398	Enovis PROCARE REDDIE BRACE,XL	\$48.29
79-82399	Enovis PROCARE REDDIE BRACE,XXL	\$48.29
79-82399-10	Enovis PROCARE REDDIE BRACE,XXXL	\$48.29
79-82399-11	Enovis PROCARE REDDIE BRACE, XXXXL	\$48.29
79-82460	Enovis PROCARE PATELLA KNEE WRAP II, UNIV	\$28.95
79-82461	Enovis PROCARE PATELLA KNEE WRAP II, XL	\$66.78
79-82612	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,XS	\$9.98
79-82613	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,S	\$9.98
79-82615	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,M	\$9.98
79-82617	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,L	\$9.98
79-82618	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,XL	\$9.98

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-82619	Enovis PROCARE CLINIC KNEE SLEEVE,OPN PAT,XXL	\$11.83
79-82633	Enovis PROCARE KNEE SUPPORT W/PATELLA,S	\$14.71
79-82635	Enovis PROCARE KNEE SUPPORT W/PATELLA,M	\$14.71
79-82637	Enovis PROCARE KNEE SUPPORT W/PATELLA,L	\$14.71
79-82638	Enovis PROCARE KNEE SUPPORT W/PATELLA,XL	\$14.71
79-82639	Enovis PROCARE KNEE SUPPORT,W/PATELLA,XXL	\$14.76
79-82702	Enovis PROCARE SPORT KNEE,UNIV BUTT,XS	\$36.43
79-82703	Enovis PROCARE SPORT KNEE,UNIV BUTT,S	\$36.43
79-82705	Enovis PROCARE SPORT KNEE,UNIV BUTT,M	\$36.43
79-82707	Enovis PROCARE SPORT KNEE,UNIV BUTT,L	\$36.43
79-82708	Enovis PROCARE SPORT KNEE,UNIV BUTT,XL	\$36.43
79-82709	Enovis PROCARE SPORT KNEE,UNIV BUTT,XXL	\$36.43
79-82709-10	Enovis PROCARE SPORT KNEE,UNIV BUTT, XXXL	\$36.43
79-82722	Enovis PROCARE STABILIZED KN W/UNIV BUTT,XS	\$34.16
79-82723	Enovis PROCARE STABILIZED KN W/UNIV BUTT,S	\$34.16
79-82725	Enovis PROCARE STABILIZED KN W/UNIV BUTT,M	\$34.16
79-82727	Enovis PROCARE STABILIZED KN W/UNIV BUTT,L	\$34.16
79-82728	Enovis PROCARE STABILIZED KN W/UNIV BUTT,XL	\$34.16
79-82729	Enovis PROCARE STABILIZED KN W/UNIV BUTT,XXL	\$34.16
79-82729-10	Enovis PROCARE STABILIZED KN W/UNIV BUTT, 3XL	\$34.86
79-82729-11	Enovis PROCARE STABILIZED KN W/UNIV BUTT, 4XL	\$34.86
79-92852	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,XS	\$37.69
79-92853	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,S	\$37.69
79-92855	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,M	\$37.69
79-92857	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,L	\$37.69
79-92858	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,XL	\$37.69
79-92859	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT,XXL	\$37.69
79-92859-10	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT 3XL	\$43.61
79-92859-11	Enovis PROCARE PATELLA STAB W/HRSSHoe BTT 4XL	\$43.61
79-94040	Enovis PROCARE KNEE RANGER - LITE	\$102.21
79-94230	Enovis PROCARE KNEE RANGER II SHORT, UNIV	\$72.97
79-94240	Enovis PROCARE KNEE RANGER II LONG,UNIV	\$107.51
79-94253	Enovis PROCARE KNEE RANGER II LONG, SMALL	\$104.52
79-94255	Enovis PROCARE KNEE RANGER II LONG, MEDIUM	\$104.52
79-94257	Enovis PROCARE KNEE RANGER II LONG, LARGE	\$104.52
79-94258	Enovis PROCARE KNEE RANGER II LONG, X-LARGE	\$104.52
79-94353	Enovis PROCARE WEEKENDER KNEE BRACE RT/SM	\$176.03
79-94355	Enovis PROCARE WEEKENDER KNEE BRACE RT/MED	\$176.03

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-94358	Enovis PROCARE WEEKENDER KNEE BRACE RT/XL	\$176.03
79-94359	Enovis PROCARE WEEKENDER KNEE BRACE RT/XXL	\$176.03
79-94363	Enovis PROCARE WEEKENDER KNEE BRACE LT/SM	\$176.03
79-94365	Enovis PROCARE WEEKENDER KNEE BRACE LT/MED	\$176.03
79-94367	Enovis PROCARE WEEKENDER KNEE BRACE LT/LG	\$176.03
79-94368	Enovis PROCARE WEEKENDER KNEE BRACE LT/XL	\$176.03
79-94369	Enovis PROCARE WEEKENDER KNEE BRACE LT/XXL	\$176.03
79-94373	Enovis PROCARE WEEKENDER STRAP SET - S	\$19.45
79-94375	Enovis PROCARE WEEKENDER STRAP SET - M	\$19.45
79-94377	Enovis PROCARE WEEKENDER STRAP SET - L	\$19.45
79-94378	Enovis PROCARE WEEKENDER STRAP SET - XL	\$19.45
79-94379	Enovis PROCARE WEEKENDER STRAP SET - XXL	\$19.45
79-94400	Enovis PROCARE HINGE,DUAL AXIS,3/16 OFFSET,PR	\$38.78
79-94402	Enovis PROCARE ACTION HINGED KNEE WRAP/SHT XS	\$71.48
79-94403	Enovis PROCARE ACTION HNGD KNE WRAP,SHT,S	\$71.48
79-94405	Enovis PROCARE ACTION HNGD KNE WRAP,SHT,M	\$71.48
79-94407	Enovis PROCARE ACTION HNGD KNE WRAP,SHT,L	\$71.48
79-94408	Enovis PROCARE ACTION HNGD KNE WRAP,SHT,XL	\$71.48
79-94409	Enovis PROCARE ACTION HNGD KNE WRAP,SH,XXL	\$71.48
79-94410	Enovis PROCARE DUAL AXIS HINGE 3/16 OFFSET,PR	\$38.78
79-94433	Enovis PROCARE PATELLA STABILIZER, REINF.S	\$65.13
79-94435	Enovis PROCARE PATELLA STABILIZER, REINF.M	\$35.95
79-94437	Enovis PROCARE PATELLA STABILIZER, REINF.L	\$35.95
79-94438	Enovis PROCARE PATELLA STABILIZER, REINF.XL	\$65.13
79-94439	Enovis PROCARE PATELLA STABILIZER, REINF. XXL	\$74.23
79-94463	Enovis PROCARE PATELLA STAB, LATERAL, RT, S	\$73.23
79-94465	Enovis PROCARE PATELLA STAB,LATERAL,RT,M	\$73.23
79-94467	Enovis PROCARE PATELLA STAB,LATERAL,RT,L	\$73.23
79-94468	Enovis PROCARE PATELLA STAB,LATERAL,RT,XL	\$73.23
79-94469	Enovis PROCARE PATELLA STAB,LATERAL,RT,XXL	\$73.23
79-94469-10	Enovis PROCARE PATELLA STAB,LATERAL,RT,XXXL	\$73.23
79-94469-11	Enovis PROCARE PATELLA STAB,LATERAL,RT,XXXXL	\$73.23
79-94483	Enovis PROCARE PATELLA STABILIZER, HINGED S	\$103.58
79-94485	Enovis PROCARE PATELLA STABILIZER, HINGED M	\$103.58
79-94487	Enovis PROCARE PATELLA STABILIZER, HINGED L	\$103.58
79-94488	Enovis PROCARE PATELLA STABILIZER, HINGED XL	\$103.58
79-94489	Enovis PROCARE PATELLA STABILIZER, HINGED XXL	\$103.58
79-96011	Enovis PROCARE 10" UNIV BASIC KNEE SPLINT	\$20.12

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-96012	Enovis PROCARE 12" UNIV BASIC KNEE SPLINT	\$23.00
79-96014	Enovis PROCARE 14" UNIV BASIC KNEE SPLINT	\$23.49
79-96016	Enovis PROCARE 16" UNIV BASIC KNEE SPLINT	\$23.79
79-96016-004	Enovis PROCARE 16" UNIV BASIC KNEE SPLINT	\$33.01
79-96018	Enovis PROCARE 18" UNIV BASIC KNEE SPLINT	\$24.60
79-96019	Enovis PROCARE 20" UNIV BASIC KNEE SPLINT	\$25.47
79-96019-004	Enovis PROCARE 20" UNIV BASIC KNEE SPLINT	\$33.01
79-96021	Enovis PROCARE 22" UNIV BASIC KNEE SPLINT	\$27.71
79-96024	Enovis PROCARE 24" UNIV BASIC KNEE SPLINT	\$29.36
79-96024-004	Enovis PROCARE 24" UNIV BASIC KNEE SPLINT	\$33.01
79-96026	Enovis PROCARE 26" UNIV BASIC KNEE SPLINT	\$34.27
79-96870	Enovis PROCARE 10" UNIV FOAM KNEE	\$35.25
79-96871	Enovis PROCARE 12" UNIV FOAM KNEE	\$36.09
79-96880	Enovis PROCARE 16" UNIV FOAM KNEE	\$44.04
79-96881	Enovis PROCARE 18" UNIV FOAM KNEE	\$47.88
79-96882	Enovis PROCARE 20" UNIV FOAM KNEE	\$49.46
79-96890	Enovis PROCARE 22" UNIV FOAM KNEE	\$52.42
79-96891	Enovis PROCARE 24" UNIV FOAM KNEE	\$54.40
79-96900	Enovis PROCARE 26" UNIV FOAM KNEE	\$62.15
79-96901	Enovis PROCARE TOE CAP PROTECTOR, 12 PK	\$35.13
79-96914	Enovis PROCARE 14" UNIVERSAL FOAM KNEE	\$39.06

ANKLE BRACES & SUPPORT – ANKLE STIRRUPS

AIRCAST®

Ankle Stirrups

Aircast pioneered functional management of ankle sprains with the Air-Stirrup® Ankle Brace. Setting the standard of care for ankle sprains, the Air-Stirrup's design features the patented Duplex™ aircell technology that enhances circulation and reduces swelling.

Air-Stirrup® Ankle Brace

Anatomically designed semi-rigid shells lined with the patented Duplex™ aircell system for support and graduated compression during ambulation. The compression promotes efficient edema reduction in addition to helping accelerate rehabilitation from ankle sprains.

Indication: Acute injury, Post op; Ankle sprains grade I, II, III; Chronic instability.

Left	Right	Height	Size
02JL	02JR	6"	Ped.
02CL	02CR	8.75"	S
02BL	02BR	9"	M
02AL	02AR	10.5"	L
<i>Maximum Ankle Circumference 16.5"</i>			
0213	Heel assembly replacement		
0130	Replacement inflation tubing		





PROCARE®

Surround® FLOAM™ Ankle

Provides stability to the ankle by minimizing ankle rotation. Comprised of two polypropylene outer shells connected with soft, adjustable foot piece. A unique FLOAM™ bladder provides uniform compression and a customized fit. FLOAM™ has no memory, and therefore prohibits pressure from being displaced to bony prominences. Fits right or left ankle.

Indication: Acute ankle sprains; chronic ankle instabilities; osteoarthritis.

	Height	Size
79-81193	6"	Ped.
79-81195	9"	M
79-81197	10"	L



PROCARE®

Surround® Gel Ankle

Rigid thermoplastic shells with adjustable heel strap to help prevent inversion/eversion. Air and gel bladder may be used for cold therapy to treat tender or swollen ankles. Fits right or left ankle.

Indication: Acute injury; Post-op; Ankle sprains grade I, II, III; Chronic instability.

	Height	Size
79-97863	8.5"	S
79-97865	9"	M
79-97867	10"	L



PRO CARE®

Excelerator™ Stirrup Ankle Splint

Constructed of a three-layer foam system for added comfort and bi-lateral rigid shells to help prevent inversion and eversion. Universal sizing fits right or left ankle. Contact closure straps for ease of application and removal.

Indication: Acute ankle sprain; chronic ankle instability; subtalar joint pain, osteoarthritis.

	Height	Size
79-81323	8"	S
79-81327	9.5"	L



PRO CARE®

Stirrup Ankle Support

Contoured, high impact molded shells help prevent inversion/eversion and provide protection for acute ankle sprains and chronic instabilities. Removable, padded nylon/neoprene liner is intended to provide superior fit and comfort around ankle joint. Contoured shells for left and right ankle.

Indication: Acute ankle sprains; ankle instability.

Right	Left	Height	Size
79-81250	79-81260	10"	Univ.
79-81251	79-81261	8.5"	Univ.



PRO CARE®

Universal Ankle Brace

Rigid thermoplastic one-piece outer shell helps prevent inversion/eversion. Cushioned liner incorporates neoprene forefoot closure. Elastic proximal strap provides additional compression.

Indication: Acute and chronic ankle instabilities, as well as, prophylactic use.

	Height	Size
79-81330	10"	Univ.



AIRCAST®

Leg Brace

Designed to provide functional management of stress fractures and the graduated management of certain stable fractures of the lower leg. Features anatomically designed shells lined with the patented Duplex aircell system. The Leg Brace can be ordered alone or with an optional Anterior Panel for additional tibial protection. A sock is included with each brace.

Indication: Stress fracture; Stable fracture (graduated fracture management).

Left	Right	Height	Size
03CL	03CR	13"	S
03DL	03DR	13"	S*
03AL	03AR	15.5"	Standard
03BL	03BR	15.5"	Standard*

*With Anterior Panel.

PLANTAR FASCIITIS

**PROCARE®****Plantar Fasciitis Night Splint**

Allows for controlled application of a prolonged stretch on the Plantar Fascia and Achilles Tendon. Lightweight, sturdy, one-piece flexible shell. Adjustable bilateral dorsi-flexion straps. Cool terry liner and padded calf and ankle straps for additional comfort. Shipped with toe wedge for additional stretch on plantar fascia. Fits right or left.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
79-97753	≤ 6	≤ 7	S
79-97755	6.5 – 8.5	7.5 – 9.5	M
79-97757	9 – 11.5	10 – 12.5	L
79-97758	12+	13+	XL

**PROCARE®****Dorsiwedge™ Night Splint**

Soft padded, contoured posterior shell provides passive dorsi-flexion to address Plantar Fasciitis and Achilles Tendonitis. Removable liner may be washed for patient hygiene. Fits right or left foot. Shipped with full-length toe wedge for additional stretch on plantar fascia.

Indication: Plantar Fasciitis.

	Shoe Size Male	Shoe Size Female	Size
79-81403	≤ 6	≤ 6.5	S
79-81405	6.5 – 9.5	7 – 10	M
79-81407	10 – 12	10.5 – 12.5	L
79-81408	12+	13+	XL

POST-OP AND CAST SHOES

PROCARE®

Squared Toe Post-Op Shoe

The square toe design acts as a bumper and provides additional room and comfort for the patient. Features high ankle strapping that maintains foot in proper location and a rigid rocker sole for ease of ambulation and pressure reduction. Design allows dressing to be kept clean and toes dry, providing more protection following post-operative procedures and forefoot trauma.

Indication: Post-op; forefoot trauma.

	Shoe Size Male	Shoe Size Female	Size
79-81232	3 – 5	4 – 6	XS
79-81233	5.5 – 7	6.5 – 8	S
79-81235	7.5 – 9	8.5 – 10	M
79-81237	9.5 – 12	10.5 – 13	L
79-81238	112.5+	13+	XL




PROCARE®

IMPAX Grid Insole

Designed for use with the ProCare Squared Toe Post-Op shoe, the dual-laminate insole is designed with Impax Grid technology to redistribute load and relieve pressure. Includes one Impax insole.

Indication: Redistribute load; relieve pressure.

	Shoe Size Male	Shoe Size Female	Size
79-81332	3 – 5	4 – 6	XS
79-81333	5.5 – 7	6.5 – 8	S
79-81335	7.5 – 9	8.5 – 10	M
79-81337	9.5 – 12	10.5 – 13	L
79-81338	12.5+	13.5+	XL





PROCARE®

Post-Op Shoe

Rigid E.V.A. rocker sole and contoured heel for added patient comfort. Low-profile foam/nylon mesh upper with tricot lining. Protective tongue provides maximum comfort. Available in loop-lock design.

Indication: Post foot surgery; metatarsal injuries; soft tissue procedures.

	Shoe Size	Size
Male Loop Lock Closure		
79-90182	5 – 7	XS
79-90183	7 – 9	S
79-90185	9 – 11	M
79-90187	11 – 13	L
79-90188	13+	XL
Female Loop Lock Closure		
79-90192	2 – 4	XS
79-90193	4 – 6	S
79-90195	6 – 8	M
79-90197	8 – 10	L



PROCARE®

Off-Loading Diabetic Shoe

Specifically designed to relieve pressure on ulcerated areas of the foot during the treatment of diabetic ulcers. Squared toe box and wide foot bed leave ample room for dressing. Removable reinforced toe cover protects from external aggression and allows easy access to the treatment area. Seamless design helps eliminating potential for skin breakdown. Impax Grid tri-laminate insole is designed to help redistribute load away from ulcers by simply removing Impax foam blocks. Fits Right or Left.

Indication: Pressure relief during treatment of diabetic ulcers.

	Shoe Size Male	Shoe Size Female	Size
79-81513	4 – 6	6.5 – 8.5	S
79-81515	6.5 – 8.5	9 – 11	M
79-81517	9 – 11	11.5 – 12.5	L
79-81518	11.5 – 14	13+	XL
Impax Diabetic Shoe Insole Accessory			
79-81523			S
79-81525			M
79-81527			L
79-81528			XL



PROCARE®

Med/Surg Shoe

Designed for a variety of post-op and post-trauma applications. Soft foam mesh upper with padded tricot lining. Semi-flexible skid resistant sole and comfortable contoured heel. Large overlapping toe area accommodates for post-operative dressing and swelling. Loop lock contact closure straps. Fits right or left foot. Latex free.

Indication: Post-op; Post-Trauma applications.

Male	Shoe Size Male	Female	Shoe Size Female	Size
79-81132	5 – 7	79-81142	2 – 4	Ped./XS
79-81133	7 – 9	79-81143	4 – 6	S
79-81135	9 – 11	79-81145	6 – 8	M
79-81137	11 – 13	79-81147	8 – 10	L
79-83118	13+			XL



PROCARE®

Rocker Cast Boot

Non-slip rocker sole to allow for more natural gait during ambulation. Sandal style with durable canvas upper. Forefoot and heel loop-lock contact closure for ease of application. Latex free.

Indication: For use with lower leg casted fractures.

	Length	Toe Width	Color	Size
79-81110	4.5"	2.5"	Ped print	3XS
79-81111	5.5"	3.0"	Ped print	2XS
79-81112	7.25"	3.5"	Black	XS
79-81113	9.25"	4.25"	Black	S
79-81114	10.5"	4.5"	Black	S/M
79-81115	11.0"	4.75"	Black	M
79-81116	12.0"	5.0"	Black	M/L
79-81117	12.5"	5.25"	Black	L
79-81118	13.5"	6.0"	Black	XL



PROCARE®

Darco™ SlimLine Cast Boot

Rocker bottom with weatherproof, stretchable upper and contact closure.

Indication: For use with lower leg casted fractures.

	Shoe Size Male	Shoe Size Female	Size
79-98432		3 – 5	XS
79-98433		5.5 – 8	S
79-98435	6 – 8	8.5 – 10	M
79-98437	10.5 – 12	10.5 – 12	L
79-98438	12.5 – 14	12+	XL

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
03AL	Enovis AIRCAST LEG BRACE, LEFT	\$126.29
03AR	Enovis AIRCAST LEG BRACE, RIGHT	\$126.29
03BL	Enovis AIRCAST LEG BRACE W/ANT PANEL, LEFT	\$147.26
03BR	Enovis AIRCAST LEG BRACE W/ANT PANEL, RIGHT	\$147.26
03CL	Enovis AIRCAST LEG BRACE, SMALL, LEFT	\$126.29
03CR	Enovis AIRCAST LEG BRACE SMALL, RIGHT	\$126.29
03DL	Enovis AIRCAST LEG BRACE WITH ANTERIOR PANEL, SMALL, LEFT	\$147.26
03DR	Enovis AIRCAST LEG BRACE WITH ANTERIOR PANEL, SMALL, RIGHT	\$147.26
79-72280	Enovis PROCARE TULI'S HEEL CUPS (REG, PAIR)	\$21.07
79-72281	Enovis PROCARE TULI'S HEEL CUPS (LARGE, PAIR)	\$21.07
79-81193	Enovis PROCARE SURROUND ANKLE,SM, 6"	\$37.62
79-81195	Enovis PROCARE SURROUND ANKLE, MED, 9"	\$37.62
79-81197	Enovis PROCARE SURROUND ANKLE, LG, 10"	\$37.62
79-81250	Enovis PROCARE STIRRUP ANKLE SPRT,UNIV,RT	\$106.53
79-81260	Enovis PROCARE STIRRUP ANKLE SPRT,UNIV,LT	\$106.53
79-81312-900	Enovis PROCARE LACE-UP ANKLE BRACE, XS	\$30.87
79-81313-900	Enovis PROCARE LACE-UP ANKLE BRACE, S	\$29.97
79-81315-900	Enovis PROCARE LACE-UP ANKLE BRACE, M	\$29.97
79-81317-900	Enovis PROCARE LACE-UP ANKLE BRACE, L	\$29.97
79-81318-900	Enovis PROCARE LACE-UP ANKLE BRACE, XL	\$30.87
79-81323	Enovis PROCARE EXCELERATOR ANKLE STIRRUP, 8"	\$57.60
79-81327	Enovis PROCARE EXCELERATOR ANKLE STIRRUP 9.5"	\$57.60
79-81330	Enovis PROCARE UNIVERSAL ANKLE BRACE	\$24.05
79-81332	Enovis PROCARE SQ. TOE POSTOP SHOE,IMPAX INSOLE, XS	\$15.55
79-81333	Enovis PROCARE SQ. TOE POSTOP SHOE,IMPAX INSOLE, S	\$15.55
79-81335	Enovis PROCARE SQ. TOE POSTOP SHOE,IMPAX INSOLE, M	\$15.55
79-81337	Enovis PROCARE SQ. TOE POSTOP SHOE,IMPAX INSOLE, L	\$15.55
79-81338	Enovis PROCARE SQ. TOE POSTOP SHOE,IMPAX INSOLE, XL	\$15.55
79-81352-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, XS	\$34.41
79-81353-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, S	\$34.41
79-81355-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, M	\$34.41
79-81357-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, L	\$34.41
79-81358-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, XL	\$34.41
79-81359-10-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, XXXL	\$36.97
79-81359-900	Enovis PROCARE STABILIZING ANKLE SUPPORT, XXL	\$35.45
79-81372-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,XS	\$25.10
79-81373-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,S	\$25.10
79-81375-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,M	\$25.10

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-81377-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,L	\$25.10
79-81378-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,XL	\$25.10
79-81379-900	Enovis PROCARE DOUBLE STRAP,ANKLE,SUPPORT,2XL	\$26.80
79-81403	Enovis PROCARE PROWEDGE NIGHT SPLINT,SM	\$65.20
79-81405	Enovis PROCARE PROWEDGE NIGHT SPLINT,MED	\$65.20
79-81407	Enovis PROCARE PROWEDGE NIGHT SPLINT,LG	\$65.20
79-81408	Enovis PROCARE PROWEDGE NIGHT SPLINT,XL	\$65.20
79-82302	Enovis PROCARE ANKLE SLEEVE,XS	\$29.72
79-82303	Enovis PROCARE ANKLE SLEEVE,S	\$29.72
79-82305	Enovis PROCARE ANKLE SLEEVE,M	\$29.72
79-82307	Enovis PROCARE ANKLE SLEEVE,L	\$29.72
79-82308	Enovis PROCARE ANKLE SLEEVE,XL	\$29.72
79-97753	Enovis PROCARE PLANTAR FASCIITIS NIGHT SPLT,S	\$70.46
79-97755	Enovis PROCARE PLANTAR FASCIITIS NIGHT SPLT,M	\$70.46
79-97757	Enovis PROCARE PLANTAR FASCIITIS NIGHT SPLT,L	\$70.46
79-97758	Enovis PROCARE PLANTAR FASCIITIS NIGHT SPL,XL	\$70.46
79-97863	Enovis PROCARE SURROUND GEL ANKLE,SMALL, 8.5"	\$28.29
79-97865	Enovis PROCARE SURROUND GEL ANKLE, ATH, 9"	\$28.29
79-97867	Enovis PROCARE SURROUND GEL ANKLE,REG, 10"	\$28.29
81-81102	Enovis PROCARE SILICONE HEEL CUPS, XS	\$61.45
81-81103	Enovis PROCARE SILICONE HEEL CUPS, S/M	\$61.45
81-81107	Enovis PROCARE SILICONE HEEL CUPS, L/XL	\$61.45

WALKING BRACES

AIRCRAFT®

Elite

The Elite has a durable, semi-rigid shell that helps support the limb while providing full-shell protection. Overlapping Duplex aircells line the shell providing intermittent pneumatic compression for 3x faster edema reduction¹. The front panel is pre-inflated and combined with the overlapping aircells, provides a 'total contact' fit. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Stable fracture of foot/ankle/lower leg; Severe ankle sprain; Post-operative immobilization; Reduction of edema (swelling) 3x faster¹.

Shoe Size			
	Male	Female	Size
01EP-XS	Up to 4	Up to 5	XS
01EP-S	4-7	5-8	S
01EP-M	7-10	8-11	M
01EP-L	10-13	11-15	L
01EP-XL	13+	15+	XL
Accessories			
01WXT-x	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E102x	Liner & Insole Replacement Kit, Elite		
01E104x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.

1. Conservative Therapy for Acute Lateral Ligament Lesions – Single Chamber vs. Two-Chamber Orthosis Systems; Schmidt, Mainers, Reintges, Lipke, Benesch, Gerngross- Surgery Dept of the Federal Arm Hospital, Ulm, Germany, 1999





AIRCAST®

Standard

The Standard has a durable, semi-rigid shell that supports the limb while providing protection. Housed within the shell are two adjustable distal aircells that provide compression and support the malleoli. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Stable fracture of foot/ankle/lower leg; Severe ankle sprain; Post-operative immobilization.

Shoe Size			
	Male	Female	Size
01EF-XS	Up to 4	Up to 5	XS
01EF-S	4-7	5-8	S
01EF-M	7-10	8-11	M
01EF-L	10-13	11-15	L
01EF-XL	13+	15+	XL
Accessories			
	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E102x	Liner & Insole Replacement Kit, Elite		
01E104x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.



AIRCAST®

Short

The Short has a durable, semi-rigid shell that helps support the limb while providing protection. The pre-inflated front panel and two adjustable distal aircells, housed within the shell, provide compression and support. SoftStrike technology absorbs and dissipates shock, while the light weight and lab-tested rocker sole combine to encourage a natural gait and a continuation of everyday activities.

Indication: Metatarsal fracture; Stress and stable fracture of the foot; Bunionectomy, hallux valgus; Soft tissue injury; Severe ankle sprain.

Shoe Size			
	Male	Female	Size
01ES-XS	Up to 4	Up to 5	XS
01ES-S	4-7	5-8	S
01ES-M	7-10	8-11	M
01ES-L	10-13	11-15	L
01ES-XL	13+	15+	XL
Accessories			
01WXT-x	Weather Cover, Tall		
01WXS-x	Weather Cover, Short		
0130A-x	Hygiene Cover		
01T-x	Toe Cover		
01P414	Brace-Lok™		
01K-x	Heel Wedges		
0129P	Tube Stretch Sock		XS
0129A	Tube Stretch Sock		Univ.
Replacement Parts			
01E100x	Liner & Insole Replacement Kit, Elite		
01E100x	Strap Replacement Kit, Elite & Standard		
01E1050	Extension Strap Kit		Univ.

AIRCAS[®]**XP Walker™ (extra pneumatic)**

The XP Walker has a lightweight, durable semi-rigid shell that helps support the limb while providing full-shell protection. Lining the shell are overlapping Duplex aircells providing intermittent pneumatic compression for edema reduction and callus formation. These aircells can be custom inflated for a "total contact" fit. Includes two socks and a hand bulb — Can be used with the Brace Lok™ system.

Indication: Stable foot and/or ankle fracture; severe ankle sprain; post-operative use.

	Shoe Size		Size
	Male	Female	
01P-P	Up to 4	Up to 5	Ped.
01P-S	4 – 7	5 – 8	S
01P-M	7 – 10	8 – 11	M
01P-L	10 – 13	11 – 15	L
01P-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
0193	Pneumatic Care Kit (2 socks, 1 hand bulb)		Univ.
0120	Replacement Hand Bulb		

AIRCAS[®]**FP Walker™ (foam pneumatic)**

The FP Walker has a lightweight, durable semi-rigid shell that supports the limb while providing protection. Housed within the shell are two adjustable distal aircells that provide compression and support the malleoli. Includes a hand bulb.

Indication: Stable foot and/or ankle fracture; severe ankle sprain; post-operative use.

	Shoe Size		Size
	Male	Female	
01F-P	Up to 4	Up to 5	Ped.
01F-S	4 – 7	5 – 8	S
01F-M	7 – 10	8 – 11	M
01F-L	10 – 13	11 – 15	L
01F-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
1120SC	Replacement Hand Bulb		Univ.



AIRCAST®

SP Walker™ (short pneumatic)

The SP Walker has a lightweight, durable semi-rigid shell with a pre-inflated anterior aircell. Two adjustable aircells located at the malleoli may be individually inflated for “total contact fit”. The generous foot base has ample room for dressings without sacrificing comfort. Includes a hand bulb.

Indication: Metatarsal fracture; forefoot and midfoot injury; acute or post-operative use; bunionectomy; soft tissue injury; severe ankle sprain.

	Shoe Size		Size
	Male	Female	
01A-P	Up to 4	Up to 5	Ped.
01A-S	4 – 7	5 – 8	S
01A-M	7 – 10	8 – 11	M
01A-L	10 – 13	11 – 15	L
01A-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
1120SC	Replacement Hand Bulb		Univ.



AIRCAST®

XP Diabetic Walker™ System

The XP Diabetic Walker incorporates an optimized rocker sole, semi-rigid shells, an innovative dual density insole and four aircells for maximum off-loading and “total contact” fit of the diabetic foot. Two extra insoles (including the dual-laminate Impax™ Grid insole), three stockings and one hand bulb with pressure gauge are included. Can be used with the Brace-Lok™ system.

Indication: Stable foot and/or ankle fracture; post-operative; hindfoot and heel foot protection, neuropathic ulcers, soft tissue injury.

	Shoe Size		Size
	Male	Female	
01PD-S	4 – 7	5 – 8	S
01PD-M	7 – 10	8 – 11	M
01PD-L	10 – 13	11 – 15	L
01PD-XL	13+	15+	XL
Accessories			
01159	Extension Straps		Univ.
01PG	Hand Bulb with Pressure Gauge		Univ.



PROCARE®

MiniTrax

MiniTrax™ is the only pediatric walking boot to feature kid-friendly design elements that make wearing a boot less intimidating. Beneath the fun exterior though, the MiniTrax is a serious boot with big benefits. Lightweight design includes anatomically shaped plastic uprights that help to facilitate a better fit. Easy to apply and adjust with a simple strapping and liner system. Low-profile rocker sole of the MiniTrax features a durable, non-slip tread tough enough to withstand the test of time—and of a kid.

Indication: Ideal following acute ankle sprains, metatarsal fractures stress and stable fractures of the lower leg.

MiniTrax Sizing Guide				
	Foot Length	Leg Length	Shoe Size	Size
79-95510	≤ 6"	6" - 7"	≤ 8.5	S (Ped)
79-95520	6" - 7"	≤ 11"	8.5 - 11.5	M (Ped)
79-955130	7" - 8"	≤ 14"	11.5 - 1.5	L (Youth)
Accessories				
79-95610	MiniTrax Replacement Liners			S (Ped)
79-95620	MiniTrax Replacement Liners			M (Ped)
79-95630	MiniTrax Replacement Liners			L (Youth)

PROCARE®

MaxTrax™ Air

Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. Adjustable medial/lateral air cell system inflates both sides equally to accommodate changes in edema during rehab and maximize patient compliance. Cushioned inner/outer sole with low profile rocker bottom provides increased stability and comfort. Soft nylon/foam liner is constructed to ensure no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walking brace to accommodate shorter tibia/fibula.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

		Shoe Size			
	Male	Female	Junior	Replacement	Liner
MaxTrax Air					
79-95412		≤4	≤5	79-95362	Ped./XS
79-95413	≤5	4.5 – 6		79-95363	S
79-95415	5.5 – 10	6.5 – 11		79-95365	M
79-95417	10.5 – 13.5	11.5 – 14.5		79-95367	L
79-95418	14 – 17			79-95368	XL
MaxTrax Air Ankle					
79-95422		≤4	≤5	79-95372	Ped./XS
79-95423	≤5	4.5 – 6		79-95373	S
79-95425	5.5 – 10	6.5 – 11		79-95375	M
79-95427	10.5 – 13.5	11.5 – 14.5		79-95377	L
79-95428	14 – 17			79-95378	XL
Anterior/Posterior and Medial/Lateral Shells					
79-95402					XS/S
79-95403					M/XL



PROCARE®**XcelTrax™ Air**

Designed for ideal support and comfort following trauma or post-operative procedure. The lightweight design and unique pre-shaped ergonomic frame help to modify the level of pressure on the injury at both the malleoli and lower leg. The soft nylon/foam liner has an integrated pump design making the pneumatic liner simple to inflate and adjust.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

Shoe Size				
	Male	Female	Replacement Liner	Size
XcelTrax Air				
79-95512	2 – 4	3.5 – 5	79-95532	XS
79-95513	4.5 – 7	6 – 8	79-95533	S
79-95515	7.5 – 10.5	8.5 – 11.5	79-95535	M
79-95517	10.5 – 12.5	11.5 – 13.5	79-95537	L
79-95518	12.5+	13.5+	79-95538	XL
XcelTrax Air Ankle				
79-95522	2 – 4	3.5 – 5	79-95342	XS
79-95523	4.5 – 7	6 – 8	79-95343	S
79-95525	7.5 – 10.5	8.5 – 11.5	79-95345	M
79-95527	10.5 – 12.5	11.5 – 13.5	79-95347	L
79-95528	12.5+	13.5+	79-95348	XL
Accessories				

**PROCARE®****Nextep™ Contour Air**

Provides protection and edema control following trauma or post-operative procedures. A medial/lateral dual air cell system inflates both sides of the liner equally and is easily adjusted with an integrated pump for proper compression, comfort, and custom fit. A secure, full wrap-around shell provides maximum support for ligament and tendon repair, osteotomy, pinning procedures, stress fractures or post-operative bunion surgery. Cushioned, non-slip rocker sole. Washable brush nylon foam liner with toe cover for maximum comfort.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

Shoe Size				
	Male	Female	Replacement Liner	Size
Nextep Contour Air				
79-95273	3.5 – 7.5	4.5 – 8.5	79-95363	S
79-95275	7 – 10	8 – 11.5	79-95365	M
79-95277	10 – 13	11 – 14	79-95367	L
Nextep Contour Air "Shortie"				
79-95443	3.5 – 7.5	4.5 – 8.5	79-95373	S
79-95445	7 – 10	8 – 11.5	79-95375	M
79-95447	10 – 13	11 – 14	79-95377	L
Accessories				
79-95440	Toe Cover 3/pk.			Univ.

PROCARE®**Nextep™ Contour 2 Air**

Provides protection and edema control following trauma or post-operative procedures. A medial/lateral dual air cell system inflates both sides of the liner equally and is easily adjusted with an integrated pump for proper compression, comfort, and custom fit. Cushioned, non-slip rocker sole. Washable brush nylon foam liner for maximum comfort.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

Shoe Size				
	Male	Female	Replacement Liner	Size
Nextep Contour 2 Air				
79-95173	3.5 – 7.5	4.5 – 8.5	79-95363	S
79-95175	7 – 10	8 – 11.5	79-95365	M
79-95177	10 – 13	11 – 14	79-95367	L
Accessories				
79-95440	Toe Cover 3/pk.			Univ.



PROCARE®

MaxTrax™

Designed with a wider foot bed for comfort and stability following lower leg, foot and ankle trauma and post-op procedures. A low profile rocker bottom helps promote natural gait and reduce plantar pressures. Cushioned inner and outer sole provides shock absorption to help improve patient comfort during ambulation. Soft nylon/foam liner is constructed to ensure no seams are in critical post-surgical areas. Malleable uprights breakdown on standard height walking brace to accommodate shorter tibia/fibula.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
MaxTrax					
79-95322		≤4	≤5	79-95192	Ped./XS
79-95323	≤5	4.5 – 6		79-95193	S
79-95325	5.5 – 10	6.5 – 11		79-95195	M
79-95327	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95328	14 – 17			79-95198	XL
MaxTrax Ankle					
79-95342		≤4	≤5	79-95182	Ped./XS
79-95343	≤5	4.5 – 6		79-95183	S
79-95345	5.5 – 10	6.5 – 11		79-95185	M
79-95347	10.5 – 13.5	11.5 – 14.5		79-95187	L
79-95348	14 – 17			79-95188	XL
Anterior/Posterior and Medial/Lateral Shells					
79-95402					XS/S
79-95403					M/XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

XcelTrax™

Designed for ideal support and comfort following trauma or post-operative procedure. The lightweight design and unique pre-shaped ergonomic frame help to modify the level of pressure on the injury at both the malleoli and lower leg.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size			Replacement Liner	Size
	Male	Female			
XcelTrax					
79-95492	2 – 4	3.5 – 5	79-95552		XS
79-95493	4.5 – 7	6 – 8	79-95553		S
79-95495	7.5 – 10.5	8.5 – 11.5	79-95555		M
79-95497	10.5 – 12.5	11.5 – 13.5	79-95557		L
79-95498	12.5+	13.5+	79-95558		XL
XcelTrax Ankle					
79-95502	2 – 4	3.5 – 5	79-95562		XS
79-95503	4.5 – 7	6 – 8	79-95563		S
79-95505	7.5 – 10.5	8.5 – 11.5	79-95565		M
79-95507	10.5 – 12.5	11.5 – 13.5	79-95567		L
79-95508	12.5+	13.5+	79-95568		XL
Accessories					
79-95210	Walker Extension Pad				Univ.

**PROCARE®****Nextep™ Contour**

Lightweight, durable, one-piece wrap-around shell, provides maximum immobilization and support following ligament and tendon repair, osteotomy, pinning procedures, stress fractures or post-operative bunion surgery. Unique design of circumferential shell also accommodates larger calf sizes. Low profile, cushioned rocker sole absorbs impact of heel strike to improve comfort during ambulation. Washable nylon foam liner with removable toe cover.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95063	3.5 – 7.5	4.5 – 8.5	79-95163	S
79-95065	7 – 10	8 – 11.5	79-95165	M
79-95067	10 – 13	11 – 14	79-95167	L
Nextep Contour Shortie				
79-95083	3.5 – 7.5	4.5 – 8.5		S
79-95085	7 – 10	8 – 11.5		M
79-95087	10 – 13	11 – 14		L
Accessories				
79-95210	Walker Extension Pad			Univ.
79-95440	Toe Cover 3/pk.			Univ.

**PROCARE®****SideKICK™**

The low profile design provides comfortable stability for acute ankle sprains, soft tissue and stress fractures to the lower leg and ankle, and stable fractures of the ankle. Other features include fixed, malleable uprights, a breathable fiber leg wrap, removable heel pads, and non-skid rocker sole.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95033	≤6	≤7	79-95123	S
79-95035	6 – 10	7 – 11	79-95125	M
79-95037	10+	11+	79-95127	L

**PROCARE®****Nextep™ Contour 2**

Lightweight, durable, one-piece wrap-around posterior shell, provides maximum immobilization and support. Brushed nylon/foam liner, low profile, cushioned non-skid sole. Unique design accommodates larger calf sizes. Rocker sole.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-95073	3.5 – 7.5	4.5 – 8.5	79-95163	S
79-95075	7 – 10	8 – 11.5	79-95165	M
79-95077	10 – 13	11 – 14	79-95167	L
Accessories				
79-95210	Walker Extension Pad			Univ.
79-95440	Toe Cover 3/pk.			Univ.

**PROCARE®****ProSTEP™**

Lightweight, sturdy one-piece design with open heel. Non-skid rocker sole and washable fiber leg wrap with removable heel pads.

Indication: Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

	Shoe Size		Replacement Liner	Size
	Male	Female		
79-98793	≤6	≤7	79-95123	S
79-98795	6 – 10	7 – 11	79-95125	M
79-98797	10+	11+	79-95127	L



PROCARE®

MaxTrax™ Air ROM

Protected range-of-motion and adjustable air cell system provides stability and helps decrease pain and edema following trauma or post-operative procedures. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion. Cushioned inner/outer sole with low profile rocker bottom provides increased stability and comfort. Soft nylon/foam liner.

Indication: Ideal following Achilles tendon repair for controlled rehab, and stable fractures of the lower leg.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
79-95432		≤4	≤5	79-95362	Ped./XS
79-95433	≤5	4.5 – 6		79-95363	S
79-95435	5.5 – 10	6.5 – 11		79-95365	M
79-95437	10.5 – 13.5	11.5 – 14.5		79-95367	L
79-95438	14 – 17			79-95368	XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

MaxTrax™ ROM

Provides protected range-of-motion following trauma or post-operative procedure. Ideal following Achilles tendon repair for controlled rehab. ROM settings between 45° plantar-flexion and 30° dorsi-flexion in 7.5° increments. Uprights lock in fixed positions of 0°, 7.5°, 15°, 22.5°, and 30° plantar- and dorsi-flexion. Malleable uprights breakdown on standard height to accommodate shorter tibia/fibula. Cushioned inner/outer sole. Low profile rocker bottom.

Indication: Ideal following Achilles tendon repair for controlled rehab, and stable fractures of the lower leg.

	Shoe Size			Replacement Liner	Size
	Male	Female	Junior		
MaxTrax ROM					
79-95332		≤4	≤5	79-95192	Ped./XS
79-95333	≤5	4.5 – 6		79-95193	S
79-95335	5.5 – 10	6.5 – 11		79-95195	M
79-95337	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95338	14 – 17			79-95198	XL
MaxTrax ROM Ankle					
79-95352		≤4	≤5	79-95192	Ped./XS
79-95353	≤5	4.5 – 6		79-95193	S
79-95355	5.5 – 10	6.5 – 11		79-95195	M
79-95357	10.5 – 13.5	11.5 – 14.5		79-95197	L
79-95358	14 – 17			79-95198	XL
Accessories					
79-95210	Walker Extension Pad				Univ.



PROCARE®

MaxTrax™ Diabetic Walker

Designed with a wider foot bed for increased stability and a soft nylon/foam liner with no seams in critical post-surgical areas and/or areas where neuropathy may occur. Tri-laminate insole is designed with Impax™ Grid technology to redistribute load away from ulcers. Low profile rocker bottom cushioned outer sole for shock absorption. Optional Toe Cover and a set of Upright Pads for additional patient comfort.

Indication: Ideal for replacement of total contact casting; treatment of pre-ulcerative or ulcerative conditions of the plantar surface of the foot; charcot foot.

Shoe Size				
	Male	Female	Replacement Liner	Size
79-95453	≤5	4.5 – 6.5	79-95473	S
79-95455	5.5 – 10	6.5 – 11	79-95475	M
79-95457	10.5 – 13.5	11.5 – 14.5	79-95477	L
79-95458	14 – 17		79-95478	XL
79-95463	Impax Foam Insole			S
79-95465	Impax Foam Insole			M
79-95467	Impax Foam Insole			L
79-95468	Impax Foam Insole			XL
Accessories				
79-95450	MaxTrax Diabetic Toe Cover			Univ.
79-95460	MaxTrax Upright Pads (1 pair)			Univ.
79-95210	Walker Extension Pad			Univ.



PROCARE®

Achilles Wedge

4" neoprene padded foam wedge fits into most closed heel walkers. Each section is 3/4" in height and may be "peeled-away" during Achilles Tendon rehab.

Size	
79-95130	Univ.



PROCARE®

Weather Cover

The ProCare Walking Brace Weather Cover is designed to repel the elements and to keep the walking brace clean and dry. The Velcro closures on the back of the cover allow for easy application. Open bottom allows sure-footing for ambulatory patients.

Size	
79-95093	S
79-95095	M
79-95097	L

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
01A-L	Enovis AIRCAST SP-WALKER, LARGE	\$147.77
01A-M	Enovis AIRCAST SP-WALKER, MEDIUM	\$147.77
01A-P	Enovis AIRCAST SP- WALKER PEDIATRIC	\$147.77
01A-S	Enovis AIRCAST SP-WALKER, SMALL	\$147.77
01A-XL	Enovis AIRCAST SP-WALKER, X-LARGE	\$147.77
01EF-L	Enovis AIRCAST AIRSELECT, STANDARD, LARGE	\$121.47
01EF-M	Enovis AIRCAST AIRSELECT, STANDARD, MEDIUM	\$121.47
01EF-S	Enovis AIRCAST AIRSELECT, STANDARD, SMALL	\$121.47
01EF-XL	Enovis AIRCAST AIRSELECT, STANDARD, XLARGE	\$121.47
01EF-XS	Enovis AIRCAST AIRSELECT, STANDARD, XSMALL	\$121.47
01EP-L	Enovis AIRCAST AIRSELECT, ELITE, LARGE	\$172.52
01EP-M	Enovis AIRCAST AIRSELECT, ELITE, MEDIUM	\$172.52
01EP-S	Enovis AIRCAST AIRSELECT, ELITE, SMALL	\$172.52
01EP-XL	Enovis AIRCAST AIRSELECT, ELITE, XLARGE	\$172.52
01EP-XS	Enovis AIRCAST AIRSELECT, ELITE, XSMALL	\$172.52
01ES-L	Enovis AIRCAST AIRSELECT, SHORT, LARGE	\$113.80
01ES-M	Enovis AIRCAST AIRSELECT, SHORT, MEDIUM	\$113.80
01ES-S	Enovis AIRCAST AIRSELECT, SHORT, SMALL	\$113.80
01ES-XL	Enovis AIRCAST AIRSELECT, SHORT, XLARGE	\$113.80
01ES-XS	Enovis AIRCAST AIRSELECT, SHORT, XSMALL	\$113.80
01F-L	Enovis AIRCAST FP WALKER, LARGE	\$158.59
01F-M	Enovis AIRCAST FP WALKER, MEDIUM	\$158.59
01F-P	Enovis AIRCAST FOAM WALKER, PEDIATRIC	\$158.59
01F-S	Enovis AIRCAST FP WALKER, SMALL	\$158.59
01F-XL	Enovis AIRCAST FOAM WALKER, X-LARGE	\$158.59
01PD-L	Enovis AIRCAST XP DIABETIC WALKER SYSTEM, L	\$240.76
01PD-M	Enovis AIRCAST XP DIABETIC WALKER SYSTEM, M	\$240.76
01PD-S	Enovis AIRCAST XP DIABETIC WALKER SYSTEM, S	\$240.76
01PD-XL	Enovis AIRCAST PNEUMATIC WLKR DIABETIC SYS,XL	\$240.76
01P-L	Enovis AIRCAST XP WALKER, LARGE	\$228.60
01PL-210	Enovis AIRCAST XP WALKER, L	\$163.73
01P-M	Enovis AIRCAST XP WALKER, MEDIUM	\$228.60
01PM-210	Enovis AIRCAST XP WALKER, M	\$152.54
01P-P	Enovis AIRCAST PNEUMATIC WALKER, PEDIATRIC	\$228.60
01P-S	Enovis AIRCAST XP WALKER, SMALL	\$228.60
01P-XL	Enovis AIRCAST PNEUMATIC WALKER, X-LARGE	\$234.53
79-95033	Enovis PROCARE SIDEKICK LOW PROFILE WALKER,S	\$131.22
79-95035	Enovis PROCARE SIDEKICK LOW PROFILE WALKER,M	\$131.22

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-95037	Enovis PROCARE SIDEKICK LOW PROFILE WALKER,L	\$131.22
79-95063	Enovis PROCARE NEXTEP CONTOUR WALKER,SMALL	\$113.14
79-95065	Enovis PROCARE NEXTEP CONTOUR WALKER,MEDIUM	\$113.14
79-95067	Enovis PROCARE NEXTEP CONTOUR WALKER,LARGE	\$113.14
79-95073	Enovis PROCARE NEXTEP CONTOUR2 WALKER,SMALL	\$59.19
79-95075	Enovis PROCARE NEXTEP CONTOUR2 WALKER,MEDIUM	\$59.19
79-95077	Enovis PROCARE NEXTEP CONTOUR2 WALKER,LARGE	\$59.19
79-95093	Enovis PROCARE WEATHER COVER,PROCARE,WALKER,S	\$32.16
79-95095	Enovis PROCARE WEATHER COVER,PROCARE,WALKER,M	\$32.16
79-95097	Enovis PROCARE WEATHER COVER,PROCARE,WALKER,L	\$32.16
79-95130	Enovis PROCARE ACHILLES WEDGE	\$21.99
79-95173	Enovis PROCARE NEXTEP CONTOUR 2 AIR WALKER, S	\$74.83
79-95175	Enovis PROCARE NEXTEP CONTOUR 2 AIR WALKER, M	\$74.83
79-95177	Enovis PROCARE NEXTEP CONTOUR 2 AIR WALKER, L	\$74.83
79-95273	Enovis PROCARE NEXTEP CONTOUR W/AIR,S	\$117.75
79-95275	Enovis PROCARE NEXTEP CONTOUR W/AIR,M	\$117.75
79-95277	Enovis PROCARE NEXTEP CONTOUR W/AIR,L	\$117.75
79-95310	Enovis PROCARE MAXTRAX STRAP EXTENSION KIT	\$8.12
79-95322	Enovis PROCARE MAXTRAX WALKER,XS	\$125.68
79-95323	Enovis PROCARE MAXTRAX WALKER,S	\$125.68
79-95325	Enovis PROCARE MAXTRAX WALKER,M	\$125.68
79-95327	Enovis PROCARE MAXTRAX WALKER,L	\$125.68
79-95328	Enovis PROCARE MAXTRAX WALKER,XL	\$125.68
79-95332	Enovis PROCARE MAXTRAX ROM WALKER,XS	\$116.17
79-95333	Enovis PROCARE MAXTRAX ROM WALKER,S	\$116.17
79-95335	Enovis PROCARE MAXTRAX ROM WALKER,M	\$116.17
79-95337	Enovis PROCARE MAXTRAX ROM WALKER,L	\$116.17
79-95338	Enovis PROCARE MAXTRAX ROM WALKER,XL	\$116.17
79-95342	Enovis PROCARE MAXTRAX ANKLE,WALKER,XS	\$120.51
79-95343	Enovis PROCARE MAXTRAX ANKLE,WALKER,S	\$120.51
79-95345	Enovis PROCARE MAXTRAX ANKLE,WALKER,M	\$120.51
79-95347	Enovis PROCARE MAXTRAX ANKLE,WALKER,L	\$120.51
79-95348	Enovis PROCARE MAXTRAX ANKLE,WALKER,XL	\$120.51
79-95352	Enovis PROCARE MAXTRAX ANKLE ROM,WALKER,XS	\$72.97
79-95353	Enovis PROCARE MAXTRAX ANKLE ROM,WALKER,S	\$72.97
79-95355	Enovis PROCARE MAXTRAX ANKLE ROM,WALKER,M	\$72.97
79-95357	Enovis PROCARE MAXTRAX ANKLE ROM,WALKER,L	\$72.97
79-95358	Enovis PROCARE MAXTRAX ANKLE ROM,WALKER,XL	\$72.97

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-95412	Enovis PROCARE MAXTRAX AIR WALKER,XS	\$134.55
79-95413	Enovis PROCARE MAXTRAX AIR WALKER,S	\$134.55
79-95413-012	Enovis PROCARE MAXTRAX AIR WALKER, S	\$54.46
79-95415	Enovis PROCARE MAXTRAX AIR WALKER,M	\$134.55
79-95415-012	Enovis PROCARE MAXTRAX AIR WALKER, M	\$54.46
79-95417	Enovis PROCARE MAXTRAX AIR WALKER,L	\$134.55
79-95417-012	Enovis PROCARE MAXTRAX AIR WALKER, L	\$54.46
79-95418	Enovis PROCARE MAXTRAX AIR WALKER,XL	\$134.55
79-95422	Enovis PROCARE MAXTRAX AIR ANKLE WALKER,XS	\$131.86
79-95423	Enovis PROCARE MAXTRAX AIR ANKLE WALKER,S	\$131.86
79-95425	Enovis PROCARE MAXTRAX AIR ANKLE WALKER,M	\$131.86
79-95427	Enovis PROCARE MAXTRAX AIR ANKLE WALKER,L	\$131.86
79-95428	Enovis PROCARE MAXTRAX AIR ANKLE WALKER,XL	\$131.86
79-95432	Enovis PROCARE MAXTRAX ROM AIR WALKER,XS	\$126.71
79-95433	Enovis PROCARE MAXTRAX ROM AIR WALKER,S	\$126.71
79-95435	Enovis PROCARE MAXTRAX ROM AIR WALKER,M	\$126.71
79-95437	Enovis PROCARE MAXTRAX ROM AIR WALKER,L	\$126.71
79-95438	Enovis PROCARE MAXTRAX ROM AIR WALKER,XL	\$126.71
79-95453	Enovis PROCARE MAXTRAX DIABETIC, S	\$121.89
79-95455	Enovis PROCARE MAXTRAX DIABETIC, M	\$121.89
79-95457	Enovis PROCARE MAXTRAX DIABETIC, L	\$121.89
79-95458	Enovis PROCARE MAXTRAX DIABETIC, XL	\$121.89
79-98793	Enovis PROCARE PRO-STEP WALKER,SMALL	\$107.29
79-98795	Enovis PROCARE PRO-STEP WALKER,MEDIUM	\$107.29
79-98797	Enovis PROCARE PRO-STEP WALKER,LARGE	\$107.29

PATIENT & PRESSURE CARE

PROCARE®**Body Wedge**

Stabilizes the patient securely in several positions without slipping or sliding. Supports arm or leg when lying on side or at foot of bed to prevent sliding in supine positions. Can also be used as a back support, footrest or blanket support. Comes compressed to reduce package sizing.

Indication: Pressure ulcer prevention, post-operative stabilization and positioning.

	Qty.	Length	Height	Width
79-90860	6/bx.	11.69"	8.375"	20"





PROCARE®

Hip Abduction Pillows

Designed for use following hip surgery when immobilization or post-operative positioning is required. All-foam construction for maximum support. Soft contact closure foam laminate straps help keep positioner in place and provide easy adjustment. Universal designs accommodate wide variety of patient sizes. 79-90180 has a 4" tear-off block on both ends of pillow for adjustment. Comes compressed to reduce package sizing.

Indication: Help keep hips in abduction (away from center of the body) following injury or surgery.

	Length	Height	Width Bottom	Size
79-90173	18"	6"	12"	S
79-90175	23"	6"	15"	M
79-90177	26"	6"	18"	L
79-90180	Max. 27" / Min. 19"	3"	6.5"	Univ.
79-90190*	18"	6"	18"	Univ.

*Foam wedge with removable, flannel lined washable cover and pile foam straps with contact closure.



PROCARE®

Leg Elevator

Full foam construction to elevate and maintain proper knee and lower leg position at 45° angle. Comes compressed to reduce package sizing.

Indication: Ideal for patients with leg casts, or those requiring stable elevation and internal/external rotation control, post-trauma or surgical procedures.

	Length	Height	Width
79-90191	31.5"	10"	10.25"



PROCARE®

Arm Elevation Pillow

Sculptured foam support aids in the reduction of pre/post-cast or post-surgical swelling. Provides comfortable, effective support with minimal restraint. Comes compressed to reduce package sizing. Excellent for use on wrist or shoulder surgery and mastectomy patients.

Indication: Intended for post-surgical or post-injury elevation of the hand and arm to promote venous return to minimize swelling.

	Height	Width	Length
79-90500	16.75"	11.5"	12.75"



PROCARE®


Donuts

Versatile positioners that provide anatomical support and comfort. May be used as a head rest; for protection of the ears and nerves of the head and face. Additional indications include protection of the elbows, knees and heels.

Indication: Cushions the head, ear or iliac crest. May be used to support elbows, heels and knees.

	Qty.	Height	Diameter
79-90880* Donuts w/Center (6", 4.5", 3" center)	36/bx.	2"	9"
79-90890 Donuts w/o (3" Center)	36/bx.	2"	7"
79-90891 Donuts w/o (4.5" Center)	36/bx.	2"	9"

*Donuts w/center (photo shown above).



PROCARE®

Foam Arm Boards

Helps maintain hemodynamic during surgery and enhances pressure redistribution. The Foam Arm Boards protect bony prominences and pressure areas such as elbows, knees, and heels. Comes compressed to reduce package sizing.

Indication: Supports entire arm during surgery. Can be used on an operating room table to enhance pressure redistribution.

	Qty.	Length	Height	Width
79-90820	12 pair/bx.	16"	1.75"	7.5"




PROCARE®

Head Cradles

Provide comfortable immobilization of head and neck. May be used for lateral, supine or prone positions. Designed to conform to sensitive facial features while maintaining head position. Comes compressed to reduce package sizing.

Indication: Provides comfortable immobilization of head and neck. May be used for lateral, supine or prone positions.

	Qty.	Length	Height	Width
79-90872 Head Cradle w/o Trach Slit	12/bx.	8"	4.5"	9"
79-90873 Head Cradle w/Trach Slit	36/bx.	8"	4"	9"



PROCARE®

Laminectomy Arm Cradles with Trough

Used to support the body off the chest while in a prone position. Offers cushioning for arms, shoulders, elbows and ulnar nerves. Design allows for varying arm positions. Comes compressed to reduce package sizing.

Indication: Offers cushioning for arms, shoulders, elbows, and ulnar nerves during laminectomy procedures or other prone positioned surgeries.

	Qty.	Length	Height	Width
79-90810	5 pair/bx.	23.5"	3"	5"



PROCARE®

Ulnar Nerve Protectors

Protects against pinching and entrapment of the ulnar nerve. Provides cushioning on bony prominences of the arm. Available in standard and premium convoluted foam density.

Indication: Ideal for protecting against pinching and entrapment of the ulnar nerve. Provides cushioning on bony prominences of the arm.

	Qty.	Length	Height	Width
Standard 1.0 Density				
79-90800	36 pair/bx.	14.75"	2"	6.4"
Premium 1.5 Density				
79-90801*	36 pair/bx.	14.75"	1.75"	6.4"

*Premium 1.5 Density (photo shown above).



PROCARE®

Prone Head Positioners

Provides comfortable immobilization of head and neck. May be used for prone, lateral or supine positions. Designed to conform to sensitive facial features while maintaining head position. The Premium Prone Head Positioner maximizes comfort by using an extra soft comfortable support, reducing the pressure on the patients eyes and face. Comes compressed to reduce package sizing.

Indication: Provides comfortable immobilization of head and neck. May be used for lateral, supine or prone positions.

	Qty.	Length	Height	Width
79-90870	12/bx.	8.25"	4.125"	8.25"
79-90871*	12/bx.	10.5"	7"	9.5"

*Premium Prone Head Positioner (photo shown above).



PROCARE®

***Non-TAA Compliant**

Flex-Master® Non Sterile Clip Closure Bandage

Latex free, 100% woven cotton/elastic minimizes moisture retention and allergic reactions. Available in standard and double lengths.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
Standard Length			
79-98843	10	5.5 yds.	2"
79-98845	10	5.5 yds.	3"
79-98847	10	5.5 yds.	4"
79-98848	10	5.5 yds.	6"
Double Length			
79-98857	12	11 yds.	4"
79-98859	12	11 yds.	6"



PROCARE®

Flex-Master® Sterile Clip Closure Bandage

Standard and double lengths. Latex free, woven cotton/elastic construction with clip closure.

Indication: Ideal for post-surgical applications such as thoracic, knee, amputations, mastectomy, vascular, shoulder and burns.

	Per Pk.	Length	Width
Standard Length			
79-98723	20	5.5 yds.	2"
79-98725	20	5.5 yds.	3"
79-98727	20	5.5 yds.	4"
79-98728	20	5.5 yds.	6"
Double Length			
79-98717	20	11 yds.	4"
79-98718	20	11 yds.	6"



PROCARE®

Elastic Bandages - Self Closure

Constructed of breathable elastic. Each bandage is individually wrapped and may be secured by hook and loop closure or with two clip fasteners. Latex free.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
79-98233	10	5 yds.	2"
79-98235	10	5 yds.	3"
79-98237	10	5 yds.	4"
79-98238	10	5 yds.	6"



PROCARE®

Elastic Bandages w/Clip Closure

Cotton/Elastic wrap stretches to a full 4.5 yards and secures with two clips. Washable and auto clavable. Latex free.

Indication: Ideal for groin wraps, ice packs, knees, and ankles.

	Per Pk.	Length	Width
79-98583	10	4.5 yds.	2"
79-98585	10	4.5 yds.	3"
79-98587	10	4.5 yds.	4"
79-98588	10	4.5 yds.	6"



PROCARE®

Latex Rubber Bandage

Constructed of natural latex and designed to force blood out of an extremity by compressing the superficial vessels. Available in 3", 4", and 6" widths.

Indication: Used in conjunction with a tourniquet cuff to help accomplish exsanguinations (removal of blood) from a limb.

	Length	Width
79-98900	144"	3"
79-98910	144"	4"
79-98920	144"	6"



PROCARE®

Heel and Elbow Protectors

Constructed of a breathable interfacing with fiberfill lining. Three-layer comfort cushion helps to eliminate pressure points. Anatomical design provides customized patient comfort. Hook and loop closure for ease of application and removal. For protection of the elbow and heel for injured and/or bedridden patients.

Indication: Ideal for pressure reduction of heel and/or elbow and pressure protection.

	Style	Qty.	Length
79-81001	Personal	Bulk Pairs, 18 Pair	Univ.
79-81002	Prevent	Indiv. Pairs, 18 Pair	Univ.



PROCARE®

Catheter Strap

2" wide flannel lined elastic strap with contact closure. Special tab closure maintains tube position without occlusion.

Indication: Commonly used for catheter stabilization.

	Qty.	Length	Size
79-89200	12/pk.	24"	Univ.



PROCARE®

Heel/Elbow Protector - Poly/Pile

Choice of polyester pile or vented foam with velfoam contact closure straps. Sold in pairs.

Indication: Ideal for pressure reduction of heel and/or elbow and pressure protection.

	Qty.	Size
79-81003	1 pair	S
79-81005	1 pair	M
79-81007	1 pair	L



PROCARE®

Podous Boot

Lightweight, quality plastic shell with poly/pile liner help treat and prevent lower extremity disorders associated with trauma or immobility including pressure necrosis, ankle contractures, and post-operative hip abduction or rotation. Adjustable toe extension helps protect toes from bedding. Anti-rotation bar may be repositioned to assist with internal or external rotation. Removable transfer pad on dorsal side for limited patient ambulating. Fits right or left foot. Latex free.

Indication: Helps treat and prevent lower extremity disorders including pressure necrosis, ankle contractures, and post-operative hip abduction or rotation.

	Shoe Size			Length
	Male	Female	Calf Circ.	
79-90550	≤ 10	≤ 11	≤ 22"	Regular
79-90551	10.5 +	11.5 +	≤ 25"	XL
Replacement Liners				
79-90560	Liner Only	3/pk.		Regular
79-90561	Liner Only	3/pk.		XL



PROCARE®

Heel/Elbow Protectors - Universal

Choice of polyester pile or vented foam with velfoam contact closure straps. One size fits all. Sold in pairs.

Indication: For protection of the elbow and heel for injured and/or bedridden patients.

	Description	Qty.	Size
79-81040	Poly/Pile	1 Pair	Univ.
79-81050	Foam	1 Pair	Univ.
79-99740	Foam	3 Pair	Univ.



PROCARE®

Drop Foot Boot

90° upright frame with internal/external rotation in 10° increments helps control foot drop. Disposable poly/pile liner provides toe and heel protection.

Indication: Prevents Achilles contracture. Ideal for ankle contractures, muscle tightening, correct plantar flexion contractures.

	Description	Size
79-90160	Drop Foot Boot, Complete	Univ.
79-90170	Disposable Liner	Univ.
79-90161	Frame Only	Univ.



PROCARE®

Pre-Vent™ Ulnar Nerve Protector

Constructed of a breathable interfacing with fiberfill lining. Provides cushioned support for patient limb. Fluid repellent outer liner. Hook and loop closure for ease of application and removal. Can be used as post-operative support, insulated positioner for cryotherapy.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91100	18 Pair	Univ.



PROCARE®

Pre-Vent™ O.R. Table Pad

Constructed of fluid-attracting face fabric with fiberfill lining. Impervious backing traps fluid, preventing pooling and expediting post-operative clean-up. Three section design can be customized to accommodate any surgical position. Does not interrupt heat transfer from warming blankets.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91110	4/pk.	Univ.



PROCARE®

Pre-Vent™ Knee Crutch Pad

Constructed of a breathable interfacing with fiberfill lining. Provides cushioned support for patient limb. Fluid repellent outer liner. Hook and loop closure for ease of application and removal. Can be used as post-operative support, insulated positioner for cryotherapy applications.

Indication: Traps fluids from patient, prevents pooling on O.R. surface area and expedites post-operative clean-up.

	Qty.	Size
79-91120	12 Pair	Univ.



PROCARE®

Adjustable Anodized Aluminum Crutches

Packaged in pairs complete with tips, grips, and pads. Adjustable by using wing nut. 300 pound maximum weight. Latex free.

Indication: Can be used for many conditions including injuries, arthritis, or problems with balance.

	Approximate User Height	Size
79-91333	4'6" – 5'2"	Youth/S
79-91335	5'1" – 5'9"	Adult/M
79-91337	5'10" – 6'6"	Tall/L

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-81001	Enovis PROCARE PERSONAL HEEL / ELBOW, 18 PR	\$157.12
79-81003	Enovis PROCARE HEEL/ELBOW PROTECTOR,S,PAIR	\$20.52
79-81040	Enovis PROCARE KODEL HEEL/ELBOW PROT.,PR	\$9.90
79-81050	Enovis PROCARE FOAM HEEL/ELBOW PROT.,PR	\$9.34
79-89200	Enovis PROCARE CATHETER STRAP, ADULT (12PK)	\$40.81
79-90160	Enovis PROCARE DROP FOOT BOOT - COMPLETE	\$65.23
79-90161	Enovis PROCARE DROP FOOT BOOT - FRAME ONLY	\$51.32
79-90170	Enovis PROCARE DROP FOOT BOOT - LINER ONLY	\$30.38
79-90173	Enovis PROCARE HIP ABDUCTION PILLOW,S	\$51.79
79-90175	Enovis PROCARE HIP ABDUCTION PILLOW,M	\$51.79
79-90177	Enovis PROCARE HIP ABDUCTION PILLOW,L	\$51.79
79-90180	Enovis PROCARE SLIMLINE HIP ABDUC.PILLOW,UNIV	\$34.37
79-90191	Enovis PROCARE LEG ELEVATOR, UNIV	\$46.51
79-90500	Enovis PROCARE ARM ELEVATION PILLOW	\$52.61
79-90550	Enovis PROCARE PODOUS BOOT, REG, UNIV	\$105.47
79-90551	Enovis PROCARE PODOUS BOOT, XL, UNIV	\$105.47
79-90560	Enovis PROCARE PODOUS BT LNR,REG,UNV,3PK	\$52.19
79-90561	Enovis PROCARE PODOUS BT LNR,XL,UNIV,3PK	\$48.79
79-90800	Enovis PROCARE ULNAR NERVE PROTECTORS 1.0 DENSITY, 36 PR/BX	\$135.82
79-90801	Enovis PROCARE ULNAR NERVE PROTECTORS 1.5 DENSITY, 36 PR/BX	\$212.26
79-90810	Enovis PROCARE LAMINECTOMY ARM CRADLE/TROUGH 5 PR/BX	\$115.49
79-90820	Enovis PROCARE FOAM ARM BOARD PADS7.5"X16"X2", 12 PR/BX	\$37.52
79-90860	Enovis PROCARE BODY WEDGE 11.5X20X7.5, 6/BX	\$140.06
79-90870	Enovis PROCARE PRONE HEAD POSITIONER 8.25"X8.25"4.125", 12/BX	\$124.71
79-90871	Enovis PROCARE PREMIUM PRONE HEAD POSITIONER 10.5"X9.5"X7", 12/BX	\$163.46
79-90872	Enovis PROCARE HEAD CRADLE W/O TRACH SLIT 9"X8"X4.25", 12/BX	\$105.38
79-90873	Enovis PROCARE HEAD CRADLE W/ TRACH SLIT 9"X8"X4", 36/BX	\$262.13
79-90880	Enovis PROCARE DONUT W/ CENTERS 6"X4.5"X3", 36/BX	\$125.79
79-90881	Enovis PROCARE DONUT WITH OUT CENTER, 9 IN X 4.5 IN X 1.5 IN, 36BX	\$63.60
79-90882	Enovis PROCARE DONUT WITH OUT CENTER, 7 IN X 3 IN 1.5 IN, 36BX	\$51.39
79-90891	Enovis PROCARE DONUT W/O CENTER 9" DIAM. (4.5" CENTER), 36/BX	\$102.54
79-91100	Enovis PROCARE PRE-VENT ULNAR NERVE PROT,18PR	\$160.80
79-91110	Enovis PROCARE PRE-VENT O.R. TABLE PAD, 4 PK	\$102.49
79-91120	Enovis PROCARE PRE-VENT KNEE CRUTCH, 12 PR	\$177.16
79-91333-900	Enovis PROCARE ADJ ALUM CRUTCH, SM, PAIR	\$56.21
79-91335-900	Enovis PROCARE ADJ ALUM CRUTCH, MED, PAIR	\$56.21
79-91337-900	Enovis PROCARE ADJ ALUM CRUTCH, LG, PAIR	\$56.21

FRACTURE BRACING & SPLINTS

PROCARE®

Colles Splint

Pre-molded, padded aluminum splint is designed for wrist support and immobilization prior to surgery or casting. Ventilation holes provide air circulation for patient comfort.

Indication: Immobilization of the wrist in the initial short and or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-72112	79-72122	5.75" x 2.5"	Ped/XS
79-72113	79-72123	7" x 3.5"	S
79-72115	79-72125	8.5" x 3.5"	M
79-72117	79-72127	9" x 4.5"	L





PROCARE®

Forearm Splint

Padded aluminum splint is ideal for support of colles injuries prior to casting or surgery. Ventilation holes facilitate air circulation and ensure patient comfort.

Indication: Immobilization of the wrist and forearm in the initial short and/or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-71972	79-71982	7" x 2.25"	XS
79-71973	79-71983	7.5" x 2.75"	S
79-71975	79-71985	9.5" x 3.75"	M
79-71977	79-71987	10.5" x 4.25"	L



PROCARE®

Metacarpal Splints

Pre-formed, padded aluminum splint. Provides proper support for phalangeal, ulna, radius, and metacarpal fractures.

Indication: Immobilization of the wrist in the initial short and/or most severe stage of injury.

Padded Left	Padded Right	Length x Width	Size
79-72203	79-72213	8" x 3"	S
79-72205	79-72215	10.5" x 3.75"	M
79-72207	79-72217	11.5" x 4.5"	L



PROCARE®

Baseball Splints

Padded aluminum splint is easily contoured for patient comfort.

Indication: Helps maintain stability by immobilizing phalangeal joints.

Padded	Length	Qty.	Size
79-72093	4"	12 pk.	S
79-72095	4.5"	12 pk.	M
79-72097	5"	12 pk.	L



PROCARE®

Finger Protectors

Lightweight aluminum, foam padded four prong splint protects injuries to fingertips. Easily molded to fit patient.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length	Qty.	Size
79-71883	1.5"	12 pk.	S
79-71884	2.5"	12 pk.	S/M
79-71885	3"	12 pk.	M
79-71887	3.5"	12 pk.	L



PROCARE®

Posterior Tibia/Fibula Splint

Aluminum pre-formed splint is designed to aid in treatment of tibia/fibula or ankle fractures. Ventilation holes facilitate air circulation to help ensure patient comfort.

Indication: For immobilization of the foot, ankle, and/or the distal 1/3 of the tibia and fibula in the acute stage.

Padded	Unpadded	Length x Width	Size
79-72283	79-72313	11" x 7"	S
79-72285	79-72315	13.5" x 8"	M
79-72287	79-72317	16" x 9.5"	L



PROCARE®

Frog Splints

Easily contours for proper fit without bulky dressing. Padded aluminum splint helps maintain phalangeal joints in position.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length x Width	Qty.	Size
79-71963	2.75" x 2.25"	12 pk.	S
79-71965	3.15" x 2.5"	12 pk.	M
79-71967	3.75" x 3"	12 pk.	L



PROCARE®

Finger Cots

Lightweight, padded aluminum protects injuries to fingertips and designed for ease of application and removal.

Indication: Ideal for use when total finger immobilization is required.

Padded	Length	Qty.	Size
79-71903	1.5"	12 pk.	S
79-71	2.25"	12 pk.	M
79-71	3.25"	12 pk.	L
79-71	3.75"	12 pk.	XL



PROCARE®

Curved Fingers

Indicated for transverse fractures of the phalanges. Lightweight aluminum is curved to conform to finger in semi-flexed position.

Indication: Ideal for dislocations; phalangeal and metacarpal fractures.

Padded	Length x Width	Qty.	Size
79-71923	1.5" x .875"	12 pk.	S
79-71925	3" x .875"	12 pk.	M
79-71924	4" x .875"	12 pk.	M/Long
79-71926	5.5" x .875"	12 pk.	M/X-Long
79-71927	6" x .875"	12 pk.	L
79-71928	9" x .875"	12 pk.	XL



PROCARE®

Plastic Finger Splint Starter Set

30-piece assorted Plastic Finger Splints in choice of clear or skin tone colors.

Indication: Fingertip and nail bed injuries; mallet finger.

	Description	Size
79-72250	30 piece set, Skin tone	Univ.
79-72260	30 piece set, Clear	Univ.
79-72240	25 pk., Finger Straps	Univ.



PROCARE®

Plastic Finger Splints

Designed to aid in the repair of ruptured tendons. Lightweight clear or skin tone.

Indication: Fingertip and nail bed injuries; mallet finger.

Skin Tone	Clear	Length	Diameter	Qty.	Size
79-72241	79-72251	2.08"	.6"	12 pk.	#1
79-72242	79-72252	2.04"	.66"	12 pk.	#2
79-72243	79-72253	2.3"	.74"	12 pk.	#3
79-72244	79-72254	2.475"	.75"	12 pk.	#4
79-72245	79-72255	2.31"	.75"	12 pk.	#5
79-72246	79-72256	2.4"	.9"	12 pk.	#5.5
79-72247	79-72257	2.53"	.9"	12 pk.	#6
79-72248	79-72258	2.47"	1"	12 pk.	#7
79-72240	Finger Straps			25 pk.	



PROCARE®

Plastalume™ Finger Splints and Kit

Coated aluminum with foam backing. 5/8" width. Kit contains six (6) of each size splint.

Indication: Applied to the extensor surface of the finger for stability.

	Type of End	Length	Qty.
79-73211	Splint Kit		48 pcs.
79-73212	Plain	1.25"	6 pk.
79-73213	Plain	2.25"	6 pk.
79-73214	Plain	3.25"	6 pk.
79-73215	Ball	3.25"	6 pk.
79-73216	Ball	4.25"	6 pk.
79-73217	Ball	5.25"	6 pk.
79-73218	Ball	6.25"	6 pk.
79-73219	Ball	7.25"	6 pk.



PROCARE®

Aluminum Finger Splint Set

Constructed of malleable aluminum foam. 48 total piece kit contains three each of each style; Sizes small, medium and large. Frog splint, Finger Cot, Curved Finger, Baseball Splint; Finger Protectors.

Indication: Ideal for finger phalange fracture, finger joint dislocation, finger tip injury, distal finger phalange fracture distal finger joint injury, and post partial finger amputation.

	Qty.	Size
79-71020	48 pcs.	Univ.



PROCARE®

Ulnar Fracture Brace

Circumferential soft tissue compression provides cast-like rigidity for optimal control of ulnar segments. Full range of wrist and elbow motion is permitted to help prevent atrophy, increase vascularity and encourage bone growth. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Lightweight, foam lined brace is designed for the management of mid-shaft and distal ulnar fractures.

	Distal Circ.	Proximal Circ.	Size
79-97943	6" – 7"	8" – 11"	S
79-97945	6.75" – 7.75"	9.5" – 12.5"	M
79-97947	7.5" – 8.5"	11" – 14"	L



PROCARE®

Tibial Fracture Brace

A bilateral upright foot component to add medial/lateral support as well as to prevent distal slippage. Each brace comes complete with two (2) double thick cotton stockinettes. Measure length from superior border of patella to mid-level malleoli.

Indication: Foam lined preformed Sarmiento style patella tendon-bearing brace provides intimate contact needed for effective support of tibial fractures.

Right	Left	Calf Circ.	Length	Size
79-97703	79-97693	11" – 13"	15" – 16.5"	S
79-97704	79-97694	13" – 15"	18" – 20"	M/Long
79-97705	79-97695	13" – 15"	16.5" – 18.5"	M
79-97707	79-97697	15" – 18"	18" – 20"	L
79-97706	79-97696	15" – 18"	19.5" – 21.5"	L/Long



PROCARE®

Humeral Cuff

Provides total contact circumferential soft tissue compression for optimal control of humeral segments while allowing range of motion of both shoulder and elbow joints. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Prefabricated foam lined orthosis designed for the management of mid-shaft and distal third humeral fractures.

	Bicep Circum.	Length	Size
79-99213	8" – 11"	6"	S
79-99215	10" – 13"	6.25"	M
79-99217	12" – 15"	7.25"	L
79-99218	14" – 17"	8"	XL



PROCARE®

Humeral Fracture Brace/Over the Shoulder

Provides total contact circumferential soft tissue compression for optimal control of humeral segments while allowing range of motion of both shoulder and elbow joints. Features a deltoid cap extension, which provides greater alignment control and prevents distal slippage. Neutral design fits right or left. Each brace comes complete with two (2) double thick cotton stockinettes.

Indication: Prefabricated foam lined orthosis designed for the management of mid-shaft and distal third humeral fractures.

	Bicep Circum.	Length	Size
79-97953	8" – 11"	10"	S
79-97955	10" – 13"	10.25"	M
79-97957	12" – 15"	11.25"	L
79-97958	14" – 17"	12.75"	XL



PROCARE®

Posterior Leaf Splint (AFO)

Durable, orthopedic grade thermoplastic may be custom trimmed and reshaped. Will fit inside most lace-up shoes.

Indication: Designed to assist ambulation for those patients suffering mild to moderate neurological deficit resulting in drop foot.

Shoe Size		Men	Women	Height	Size
Rigth	Left				
79-97743	79-97733	6-8	≤6	12.0"	S
79-97745	79-97735	8.5-10	6.5-8	13.5"	M
79-97747	79-97737	10.5-12	8.5-10	13.5"	L
79-97748	79-97738	12+	10+	14.0"	XL

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-71020	Enovis PROCARE FINGER SPLINT KIT, 48 PCS	\$149.56
79-71883	Enovis PROCARE FINGER PROT/PADDED,S,12 PK	\$31.12
79-71884	Enovis PROCARE FINGER PROT/PADDED,S/M,12 PK	\$31.12
79-71885	Enovis PROCARE FINGER PROT/PADDED,M, 12 PK	\$31.12
79-71887	Enovis PROCARE FINGER PROT/PADDED,L,12 PK	\$31.12
79-71903	Enovis PROCARE FINGER COT/PADDED,S,12 PK	\$17.98
79-71905	Enovis PROCARE FINGER COT/PADDED,M,12 PK	\$17.98
79-71907	Enovis PROCARE FINGER COT/PADDED,L,12 PK	\$17.98
79-71908	Enovis PROCARE FINGER COT/PADDED,XL,12 PK	\$17.98
79-71923	Enovis PROCARE CRVD FNGR/PADDED,1 1/2",S,12PK	\$17.00
79-71924	Enovis PROCARE CRVD FNGR/PADDED,4",M/L,12 PK	\$17.00
79-71925	Enovis PROCARE CRVD FNGR/PADDED,3",M,12 PK	\$17.00
79-71926	Enovis PROCARE CRVD FNGR/PADDED,5 1/2",12PK	\$17.00
79-71927	Enovis PROCARE CRVD FNGR/PADDED,6",L,12 PK	\$17.00
79-71928	Enovis PROCARE CURVED FNGR/PADDED,9",XL,12 PK	\$17.00
79-71963	Enovis PROCARE FROG/PADDED,S,12 PK	\$19.81
79-71965	Enovis PROCARE FROG/PADDED,M,12 PK	\$19.81
79-71967	Enovis PROCARE FROG/PADDED,L,12 PK	\$19.81
79-71972	Enovis PROCARE FOREARM,PADDED,LT,XS	\$8.13
79-71973	Enovis PROCARE FOREARM,PADDED,LT,SM	\$8.13
79-71975	Enovis PROCARE FOREARM,PADDED,LT,MD	\$8.13
79-71977	Enovis PROCARE FOREARM,PADDED,LT,LG	\$8.13
79-71982	Enovis PROCARE FOREARM,PADDED,RT,XS	\$8.13
79-71983	Enovis PROCARE FOREARM,PADDED,RT,SM	\$8.13
79-71985	Enovis PROCARE FOREARM,PADDED,RT,MD	\$8.13
79-71987	Enovis PROCARE FOREARM,PADDED,RT,LG	\$8.13
79-72093	Enovis PROCARE BASEBALL SPL/PADDED,S,12 PK	\$28.18
79-72095	Enovis PROCARE BASEBALL SPL/PADDED,M,12 PK	\$28.18
79-72097	Enovis PROCARE BASEBALL SPL/PADDED,L,12 PK	\$28.18
79-72112	Enovis PROCARE COLLES - XS/LT. -PADDED	\$8.46
79-72113	Enovis PROCARE COLLES - SM/LT. -PADDED	\$8.46
79-72115	Enovis PROCARE COLLES - MED/LT -PADDED	\$8.46
79-72117	Enovis PROCARE COLLES - LG/LT. - PADDED	\$8.46
79-72122	Enovis PROCARE COLLES - XS/RT. - PADDED	\$8.46
79-72123	Enovis PROCARE COLLES - SM/RT. - PADDED	\$8.46
79-72125	Enovis PROCARE COLLES - MED/RT -PADDED	\$8.46
79-72127	Enovis PROCARE COLLES - LG/RT. -PADDED	\$8.46
79-72203	Enovis PROCARE METACARPAL,PADDED,LT,SM	\$9.70

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-72205	Enovis PROCARE METACARPAL,PADDED,LT,MD	\$9.70
79-72207	Enovis PROCARE METACARPAL,PADDED,LT,LG	\$9.70
79-72213	Enovis PROCARE METACARPAL - SM/RT PADDED	\$9.70
79-72215	Enovis PROCARE METACARPAL - MED/RT PADDED	\$9.70
79-72217	Enovis PROCARE METACARPAL - LG/RT PADDED	\$9.70
79-72240	Enovis PROCARE STAXX STRAPS (25-PK)	\$12.59
79-72241	Enovis PROCARE STAXX SPLINTS #1/12PK	\$27.82
79-72242	Enovis PROCARE STAXX SPLINTS #2/12PK	\$27.82
79-72243	Enovis PROCARE STAXX SPLINTS #3/12PK	\$27.82
79-72244	Enovis PROCARE STAXX SPLINTS #4/12PK	\$27.82
79-72245	Enovis PROCARE STAXX SPLINTS #5/12PK	\$27.82
79-72246	Enovis PROCARE STAXX SPLINTS #5.5/12PK	\$27.82
79-72247	Enovis PROCARE STAXX SPLINTS #6/12PK	\$27.82
79-72248	Enovis PROCARE STAXX SPLINTS #7/12PK	\$27.82
79-72250	Enovis PROCARE STAXX FINGER SPLINT KIT	\$69.99
79-72251	Enovis PROCARE CLEAR STAXX SPLINTS#1/12PK IMMOBILIZATION FINGER	\$27.82
79-72252	Enovis PROCARE CLEAR STAXX SPLINTS #2/12-PK.	\$27.82
79-72253	Enovis PROCARE CLEAR STAXX SPLINTS #3/12-PK.	\$27.82
79-72254	Enovis PROCARE CLEAR STAXX SPLINTS #4/12-PK.	\$27.82
79-72255	Enovis PROCARE CLEAR STAXX SPLINTS #5/12-PK.	\$27.82
79-72256	Enovis PROCARE CLEAR STAXX SPLINTS #5.5/12-PK	\$27.82
79-72257	Enovis PROCARE CLEAR STAXX SPLINTS #6/12-PK.	\$27.82
79-72258	Enovis PROCARE CLEAR STAXX SPLINTS #7/12-PK.	\$27.82
79-72260	Enovis PROCARE CLEAR STAXX FINGER SPLINT KIT	\$69.65
79-72283	Enovis PROCARE POST,TIBIA/FIBULA S,PADDED	\$68.66
79-72285	Enovis PROCARE POST. TIBIA/FIBULA MED/PADDED	\$68.66
79-72287	Enovis PROCARE POST. TIBIA/FIBULA LG/PADDED	\$68.66
79-73211	Enovis PROCARE PLASTALUME,SPLINTS,48 PK	\$113.40
79-73212	Enovis PROCARE PLASTALUME,1 1/4",PLN END,6 PK	\$8.76
79-73213	Enovis PROCARE PLASTALUME,2 1/4",PLN END,6 PK	\$8.76
79-73214	Enovis PROCARE PLASTALUME,3 1/4",PLN END,6 PK	\$9.70
79-73215	Enovis PROCARE PLASTALUME,3 1/4",BALL END,6PK	\$11.54
79-73216	Enovis PROCARE PLASTALUME,4 1/4",BALL END,6PK	\$13.54
79-73217	Enovis PROCARE PLASTALUME,5 1/4",BALL END,6PK	\$16.51
79-73218	Enovis PROCARE PLASTALUME,6 1/4",BALL END,6PK	\$18.52
79-73219	Enovis PROCARE PLASTALUME,7 1/4",BALL END,6PK	\$20.29
79-73220	Enovis PROCARE PLASTALUME, 4.25",PLN END, 6PK	\$15.04
79-97693	Enovis PROCARE TIBIAL FRACTURE BRACE,LT,S	\$573.52

FSS Contract # V797D-50479 | SAM UEI: SN3TFMPA2TH6

PART#	DESCRIPTION	PRICE
79-97694	Enovis PROCARE TIBIAL FRACTURE BRACE,LT,M-LG	\$573.52
79-97695	Enovis PROCARE TIBIAL FRACTURE BRACE,LT,M	\$573.52
79-97696	Enovis PROCARE TIBIAL FRACTURE BRCE,LT,LG/LNG	\$573.52
79-97697	Enovis PROCARE TIBIAL FRACTURE BRACE,LT,L	\$573.52
79-97703	Enovis PROCARE TIBIAL FRACTURE BRACE,RT,S	\$573.52
79-97704	Enovis PROCARE TIBIAL FRACTURE BRACE,RT,M-LNG	\$573.52
79-97705	Enovis PROCARE TIBIAL FRACTURE BRACE,RT,M	\$573.52
79-97706	Enovis PROCARE TIBIAL FRACTURE BRCE,RT,LG/LNG	\$573.52
79-97707	Enovis PROCARE TIBIAL FRACTURE BRACE,RT,L	\$573.52
79-97733	Enovis PROCARE POSTERIOR LEAF SPLINT,LT,S	\$126.36
79-97735	Enovis PROCARE POSTERIOR LEAF SPLINT,LT,M	\$126.36
79-97737	Enovis PROCARE POSTERIOR LEAF SPLINT,LT,L	\$126.36
79-97738	Enovis PROCARE POSTERIOR LEAF SPLINT,LT,XL	\$126.36
79-97743	Enovis PROCARE POSTERIOR LEAF SPLINT,RT,S	\$126.36
79-97745	Enovis PROCARE POSTERIOR LEAF SPLINT,RT,M	\$126.36
79-97747	Enovis PROCARE POSTERIOR LEAF SPLINT,RT,L	\$126.36
79-97748	Enovis PROCARE POSTERIOR LEAF SPLINT,RT,XL	\$126.36
79-97943	Enovis PROCARE ULNAR FRACTURE BRACE, SMALL	\$75.07
79-97945	Enovis PROCARE ULNAR FRACTURE BRACE, MEDIUM	\$75.07
79-97947	Enovis PROCARE ULNAR FRACTURE BRACE, LARGE	\$75.07
79-97953	Enovis PROCARE OVER-THE-SHLDR HUMERAL CUFF,S	\$115.28
79-97955	Enovis PROCARE OVER-THE-SHLDR HUMERAL CUFF,M	\$115.28
79-97957	Enovis PROCARE OVER-THE-SHLDR HUMERAL CUFF,L	\$115.28
79-97958	Enovis PROCARE OVER-THE-SHLDR HUMERAL CUF,XL	\$115.28
79-99213	Enovis PROCARE HUMERAL CUFF,S	\$51.80
79-99215	Enovis PROCARE HUMERAL CUFF,M	\$51.80
79-99217	Enovis PROCARE HUMERAL CUFF,L	\$51.80
79-99218	Enovis PROCARE HUMERAL CUFF,XL	\$51.80